

TECHNICAL DATA SHEET

Part Number	38310-1EA-HON
Description	Silver nitrate concentrate

Technical specifications :

- ✓ silver nitrate concentrate, for 1l standard solution, 0,1 m agno₃ (0,1n)
- ✓ CAS N°7761-88-8
- ✓ Molecular weight 169.87 g/mol
- ✓ Molecular formula AgNO₃
- ✓ no data available
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Ampoule
- ✓ Package size 1 Each

- ✓ GHS pictograms GHS09 ; GHS05
- ✓ H Statements H290 ; H314 ; H410
- ✓ P Statements P260 ; P280

Sales unit :

1 Each Silver nitrate concentrate
