

TECHNICAL DATA SHEET

Part Number	38741-1EA-HON
Description	Buffer concentrate pH 2.00

Technical specifications :

- ✓ buffer concentrate ph 2,00, for 500ml buffer sol. citric acid/hcl/nacl
- ✓ CAS N°5949-29-1
- ✓ no data available
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Ampoule
- ✓ Package size 1 Each

- ✓ GHS pictograms GHS05
- ✓ H Statements H290
- ✓ P Statements P234 ; P280

Sales unit :

1 Each Buffer concentrate pH 2.00
