

TECHNICAL DATA SHEET

Part Number	38742-1EA-HON
Description	Buffer concentrate pH 3.00

Technical specifications :

- ✓ buffer concentrate ph 3,00, for 500ml buffer sol. citric acid/naoh/nacl
- ✓ CAS N°100-56-1
- ✓ Molecular weight no data available
- ✓ no data available
- ✓ Flash point no data available
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Ampoule
- ✓ Package size 1 Each

- ✓ P Statements P280

Sales unit :

1 Each Buffer concentrate pH 3.00
