

# TECHNICAL DATA SHEET

---

<b>Part Number</b>	38744-1EA-HON
<b>Description</b>	Buffer concentrate pH 5.00

---

## Technical specifications :

- ✓ buffer concentrate pH 5,00, for 500ml buffer sol. citric acid/naoh
- ✓ CAS N°144-33-2
- ✓ Molecular weight no data available
- ✓ no data available
- ✓ Flash point Not applicable
- ✓ Material form liquid
  
- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Ampoule
- ✓ Package size 1 Each
  
- ✓ P Statements P280

## Sales unit :

1 Each Buffer concentrate pH 5.00

---