

TECHNICAL DATA SHEET

| | |
|--------------------|----------------------------|
| Part Number | 38746-1EA-HON |
| Description | Buffer concentrate pH 7.00 |

Technical specifications :

- ✓ buffer concentrate pH 7,00, for 500ml buffer sol., KH_2PO_4 / Na_2HPO_4
- ✓ no data available
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Ampoule
- ✓ Package size 1 Each

- ✓ P Statements P280

Sales unit :

1 Each Buffer concentrate pH 7.00
