

TECHNICAL DATA SHEET

Part Number	38748-1EA-HON
Description	Buffer concentrate pH 9.00

Technical specifications :

- ✓ buffer concentrate ph 9,00, for 500ml buffer sol. sodium tetraborate/hcl
- ✓ CAS N°7647-01-0
- ✓ no data available
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Ampoule
- ✓ Package size 1 Each

- ✓ P Statements P280

Sales unit :

1 Each Buffer concentrate pH 9.00
