

TECHNICAL DATA SHEET

Part Number	38750-1EA-HON
Description	Buffer concentrate pH 11.00

Technical specifications :

- ✓ buffer concentrate ph 11,00, for 500ml buffer sol. h3bo3/naoh sol./kcl
- ✓ CAS N°7447-40-7
- ✓ no data available
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Ampoule
- ✓ Package size 1 Each

- ✓ GHS pictograms GHS05 ; GHS08 ; GHS05 ; GHS08
- ✓ H Statements H360FD ; H314 ; H314 ; H360FD
- ✓ P Statements P201 ; P280 ; P201 ; P260 ; P280 ; P284

Sales unit :

1 Each Buffer concentrate pH 11.00
