

TECHNICAL DATA SHEET

Part Number	38751-1EA-HON
Description	Buffer concentrate pH 12.00

Technical specifications :

- ✓ buffer concentrate ph 12,00, for 500ml buffer sol., na2hpo4 / naoh sol.
- ✓ CAS N°1310-73-2
- ✓ no data available
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Ampoule
- ✓ Package size 1 Each

- ✓ GHS pictograms GHS07
- ✓ H Statements H315 ; H319
- ✓ P Statements P280

Sales unit :

1 Each Buffer concentrate pH 12.00
