

TECHNICAL DATA SHEET

Part Number	38752-1EA-HON
Description	Buffer concentrate pH 13.00

Technical specifications :

- ✓ buffer concentrate ph 13,00, for 500ml buffer sol. glycine/naoh/nacl
- ✓ CAS N°7647-14-5
- ✓ no data available
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Ampoule
- ✓ Package size 1 Each

- ✓ GHS pictograms GHS05
- ✓ H Statements H290 ; H314
- ✓ P Statements P234 ; P260 ; P280 ; P284

Sales unit :

1 Each Buffer concentrate pH 13.00
