



PETRI DISHES

Placas de petri / Boîtes de petri / Piastre di petri

PETRI DISH TRIPLE VENTS

placa de petri con tres vientos / boîte de petri avec trois ergots / piastra di petri tripla ventilazione

MAXIMUM TRANSPARENCY

máxima transparencia / transparence maximale / massima trasparenza


PERFECT PLANARITY

planaridad perfecta / parfaite planéité / planarità perfetta



POLYSTYRENE

poliestireno / polystyrène / polistirolo



	CODE	TYPE	STERILITY	PACKAGING
	29042	ø 90 mm, H 14,2 mm	PRODUCED IN CLEANROOM ISO 6 producido en cleanroom ISO 6 / produit en cleanroom ISO 6 / prodotta in camera bianca ISO 6	BAG OF 20 PIECES bolsa de 20 piezas / en sachet de 20 pièces / in sacchetto da 20 pezzi
	29082	ø 90 mm, H 14,2 mm	STERILE estéril / stérile / sterile	BAG OF 20 PIECES bolsa de 20 piezas / en sachet de 20 pièces / in sacchetto da 20 pezzi
	29052	ø 90 mm, H 16,2 mm	PRODUCED IN CLEANROOM ISO 6 producido en cleanroom ISO 6 / produit en cleanroom ISO 6 / prodotta in camera bianca ISO 6	BAG OF 20 PIECES bolsa de 20 piezas / en sachet de 20 pièces / in sacchetto da 20 pezzi
	29062	ø 90 mm, H 16,2 mm	STERILE estéril / stérile / sterile	BAG OF 20 PIECES bolsa de 20 piezas / en sachet de 20 pièces / in sacchetto da 20 pezzi