



DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY DECLARATION DE CONFORMITÈ/ KONFORMITÄTSEKKLÄRUNG

Nome e indirizzo della ditta
Name and address of the firm
Name und Adresse der Firma
Nom et adresse de l'entreprise

BIOSIGMA S.r.l.
Via Valletta, 6 – Zona Ind. Cantarana
30010 - Cona (VE) - ITALY
Tel. +39.0426.302224 - Fax +39.0426.302228
E-mail: info@biosigmaeu.com
<http://www.biosigma.com>

Dichiaro sotto nostra responsabilità che il dispositivo medico-diagnostico in vitro:

We declare under our sole responsibility that the in vitro diagnostic medical device:
Wir erklären in alleiniger Verantwortung, dass das Medizinprodukt für die In-vitro-Diagnostik:
Nous déclarons sous notre propre responsabilité que le dispositif médical de diagnostique in vitro:

Vedi elenco allegato / *See the enclosed list*
Produkt der Liste beiliegend / Voir liste ci-joint /

Classificazione secondo la direttiva relativa ai dispositivi medico-diagnostici in vitro 98/79/CE.
Classification according to the directive on in vitro diagnostic medical devices 98/79/EC.
Klassifizierung nach der Richtlinie über In-vitro-Diagnostika 98/79/EC.
Classification selon la directive relative aux dispositifs médicaux de diagnostic in vitro 98/79/CE.

- Dispositivo dell'elenco A, allegato II / *Produkt der Liste A, Anhang II / Dispositif de la liste B, annexe II / Device of List A, Annex II*
- Dispositivo dell'elenco B, allegato II / *Device of List B, Annex II / Produkt der Liste B, Anhang II / Dispositif de la liste B, annexe II*
- Dispositivo per test autodiagnostico non elencato nell'allegato II / *Device for self-testing not listed in Annex II / Produkt zur Eigenanwendung, das nicht in Anhang II genannt ist / Dispositif destiné à l'autodiagnostic non listé dans l'annexe II*
- Altro dispositivo / Other device / Sonstiges Produkt / Autre dispositif**

- **soddisfa tutte le disposizioni della direttiva relativa ai dispositivi medico-diagnostici in vitro 98/79/CE che lo riguardano e risulta conforme ai requisiti essenziali di cui all'allegato I della direttiva 98/79/CE recepita con il D.Lgs 332 del 08/09/2000.**
- *meets all the provisions of the directive on in vitro diagnostic medical devices 98/79/EC which apply to it and is manufactured in compliance with essential requirements of Annex 1 of the 98/79/CE Directive.*
- *allen Anforderungen der Richtlinie über In-vitro-Diagnostika 98/79/EC entspricht, die anwendbar sind und es erfüllt die grundlegenden Anforderungen des Anhangs I 98/79 / EG Richtlinie.*
- *remplit toutes les exigences de la directive relative aux dispositifs médicaux de diagnostic in vitro 98/79/CE qui le concernent et il est conforme aux exigences essentielles de l'annexe I de la directive 98/79/CE.*

Procedimento di valutazione della conformità:
Conformity assessment procedure:
Konformitätsbewertungsverfahren:
Procédure d'évaluation de la conformité:

Allegato III (escluso punto 6)
Annex III (excluding point 6)
Anhang III (mit Ausnahme von Punkt 6)
Annexe III (à l'exclusion, point 6)

Organo incaricato della valutazione della conformità (se consultato):
Notified Body (if consulted):
Konformitätsbewertungsstelle (falls beigezogen) :
Organe resp. de l'évaluat. de la conform.(si consult) :

N.A.

CONA 22/06//2018
Luogo, data / Place, date / Ort, Datum / Lieu, date

MARCOLIN MARTINO – Consigliere Delegato

Nome e funzione / Name and function / Name und Funktion / Nom et fonction



Biosigma s.r.l.
a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY DECLARATION DE CONFORMITÉ/ KONFORMITÄTSEKTLÄRUNG

Le DICHIARAZIONI DI CONFORMITA' alla Direttiva 98/79/CE relative alle varie famiglie di prodotti si riferiscono ai seguenti codici in dettaglio:

The DECLARATIONS OF CONFORMITY to Directive 98/79/EC concerning the various families of products make reference to following code number

REF	ANSE	MICROLOOPS
BSV120	Ansa 10 ul in PP STERILE –in Peel Pack da 20 pezzi	Microloop 10 µl in PP STERILE in Peel Pack bag of 20 Pcs.
BSV1200	Ansa 10 ul in PP STERILE –in Peel Pack da 10 pezzi	Microloop 10 µl in PP STERILE in Peel Pack bag of 10 Pcs.
BSV1201	Ansa 10 ul in PP STERILE –in confezione singola	Microloop 10 µl in PP STERILE individually wrapped.
BSV1205	Ansa 10 ul in PP STERILE –in Peel Pack da 5 pezzi	Microloop 10 µl in PP STERILE in Peel Pack bag of 5 Pcs
BSV121	Ansa 1 ul in PP STERILE in Peel Pack da 20 pezzi	Microloop 1 µl in PP STERILE in Peel Pack bag of 20 Pcs.
BSV1210	Ansa 1 ul in PP STERILE in Peel Pack da 10 pezzi	Microloop 1 µl in PP STERILE in Peel Pack bag of 10 Pcs.
BSV1211	Ansa 1 ul in PP STERILE in confezione singola	Microloop 1 µl in PP STERILE individually wrapped.
BSV1215	Ansa 1 ul in PP STERILE in Peel Pack da 5 pezzi	Microloop 1 µl in PP STERILE in Peel Pack bag of 5 Pcs.
REF	CASSETTE BIOPSIA /ISTOLOGIA	BIOPSY/HISTOLOGY CASSETTES
BSI4401	Anello per inclusione misura piccola BIANCO	Embedding ring small size WHITE
BSI102	Cassetta per Istologia BIANCA	Histology cassette WHITE
BSI1021	Cassetta per Istologia GRIGIA	Histology cassette GRAY
BSI1022	Cassetta per Istologia GIALLA	Histology cassette YELLOW
BSI1024	Cassetta per Istologia ROSA	Histology cassette PINK
BSI1026	Cassetta per Istologia VERDE	Histology cassette GREEN
BSI1027	Cassetta per Istologia ARANCIO	Histology cassette ORANGE
BSI1028	Cassetta per Istologia AZZURRA	Histology cassette LIGHT BLUE
BSI103	Cassette per Biopsia BIANCO	Biopsy cassette WHITE
BSI1031	Cassetta per Biopsia GRIGIA	Biopsy cassette GRAY
BSI1032	Cassette per Biopsia GIALLA	Biopsy cassette YELLOW
BSI1034	Cassette per Biopsia ROSA	Biopsy cassette PINK
BSI1036	Cassette per Biopsia VERDE	Biopsy cassette GREEN
BSI1037	Cassette per Biopsia ARANCIO	Biopsy cassette ORANGE
BSI1038	Cassette per Biopsia AZZURRA	Biopsy cassette LIGHT BLUE
BSI104	Cassetta per inclusione senza coperchio	Embedding cassette without lid
REF	CONTENITORI	CONTAINERS
BSC253	CONTENITORE 120ml PER URINA IN PP TAPPO A VITE BIANCO	URINE CONTAINER 120ml IN PP WITH WHITE SCREW CAP
BSC258	CONTENITORE 120ml PER URINA IN PP TAPPO A VITE ROSSO CONFEZIONE SINGOLA STERILE	URINE CONTAINER 120ml IN PP WITH RED SCREW CAP INDIVIDUALLY WRAPPED STERILE



Biosigma s.r.l.
a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



**DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY
DECLARATION DE CONFORMITÉ/ KONFORMITÄTSEKTLÄRUNG**

BSC660	CONTENITORE 60 ml PER URINA in PP TAPPO A VITE BIANCO	<i>URINA CONTAINER 60ml in PP WITH WHITE SCREW CAP</i>
BSC665	CONTENITORE 60 ml, TAPPO A VITE ROSSO, CON AREA SCRIVIBILE – STERILE in confezione singola	<i>CONTAINER 60ml , RED SCREW CAP WITH FROSTED SIDE INDIVIDUALLY WRAPPED - STERILE</i>
REF	PRODOTTI PER CRIOGENIA	CRYOGENICS' PRODUCTS
10001TS	Tappo esterno MARIA STELLA	<i>MARIA STELLA external cap</i>
10011TS	Tappo interno MARIO STELLA	<i>MARIO STELLA internal cap</i>
1BETS	Provetta CryoGen da 1,2 ml con base, per tappo Stella esterno	<i>CryoGen Tube 1.2 ml self standing, for external Stella cap</i>
2BETS	Provetta CryoGen da 2 ml con base, per tappo Stella esterno	<i>CryoGen Tube 2 ml self standing, for external Stella cap</i>
2BI	Provetta CryoGen 2ml con base, filetto interno	<i>CryoGen Tube 2ml self standing, internal thread</i>
2ETS	Provetta CryoGen da 2 ml, per tappo Stella esterno	<i>CryoGen Tube 2 ml, for external Stella cap</i>
2I	Provetta CryoGen 2ml filetto interno	<i>CryoGen Tube 2ml internal thread</i>
3BETS	Provetta CryoGen da 3 ml con base, per tappo Stella esterno	<i>CryoGen Tube 3 ml self standing, for external Stella cap</i>
4BETS	Provetta CryoGen da 4 ml con base, per tappo Stella esterno	<i>CryoGen Tube 4 ml self standing, for external Stella cap</i>
4BI	Provetta CryoGen da 4 ml con base, filetto interno	<i>CryoGen Tube 4ml self standing, internal thread</i>
4I	Provetta CryoGen da 4ml, filetto interno	<i>CryoGen Tube 4ml, internal thread</i>
5BETS	Provetta CryoGen da 5 ml con base, per tappo Stella esterno	<i>CryoGen Tube 5 ml self standing, for external Stella cap</i>
5I	Provetta CryoGen da 5ml filetto interno	<i>CryoGen Tube 5ml internal thread</i>
BSM58025B/8X	CRYOGEN® BOX25 per provette da 1,2 e 2 ml, fondo bianco, griglia vari colori - conf.8 pz	<i>CRYOGEN® BOX25 for 1.2 and 2 ml tubes, white bottom, coloured grid - conf.8 pcs</i>
BSM58025B/X	CRYOGEN® BOX25 per provette da 1,2 e 2 ml, fondo bianco, griglia vari colori	<i>CRYOGEN® BOX25 for 1.2 and 2 ml tubes, white bottom, coloured grid</i>
BSM58025B1D/8X	CRYOGEN® BOX25 per provette da 1,2 e 2 ml, fondo bianco con barcode, griglia vari colori - conf.8 pz	<i>CRYOGEN® BOX25 for 1.2 and 2 ml tubes, white bottom with barcode, coloured grid - conf.8 pcs</i>
BSM58025B1D/X	CRYOGEN® BOX25 per provette da 1,2 e 2 ml, fondo bianco con barcode, griglia vari colori	<i>CRYOGEN® BOX25 for 1.2 and 2 ml tubes, white bottom with barcode, coloured grid</i>
BSM580612D/5X	CRYOGEN® BOX81-2D per provette da 2 ml 2D, griglia vari colori - conf.5 pz	<i>CRYOGEN® BOX81-2D for 2 ml tubes 2D, coloured grid - conf.5 pcs</i>
BSM580612D/X	CRYOGEN® BOX81-2D per provette da 2 ml 2D, griglia vari colori	<i>CRYOGEN® BOX81-2D for 2 ml tubes 2D, coloured grid</i>
BSM58061B/5X	CRYOGEN® BOX81 per provette da 1,2 e 2 ml, fondo bianco, griglia vari colori - conf.5 pz	<i>CRYOGEN® BOX81 for 1,2 and 2 ml tubes, white bottom, coloured grid - conf.5 pcs</i>
BSM58061B/X	CRYOGEN® BOX81 per provette da 1,2 e 2 ml, fondo bianco, griglia vari colori	<i>CRYOGEN® BOX81 for 1,2 and 2 ml tubes, white bottom, coloured grid</i>
BSM58061B1D/5X	CRYOGEN® BOX81 per provette da 1,2 e 2 ml, fondo bianco con barcode, griglia vari colori - conf.5 pz	<i>CRYOGEN® BOX81 for 1,2 and 2 ml tubes, white bottom with barcode, coloured grid - conf.5 pcs</i>



Biosigma s.r.l.
a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



**DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY
DECLARATION DE CONFORMITÉ/ KONFORMITÄTSEKTLÄRUNG**

BSM58061B1D/X	CRYOGEN® BOX81 per provette da 1,2 e 2 ml, fondo bianco con barcode, griglia vari colori	<i>CRYOGEN® BOX81 for 1,2 and 2 ml tubes, white bottom with barcode, coloured grid</i>
BSM580712D/4X	CRYOGEN® BOX81-2D per provette da 5 ml 2D, griglia vari colori - conf.4 pz	<i>CRYOGEN® BOX81-2D for 5 ml tubes 2D, coloured grid - conf.4 pcs</i>
BSM580712D/X	CRYOGEN® BOX81-2D per provette da 5 ml 2D, griglia vari colori	<i>CRYOGEN® BOX81-2D for 5 ml tubes 2D, coloured grid</i>
BSM58071B/4X	CRYOGEN® BOX81 per provette da 3, 4 e 5 ml, fondo bianco, griglia vari colori - conf.4 pz	<i>CRYOGEN® BOX81 for 3, 4 and 5 ml tubes, white bottom, coloured grid - conf.4 pcs</i>
BSM58071B/X	CRYOGEN® BOX81 per provette da 3, 4 e 5 ml, fondo bianco, griglia vari colori	<i>CRYOGEN® BOX81 for 3, 4 and 5 ml tubes, white bottom, coloured grid</i>
BSM58071B1D/4X	CRYOGEN® BOX81 per provette da 3, 4 e 5 ml, fondo bianco con barcode, griglia vari colori - conf.4 pz	<i>CRYOGEN® BOX81 for 3, 4 and 5 ml tubes, white bottom with barcode, coloured grid - conf.4 pcs</i>
BSM58071B1D/X	CRYOGEN® BOX81 per provette da 3, 4 e 5 ml, fondo bianco con barcode, griglia vari colori	<i>CRYOGEN® BOX81 for 3, 4 and 5 ml tubes, white bottom with barcode, coloured grid</i>
BSM58081B/X	CRYOGEN® BOX81 per provette da 1,2 e 2 ml, griglia vari colori, coperchio inciso	<i>CRYOGEN® BOX81 for 1,2 and 2 ml tubes, coloured grid, engraved cover</i>
BSM58091B/X	CRYOGEN® BOX81 per provette da 3, 4 e 5 ml, griglia vari colori, coperchio inciso	<i>CRYOGEN® BOX81 for 3, 4 and 5 ml tubes, coloured grid, engraved cover</i>
BSM5810/X	WORKSTATION per provette CryoGen a 40 posti, vari colori	<i>WORKSTATION for CryoGen tubes 40 holes, coloured</i>
BSM58100/5X	CRYOGEN® BOX100 per provette da 1,2 e 2 ml, fondo vari colori - conf.5 pz	<i>CRYOGEN® BOX100 for 1,2 and 2 ml tubes, coloured bottom - conf.5 pcs</i>
BSM58100/X	CRYOGEN® BOX100 per provette da 1,2 e 2 ml, fondo vari colori	<i>CRYOGEN® BOX100 for 1,2 and 2 ml tubes, coloured bottom</i>
BSM581001D/B	CRYOGEN® BOX100 per provette da 1,2 e 2 ml, fondo vari colori con barcode	<i>CRYOGEN® BOX100 for 1,2 and 2 ml tubes, coloured bottom with barcode</i>
BSM581001D/5B	CRYOGEN® BOX100 per provette da 1,2 e 2 ml, fondo vari colori con barcode- conf.5 pz	<i>CRYOGEN® BOX100 for 1,2 and 2 ml tubes, coloured bottom with barcode- conf.5 pcs</i>
BSM581002D/5X	CRYOGEN® BOX100-2D per provette da 2 ml 2D, fondo vari colori - conf.5 pz	<i>CRYOGEN® BOX100-2D for 2 ml tubes 2D, coloured bottom - conf.5 pcs</i>
BSM581002D/X	CRYOGEN® BOX100-2D per provette da 2 ml 2D, fondo vari colori	<i>CRYOGEN® BOX100-2D for 2 ml tubes 2D, coloured bottom</i>
BSM581052D/4X	CRYOGEN® BOX100-2D per provette da 5 ml 2D, fondo vari colori - conf.4 pz	<i>CRYOGEN® BOX100-2D for 5 ml tubes 2D, coloured bottom - conf.4 pcs</i>
BSM581052D/X	CRYOGEN® BOX100-2D per provette da 5 ml 2D, fondo vari colori	<i>CRYOGEN® BOX100-2D for 5 ml tubes 2D, coloured bottom</i>
CL1ARBEPSTS	Provetta CryoGen® da 1,2 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 1,2 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL1ARBEPSTS/B	Provetta CryoGen® da 1,2 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 1,2 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL1ARBEPSTS/BX	Provetta CryoGen® da 1,2 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 1,2 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>



Biosigma s.r.l.
a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



**DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY
DECLARATION DE CONFORMITÉ/ KONFORMITÄTSEKTLÄRUNG**

CL1ARBEPSTS/X	Provetta CryoGen® da 1,2 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 1,2 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL1ARBIPSTS	Provetta CryoGen® da 1,2 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 1,2 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL1ARBIPSTS/B	Provetta CryoGen® da 1,2 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 1,2 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®-- IN BULK</i>
CL1ARBIPSTS/BX	Provetta CryoGen® da 1,2 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 1,2 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®-- IN BULK</i>
CL1ARBIPSTS/X	Provetta CryoGen® da 1,2 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 1,2 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARBEPS2D	Provetta CryoGen® da 2 ml con base, tappo Stella esterno, graduata, bar-code + 2D, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml self standing, external Stella cap, graduated, with bar-code + 2D, sterile - CLEARLINE®</i>
CL2ARBEPS2D/B	Provetta CryoGen® da 2 ml con base, tappo Stella esterno, graduata, bar-code + 2D, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml self standing, external Stella cap, graduated, with bar-code + 2D, sterile - CLEARLINE® - IN BULK</i>
CL2ARBEPSTS	Provetta CryoGen® da 2 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARBEPSTS/B	Provetta CryoGen® da 2 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®-- IN BULK</i>
CL2ARBEPSTS/BX	Provetta CryoGen® da 2 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®-- IN BULK</i>
CL2ARBEPSTS/X	Provetta CryoGen® da 2 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARBIPS2D	Provetta CryoGen® da 2 ml con base, tappo Stella interno, graduata, bar-code + 2D, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml self standing, internal Stella cap, graduated, with bar-code + 2D, sterile - CLEARLINE®</i>
CL2ARBIPS2D/B	Provetta CryoGen® da 2 ml con base, tappo Stella interno, graduata, bar-code + 2D, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml self standing, internal Stella cap, graduated, with bar-code + 2D, sterile - CLEARLINE® - IN BULK</i>
CL2ARBIPSTS	Provetta CryoGen® da 2 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARBIPSTS/B	Provetta CryoGen® da 2 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2ARBIPSTS/BX	Provetta CryoGen® da 2 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>



Biosigma s.r.l.
a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



**DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY
DECLARATION DE CONFORMITÉ/ KONFORMITÄTSEKTLÄRUNG**

CL2ARBIPSTS/X	Provetta CryoGen® da 2 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2AREPSTS	Provetta CryoGen® da 2 ml , tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2AREPSTS/B	Provetta CryoGen® da 2 ml , tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2AREPSTS/BX	Provetta CryoGen® da 2 ml , tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2AREPSTS/X	Provetta CryoGen® da 2 ml , tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARIPSTS	Provetta CryoGen® da 2 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARIPSTS/B	Provetta CryoGen® da 2 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®- IN BULK</i>
CL2ARIPSTS/BX	Provetta CryoGen® da 2 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®- IN BULK</i>
CL2ARIPSTS/X	Provetta CryoGen® da 2 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL3ARBEPSSN	Provetta CryoGen® da 3 ml con base, tappo a vite esterno bi-iniettato, graduata, con area di scrittura, sterile - CLEARLINE® - con bar-code a sequenza personalizzata	<i>CryoGen® Tube 3ml self standing, external bi-injected screw cap, graduated, with writing area, sterile - CLEARLINE® - with bar-code custom sequence</i>
CL3ARBEPSTS	Provetta CryoGen® da 3 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 3 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL3ARBEPSTS/B	Provetta CryoGen® da 3 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 3 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL3ARBEPSTS/BX	Provetta CryoGen® da 3 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 3 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL3ARBEPSTS/X	Provetta CryoGen® da 3 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 3 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARBEPSTS	Provetta CryoGen® da 4 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARBEPSTS/B	Provetta CryoGen® da 4 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>



Biosigma s.r.l.
a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



**DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY
DECLARATION DE CONFORMITÉ/ KONFORMITÄTSEKTLÄRUNG**

CL4ARBEPSTS/BX	Provetta CryoGen® da 4 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARBEPSTS/X	Provetta CryoGen® da 4 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARBIPSTS	Provetta CryoGen® da 4 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARBIPSTS/B	Provetta CryoGen® da 4 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARBIPSTS/BX	Provetta CryoGen® da 4 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARBIPSTS/X	Provetta CryoGen® da 4 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARIPSTS	Provetta CryoGen® da 4 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARIPSTS/B	Provetta CryoGen® da 4 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARIPSTS/BX	Provetta CryoGen® da 4 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARIPSTS/X	Provetta CryoGen® da 4 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARBEP2D	Provetta CryoGen® da 5 ml con base, tappo Stella esterno, graduata, bar-code + 2D, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, external Stella cap, graduated, with bar-code + 2D, sterile - CLEARLINE®</i>
CL5ARBEP2D/B	Provetta CryoGen® da 5 ml con base, tappo Stella esterno, graduata, bar-code + 2D, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, external Stella cap, graduated, with bar-code + 2D, sterile - CLEARLINE® - IN BULK</i>
CL5ARBEPSTS	Provetta CryoGen® da 5 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARBEPSTS/B	Provetta CryoGen® da 5 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARBEPSTS/BX	Provetta CryoGen® da 5 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARBEPSTS/X	Provetta CryoGen® da 5 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>



Biosigma s.r.l.

a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 | Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 | P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



**DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY
DECLARATION DE CONFORMITÉ/ KONFORMITÄTSEKTLÄRUNG**

CL5ARBIPS2D	Provetta CryoGen® da 5 ml con base, tappo Stella interno, graduata, bar-code + 2D, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, internal Stella cap, graduated, with bar-code + 2D, sterile - CLEARLINE®</i>
CL5ARBIPS2D/B	Provetta CryoGen® da 5 ml con base, tappo Stella interno, graduata, bar-code + 2D, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, internal Stella cap, graduated, with bar-code + 2D, sterile - CLEARLINE® - IN BULK</i>
CL5ARBIPSSN	Provetta CryoGen® da 5ml con base, tappo a vite interno bi-iniettato, graduata, con area di scrittura, sterile - CLEARLINE® - con bar-code a sequenza personalizzata	<i>CryoGen® Tube 5ml self standing, internal bi-injected screw cap, graduated, with writing area, sterile - CLEARLINE® - with bar-code custom sequence</i>
CL5ARBIPSTS	Provetta CryoGen® da 5 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARBIPSTS/B	Provetta CryoGen® da 5 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARBIPSTS/BX	Provetta CryoGen® da 5 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARBIPSTS/X	Provetta CryoGen® da 5 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARIPSTS	Provetta CryoGen® da 5 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARIPSTS/B	Provetta CryoGen® da 5 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARIPSTS/BX	Provetta CryoGen® da 5 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARIPSTS/X	Provetta CryoGen® da 5 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CRY910TS/X	CAP DISK per tappo STELLA, vari colori	<i>CAPDISK for STELLA cap / colour</i>
REF	CUVETTE	CUVETTES
BSA001	Cuvetta Macro monouso per spettrofotometria in PS	<i>Disposable spectrophotometer Macro cuvettes in PS</i>
BSA002	Cuvetta Semi-Micro monouso per spettrofotometria in PS	<i>Disposable spectrophotometer Semi-Micro cuvettes in PS</i>
BSA001M	Cuvetta Macro monouso per spettrofotometria in PMMA	<i>Disposable spectrophotometer Macro cuvettes in PMMA</i>
BSA002M	Cuvetta Semi-Micro monouso per spettrofotometria in PMMA	<i>Disposable spectrophotometer Semi-Micro cuvettes in PMMA</i>
BSA094	Cuvetta tipo Hitachi 902 Bionex	<i>Cuvette Hitachi 902 type in Bionex</i>
BSA096	Cuvetta tipo Hitachi 717 in PMMA	<i>Cuvette Hitachi 717 type in PMMA</i>
BSA097	Cuvetta tipo Hitachi 911/912 in PMMA	<i>Cuvette Hitachi 911/912 type in PMMA</i>
BSA098	Cuvetta tipo Hitachi 911/912 Bionex	<i>Cuvette Hitachi 911/912 type in Bionex</i>
BSA099	Cuvetta tipo Hitachi 704 in PMMA	<i>Cuvette Hitachi 704 type in PMMA</i>



Biosigma s.r.l.
a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



**DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY
DECLARATION DE CONFORMITÈ/ KONFORMITÄTSEKTLÄRUNG**

BSA107	Cuvetta singola in PMMA	Single cuvette in PMMA.
BSA108	Cuvetta singola in Bionex	Single cuvette in Bionex
BSA655	Cuvetta tipo Express 550 and Plus	Cuvette Express 550 and Plus type
BSA656	Cuvetta tipo Express 550 and Plus in PMMA con sacchetti di recupero	Cuvette Express 550 and Plus type in PMMA with bag
BSA84	Cuvetta multipla in PMMA	Multicells cuvette in PMMA
BSA85	Cuvetta multipla in PMMA in rack	Multicells cuvette in PMMA in rack
BSP300	Cuvetta a 12 cellette con pinze	12 cells cuvette with 2 pliers.
REF	CONTENITORI PER ESPETTORATO E LIQUIDO SEMINALE	SEMINAL FLUID & SPUTUM COLLECTORS
BSC740	FERTILITY COLLECTOR KIT- STERILE	FERTILITY COLLECTOR KIT- STERILE
REF	KITS PER L'ANALISI DELL'URINA	KITS FOR URINALYSIS
BVS101	KIT: FAST READ 102 + PIPETTE	KIT: FAST READ 102 + PIPETTES
BVS102	KIT: FAST READ 102 + PIPETTE + STAIN	KIT: FAST READ 102 + PIPETTES + STAIN
BVS171	KIT: FAST READ 102 + MT MICROTIPS	KIT: FAST READ 102 + MT MICROTIPS
BVS1715	KIT: FAST READ 102 + MT SYSTEM	KIT: FAST READ 102 + MT SYSTEM
BVS1717	KIT: FAST READ 102 + MT SYSTEM + STAIN	KIT: FAST READ 102 + MT SYSTEM + STAIN
BVS1719	KIT: FAST READ 102 + MT SYSTEM + TAPPO	KIT: FAST READ 102 + MT SYSTEM + STOPPER
BVS1720	KIT: FAST READ 102 + MT SYSTEM + TAPPO + STAIN	KIT: FAST READ 102 + MT SYSTEM + STOPPER + STAIN
BVS1721	KIT: FAST READ 102 + PROVETTA ETICHETTATA TAPPATA + MT MICROTIP	KIT: FAST READ 102 + TUBE LABELLED STOPPERED + MT MICROTIP
BVS190	MT SYSTEM	MT SYSTEM
U9028	URIGLASS KIT	URIGLASS KIT
U9029	URIGLASS KIT + STAIN	URIGLASS KIT + STAIN
REF	MICROPIPETTE	MICROPIPETTES
ABS119HP	Hp-pette volume fisso 5 µl	Hp-pette fixed volume 5 µl
ABS120HP	Hp-pette volume fisso 10 µl	Hp-pette fixed volume 10 µl
ABS1210HP	Hp-pette volume fisso 25 µl	Hp-pette fixed volume 25 µl
ABS121HP	Hp-pette volume fisso 20 µl	Hp-pette fixed volume 20 µl
ABS122HP	Hp-pette volume fisso 50 µl	Hp-pette fixed volume 50 µl
ABS123HP	Hp-pette volume fisso 100 µl	Hp-pette fixed volume 100 µl
ABS1240HP	Hp-pette volume fisso 250 µl	Hp-pette fixed volume 250 µl
ABS1241HP	Hp-pette volume fisso 500 µl	Hp-pette fixed volume 500 µl
ABS124HP	Hp-pette volume fisso 200 µl	Hp-pette fixed volume 200 µl
ABS125HP	Hp-pette volume fisso 1000 µl	Hp-pette fixed volume 1000 µl
ABS1261HP	Hp-pette volume variabile 0,1 – 2,5 µl	Hp-pette variable volume 0,1 – 2,5 µl
ABS126HP	Hp-pette volume variabile 0,5 – 10 µl	Hp-pette variable volume 0,5 – 10 µl
ABS127HP	Hp-pette volume variabile 1000 - 5000 µl	Hp-pette variable volume 1000 - 5000 µl
ABS128HP	Hp-pette volume variabile 2 -20 µl	Hp-pette variable volume 2 -20 µl
ABS1290HP	Hp-pette volume variabile 10 - 100 µl	Hp-pette variable volume 10 - 100 µl



Biosigma s.r.l.
a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



**DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY
DECLARATION DE CONFORMITÉ/ KONFORMITÄTSERKLÄRUNG**

ABS129HP	Hp-pette volume variabile 5 - 50 µl	<i>Hp-pette variable volume 5 - 50 µl</i>
ABS1310HP	Hp-pette volume variabile 50 - 200 µl	<i>Hp-pette variable volume 50 - 200 µl</i>
ABS1311HP	Hp-pette volume variabile 100 – 1000 µl	<i>Hp-pette variable volume 100 – 1000 µl</i>
ABS131HP	Hp-pette volume variabile 20 - 200 µl	<i>Hp-pette variable volume 20 - 200 µl</i>
ABS132HP	Hp-pette volume variabile 200 – 1000 µl	<i>Hp-pette variable volume 200 – 1000 µl</i>
ABS133HP	Pipette MULTICANALE a volume variabile 5-50 µl	<i>Variable volumes MULTICHANNEL Pipettes 5-50 µl</i>
ABS135HP	Pipette MULTICANALE a volume variabile 50-300 µl	<i>Variable volumes MULTICHANNEL Pipettes 50-300 µl</i>
ABS136HP	Pipette MULTICANALE a volume variabile 0,5-10 µl	<i>Variable volumes MULTICHANNEL Pipettes 0,5-10 µl</i>
ABS119HPA	Hp-pette Plus completamente autoclavabili volume fisso 5 µl	<i>Hp-pette Plus Fixed volume Fully Autoclavables 5 µl</i>
ABS120HPA	Hp-pette Plus completamente autoclavabili volume fisso 10 µl	<i>Hp-pette Plus Fixed volume Fully Autoclavables 10 µl</i>
ABS121HPA	Hp-pette Plus completamente autoclavabili volume fisso 20 µl	<i>Hp-pette Plus Fixed volume Fully Autoclavables 20 µl</i>
ABS1210HPA	Hp-pette Plus completamente autoclavabili volume fisso 25 µl	<i>Hp-pette Plus Fixed volume Fully Autoclavables 25 µl</i>
ABS122HPA	Hp-pette Plus completamente autoclavabili volume fisso 50 µl	<i>Hp-pette Plus Fixed volume Fully Autoclavables 50 µl</i>
ABS123HPA	Hp-pette Plus completamente autoclavabili volume fisso 100 µl	<i>Hp-pette Plus Fixed volume Fully Autoclavables 100 µl</i>
ABS124HPA	Hp-pette Plus completamente autoclavabili volume fisso 200 µl	<i>Hp-pette Plus Fixed volume Fully Autoclavables 200 µl</i>
ABS1240HPA	Hp-pette Plus completamente autoclavabili volume fisso 250 µl	<i>Hp-pette Plus Fixed volume Fully Autoclavables 250 µl</i>
ABS1241HPA	Hp-pette Plus completamente autoclavabili volume fisso 500 µl	<i>Hp-pette Plus Fixed volume Fully Autoclavables 500 µl</i>
ABS125HPA	Hp-pette Plus completamente autoclavabili volume fisso 1000 µl	<i>Hp-pette Plus Fixed volume Fully Autoclavables 1000 µl</i>
ABS1261HPA	Hp-pette Plus completamente autoclavabili volume variabile 0,1 – 2,5 µl	<i>Hp-pette Plus Variable volumes Fully Autoclavables 0,1 – 2,5 µl</i>
ABS126HPA	Hp-pette Plus completamente autoclavabili volume variabile 0,5 – 10 µl	<i>Hp-pette Plus Variable volumes Fully Autoclavables 0,5 – 10 µl</i>
ABS128HPA	Hp-pette Plus completamente autoclavabili volume variabile 2 - 20 µl	<i>Hp-pette Plus Variable volumes Fully Autoclavables 2 - 20 µl</i>
ABS129HPA	Hp-pette Plus completamente autoclavabili volume variabile 5 -50 µl	<i>Hp-pette Plus Variable volumes Fully Autoclavables 5 -50 µl</i>
ABS1290HPA	Hp-pette Plus completamente autoclavabili volume variabile 10 - 100 µl	<i>Hp-pette Plus Variable volumes Fully Autoclavables 10 - 100 µl</i>
ABS131HPA	Hp-pette Plus completamente autoclavabili volume variabile 20 - 200 µl	<i>Hp-pette Plus Variable volumes Fully Autoclavables 20 - 200 µl</i>
ABS1310HPA	Hp-pette Plus completamente autoclavabili volume variabile 50 - 200 µl	<i>Hp-pette Plus Variable volumes Fully Autoclavables 50 - 200 µl</i>
ABS1311HPA	Hp-pette Plus completamente autoclavabili volume variabile 100 – 1000 µl	<i>Hp-pette Plus Variable volumes Fully Autoclavables 100 – 1000 µl</i>
ABS132HPA	Hp-pette Plus completamente autoclavabili volume variabile 200 – 1000 µl	<i>Hp-pette Plus Variable volumes Fully Autoclavables 200 – 1000 µl</i>



Biosigma s.r.l.
a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY DECLARATION DE CONFORMITÉ/ KONFORMITÄTSERKLÄRUNG

ABS127HPA	Hp-pette volume Plus completamente autoclavabili variabile 1000 - 5000 µl	<i>Hp-pette Plus Variable volumes Fully Autoclavables 1000 - 5000 µl</i>
ABS133HPA	Pipette MULTICANALE Plus completamente autoclavabili a volume variabile 5-50 µl	<i>Hp-pette Plus Variable volumes MULTICHANNEL Fully Autoclavables 5-50 µl</i>
ABS135HPA	Pipette MULTICANALE Plus completamente autoclavabili a volume variabile 50-300 µl	<i>Hp-pette Plus Variable volumes MULTICHANNEL Fully Autoclavables 50-300 µl</i>
ABS136HPA	Pipette MULTICANALE Plus completamente autoclavabili a volume variabile 0,5-10 µl	<i>Hp-pette Plus Variable volumes MULTICHANNEL Fully Autoclavables 0,5-10 µl</i>
REF	PIASTRE PER MICROMETODICHE	MICROTITER PLATE
BSD001	PIASTRA PER MICROMETODI, 96 POZZETTI FONDO A "U"	<i>MICRO-METHOD PLATES, 96 WELLS "U" BOTTOM</i>
BSD010	PIASTRA PER MICROMETODI, 96 POZZETTI FONDO A "U" STERILE	<i>MICRO-METHOD PLATE, 96 WELLS "U" BOTTOM, STERILE</i>
BSD002	PIASTRA PER MICROMETODI, 96 POZZETTI FONDO A "V"	<i>MICRO-METHOD PLATE, 96 WELLS "V" BOTTOM</i>
BSD011	PIASTRA PER MICROMETODI, 96 POZZETTI FONDO A "V"	<i>MICRO-METHOD PLATE, 96 WELLS "V" BOTTOM, STERILE</i>
BSD003	PIASTRA PER MICROMETODI, 96 POZZETTI FONDO PIATTO	<i>MICRO-METHOD PLATE, 96 WELLS FLAT BOTTOM</i>
BSD012	PIASTRA PER MICROMETODI, 96 POZZETTI FONDO PIATTO STERILE	<i>MICRO-METHOD PLATE, 96 WELLS FLAT BOTTOM, STERILE</i>
BSD013	COPERCHIO PER PIASTRA MICROMETODI	<i>CAPS FOR MICRO-METHOD PLATES</i>
REF	PIASTRE PETRI	PETRI PLATES
BSM200	PIASTRE PETRI DIAM. 60MM, STERILI, SENZA VENTILAZIONE	<i>PETRI DISHES DIAM.60MM STERILE WITHOUT VENTS</i>
BSM212	PIASTRE PETRI DIAM. 90MM, CON VENTILAZIONE	<i>PETRI DISHES DIAM.90 MM WITH VENTS</i>
BSM2120	PIASTRE PETRI DIAM. 90MM, STERILI, CON VENTILAZIONE	<i>PETRI DISHES DIAM.90 MM STERILE WITH VENTS</i>
BSM230	PIASTRE PETRI DIAM. 150MM, STERILI, CON VENTILAZIONE	<i>PETRI DISHES DIAM.150 MM STERILE WITH VENTS</i>
BSM239	PIASTRE PETRI QUADRATE 120x120MM, STERILI, CON VENTILAZIONE	<i>SQUARE PETRI DISHES 120X120MM STERILE, WITH VENTS</i>
BSM235	PIASTRE PETRI A 2 SETTORI DIAM.90MM, STERILI, CON VENTILAZIONE	<i>PETRI DISHES 2 SECTORS,DIAM.90MM STERILE, WITH VENTS</i>
BSM236	PIASTRE PETRI A 3 SETTORI DIAM.90MM, STERILI, CON VENTILAZIONE	<i>PETRI DISHES 3 SECTORS, DIAM.90MM STERILE, WITH VENTS</i>
BSM237	PIASTRE PETRI A 4 SETTORI DIAM.90MM, STERILI, CON VENTILAZIONE	<i>PETRI DISHES 4 SECTORS,DIAM.90MM STERILE, WITH VENTS</i>
BSM190	PIASTRE PETRI TIPO CONTACT DIAM.55MM, STERILI CON VENTILAZIONE	<i>PETRI DISHES CONTACT TYPE DIAM. 55MM STERILE, WITH VENTS</i>
BSM191	PIASTRE PETRI TIPO CONTACT DIAM.90MM, STERILI CON VENTILAZIONE	<i>PETRI DISHES CONTACT TYPE DIAM. 90MM STERILE, WITH VENTS</i>



Biosigma s.r.l.
a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY DECLARATION DE CONFORMITÉ/ KONFORMITÄTSERKLÄRUNG

REF	PIPETTE	PIPETTES
BSV134	Pipetta Pasteur 1 ml in PE	<i>Pasteur pipettes 1 ml in PE</i>
390509	Pipetta Pasteur 1 ml in PE – conf. 500 pezzi	<i>Pasteur pipettes 1 ml in PE – conf. 500 pcs</i>
BSV135	Pipetta Pasteur 1 ml in PE conf. 5 pz. STERILE	<i>Pasteur pipettes 1ml in PE bags 5 pcs. STERILE</i>
BSV1351	Pipetta Pasteur 1 ml in PE conf. singola STERILE	<i>Pasteur pipettes 1ml in PE individually wrapped STERILE</i>
BSV136	Pipetta Pasteur 3 ml in PE	<i>Pasteur pipettes 3 ml in PE</i>
390512	Pipetta Pasteur 3 ml in PE – conf. 500 pezzi	<i>Pasteur pipettes 3 ml in PE – conf. 500 pcs</i>
BSV137	Pipetta Pasteur 3 ml in PE conf. 5 pz. STERILE	<i>Pasteur pipettes 3ml in PE bags 5 pcs. STERILE</i>
BSV1371	Pipetta Pasteur 3 ml in PE conf. singola STERILE	<i>Pasteur pipettes 3ml in PE individually wrapped STERILE</i>
BSV1344	Pipetta Pasteur 1 ml in PE	<i>Pasteur pipettes 1 ml in PE</i>
BSV13514	Pipetta Pasteur 1 ml in PE conf. singola STERILE	<i>Pasteur pipettes 1ml in PE individually wrapped STERILE</i>
BSV1364	Pipetta Pasteur 3 ml in PE	<i>Pasteur pipettes 3 ml in PE</i>
BSV1364/20	Pipetta Pasteur 3 ml in PE in sacchetti da 20 pezzi	<i>Pasteur pipettes 3 ml in PE in bags of 20 pcs.</i>
BSV140	Pipetta Pasteur per Fast Read in PE	<i>Fast Read Pasteur pipettes in PE</i>
BSV143	Pipetta Pasteur Triplo Bulbo	<i>Triple bulbs Pasteur Pipette</i>
4249	Pipetta Pasteur Triplo Bulbo	<i>Triple bulbs Pasteur Pipette</i>
BSV150	Pipetta Pasteur per test in card in PE	<i>Pasteur pipettes for card-tests in PE</i>
BSV180	Pipetta Pasteur 10 ml in PE	<i>Pasteur pipettes 10 ml in PE</i>
REF	PROVETTE E TAPPI	TUBES AND STOPPERS
4175948001	Provetta conica 10ml - in PS- tappo giallo + etichetta+ rack	<i>Conical tube 10ml in PS Yellow stopper labelled in rack</i>
BSA059	Provetta cilindrica senza bordo 5ml, 12x75mm in PS	<i>Cylindrical tube rimless 5ml, 12x75mm in PS</i>
BSA059V	Provetta cilindrica senza bordo 5ml, 12x75mm in PS	<i>Cylindrical tube rimless 5ml, 12x75mm in PS</i>
BSA059/BLU	Provetta cilindrica senza bordo 5ml, 12x75mm in PP BLU	<i>Cylindrical tube rimless 5ml, 12x75mm in PP BLUE</i>
BSA059/G	Provetta cilindrica senza bordo 5ml, 12x75mm in PP GIALLA	<i>Cylindrical tube rimless 5ml, 12x75ml in PP YELLOW</i>
BSA059V/P	Provetta cilindrica senza bordo 5ml, 12x75mm in PP neutra	<i>Cylindrical tube rimless 5ml, 12x75mm in PP neutral</i>
BSA059/RA	Provetta cilindrica senza bordo 5ml, 12x75mm in PP ROSA	<i>Cylindrical tube rimless 5ml, 12x75mm in PP PINK</i>
BSA059/V	Provetta cilindrica senza bordo 5ml, 12x75mm in PP VERDE	<i>Cylindrical tube rimless 5ml, 12x75mm in PP GREEN</i>
BSA0592	Provetta cilindrica senza bordo 5ml, 12x75mm in PS BLU	<i>Cylindrical tube rimless 5ml, 12x75mm in PS BLUE</i>
BSA0594	Provetta cilindrica senza bordo 6 ml, 13x75mm in PS	<i>Cylindrical tube rimless 6ml, 13x75mm in PS</i>
BSA0594V	Provetta cilindrica senza bordo 6 ml, 13x75mm in PS	<i>Cylindrical tube rimless 6ml, 13x75mm in PS</i>
BSA0594/P	Provetta cilindrica senza bordo 6 ml, 13x75mm in PP	<i>Cylindrical tube rimless 6ml, 13x75mm in PP</i>
BSM398	Provetta in PS 7ml 13x100mm	<i>Tube in PS 7ml 13x100mm</i>
BSM402	Provetta 16x100 mm. in PS etichettata Tappo Rosso Trattate a Raggi Beta	<i>Biotube 16x100 mm. In PS label RED Stopper Treated β rays</i>



Biosigma s.r.l.

a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY DECLARATION DE CONFORMITÉ/ KONFORMITÄTSEKTLÄRUNG

BSM410	Provetta 16x103 mm. in PS etichettata Tappo Rosso Trattate a Raggi Beta	<i>Biotube 16x103 mm. In PS label RED Stopper Treated β rays</i>
BSM440	Provetta fondo piatto vol. 10 ml, in PS, tappo a vite	<i>Flat bottom test-tube, vol. 10 ml in PS screw- cap</i>
BSM441	Provetta fondo piatto vol. 5 ml, in PP, tappo a vite	<i>Flat bottom test-tube, vol. 5 ml in PP screw-cap</i>
BSM474	Provetta conica in PP con tappo a vite, 50ml, graduata, Sterile	<i>Conical tube in PP, screw cap, 50 ml, graduated, Sterile</i>
BSM475	Provetta conica in PP con tappo a vite, 50ml, graduata, Sterile – In Rack da 25 pz.	<i>Conical tube in PP, screw cap, 50 ml, graduated, Sterile – In Racks da 25 pcs.</i>
BSM477	Provetta conica con base d'appoggio, in PP con tappo a vite, 50ml, graduata, Sterile	<i>Conical tube self-standing in PP, screw cap, 50 ml, graduated, Sterile</i>
BSM4811	Provetta conica graduata 15 ml in PS con tappo a vite separato	<i>Graduated Conical Tube 15ml in PS, separate screw cap</i>
BSM482	Provetta conica in PP con tappo a vite, 15ml, graduata, Sterile	<i>Conical tube in PP, screw cap, 15 ml, graduated, Sterile</i>
BSM483	Provetta conica in PS con tappo a vite, 15ml, graduata, Sterile	<i>Conical tube in PS, screw cap, 15 ml, graduated, Sterile</i>
BSM490	Provetta conica in PP con tappo a vite, 15ml, graduata, Sterile – In Rack da 50 pz	<i>Conical tube in PP, screw cap, 15 ml, graduated, Sterile – In Racks of 50 pcs.</i>
BSM491	Provetta conica in PS con tappo a vite, 15ml, graduata, Sterile - In rack da 50 pz.	<i>Conical tube in PS, screw cap, 15 ml, graduated, Sterile – In Racks of 50 pcs.</i>
BSP048	Provetta 10 ml 16x100 mm in PP	<i>Tube 10 ml 16x100 mm in PP</i>
BSP049	Provetta 16x100 mm in PP etichettata	<i>Tube 16x100 in PP labelled</i>
BSP054	Provetta 10 ml 16x100 mm in PMMA	<i>Tube 10 ml 16x100 mm in PMMA</i>
BSP055	Provetta 16x100 mm in PMMA etichettata	<i>Tube 16x100 in PMMA labelled</i>
BSP058	Provetta 10 ml 16x100 mm graduata in PP	<i>Tube 10 ml 16x100 mm graduated in PP</i>
BSP0584	Provetta 10 ml 16x100 mm graduata in PP tappata etichettata	<i>Tube 10 ml 16x100 mm graduated in PP capped and labelled</i>
BSP059	Provetta 10 ml graduata in PP etichettata	<i>Tube 10 ml in PP graduated labelled</i>
BSP060	Provetta conica 16x103mm graduata in PP	<i>Conical tube 16x103mm graduated in PP</i>
BSP060/100	Provetta conica 16x103mm graduata in PP	<i>Conical tube 16x103mm graduated in PP</i>
BSP060/200	Provetta conica 16x103mm graduata in PP	<i>Conical tube 16x103mm graduated in PP</i>
BSP0604	Provetta conica 10 ml graduata in PP tappo BIANCO etichettata	<i>Conical tube 10ml graduated labelled in PP WHITE stopper</i>
BSP061	Provetta conica 10 ml graduata in PP etichettata	<i>Conical Tube 10ml in PP graduated labelled</i>
BSP072	Provetta cilindrica 10 ml graduata in PS	<i>Cylindrical tube 10 ml graduated in PS</i>
BSP0721	Provetta 10ml 10x100mm in PS tappate sfuse	<i>Tube 10ml 10x100mm in PS capped in bulk</i>
BSP073	Provetta 10 ml graduata in PS etichettata	<i>Tube 10ml in PS graduated labelled</i>
BSP074	Provetta conica 10 ml graduata in PS	<i>Conical tube 10ml graduated in PS</i>
BSP074/100	Provetta conica 10 ml graduata in PS	<i>Conical tube 10ml graduated in PS</i>
BSP074/200	Provetta conica 10 ml graduata in PS	<i>Conical tube 10ml graduated in PS</i>
BSP0745	Provetta conica 16X103 mm 10 ml graduata in PS	<i>Conical tube 16x103 mm 10ml graduated in PS</i>
BSP0741	Provetta conica 10 ml in PS tappo a parte	<i>Conical Tube 10ml in PS stopper aside</i>



Biosigma s.r.l.

a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY DECLARATION DE CONFORMITÉ/ KONFORMITÄTSERKLÄRUNG

BSP075	Provetta conica 10 ml graduata in PS etichettata	<i>Conical Tube 10ml in PS graduated labelled</i>
BSP0751	Provetta conica 10 ml in PS etichettata tappo a parte	<i>Conical Tube 10ml in PS labelled stopper aside</i>
BSP076	Provetta 12x86mm 5ml cilindrica in polistirolo	<i>Test Tube 12x86mm 5ml cylindrical in polystyrene</i>
BSP086	Provetta 7 ml 16x65 mm in PP	<i>Test Tube 7 ml 16x65 mm in PP</i>
BSP111	Provetta conica 10 ml graduata in PS tappo BIANCO	<i>Conical tube 10ml graduated in PS WHITE stopper</i>
BSP111/BLU	Provetta conica 10 ml graduata in PS tappo BLU	<i>Conical tube 10ml graduated in PS BLUE stopper</i>
BSP111/G	Provetta conica 10 ml in PS tappo GIALLO	<i>Conical tube 10ml graduated in PS YELLOW stopper</i>
BSP111/RO	Provetta conica 10 ml in PS tappo ROSSO	<i>Conical tube 10ml graduated in PS RED stopper</i>
BSP112	Provetta conica 10 ml in PS tappo BIANCO etichettata	<i>Conical tube 10ml labelled in PS WHITE stopper</i>
BSP112/AR	Provetta conica 10 ml in PS tappo ARANCIO etichettata	<i>Conical tube 10ml labelled in PS ORANGE stopper</i>
BSP112/BLU	Provetta conica 10 ml in PS tappo BLU etichettata	<i>Conical tube 10ml labelled in PS BLUE stopper</i>
BSP112/G	Provetta conica 10 ml in PS tappo GIALLO etichettata	<i>Conical tube 10ml labelled in PS YELLOW stopper</i>
BSP112/RO	Provetta conica 10 ml in PS tappo ROSSO etichettata	<i>Conical tube 10ml labelled in PS RED stopper</i>
BSP112/V	Provetta conica 10 ml in PS tappo VERDE etichettata	<i>Conical tube 10ml labelled in PS GREEN stopper</i>
BSP113	Provetta conica 10 ml in PS tappo BIANCO etichettata in rack	<i>Conical tube 10ml labelled in PS WHITE stopper in rack</i>
BSP113/BLU	Provetta conica 10 ml in PS tappo BLU etichettata in rack	<i>Conical tube 10ml labelled in PS BLUE stopper in rack</i>
BSP113/G	Provetta conica 10 ml in PS tappo GIALLO etichettata in rack	<i>Conical tube 10ml labelled in PS YELLOW stopper in rack</i>
BSP113/NE	Provetta conica 10 ml in PS tappo NERO etichettata in rack	<i>Conical tube 10ml labelled in PS BLACK stopper in rack</i>
BSP113/RO	Provetta conica 10ml in PS tappo ROSSO +etichetta in rack	<i>Conical tube 10ml labelled in PS RED stopper in rack</i>
BSP1131	Provetta conica 10 ml in PS tappo NERO in rack	<i>Conical tube 10ml labelled in PS BLACK stopper</i>
BSP116/G	Provetta cilindrica 10ml in PS tappo GIALLO etichettata in rack	<i>Cylindrical tube 10ml labelled in PS YELLOW stopper in rack</i>
BSP116/VS	Provetta cilindrica 10ml in PS tappo VERDE SCURO etichettata in rack	<i>Cylindrical tube 10ml labelled in PS DARK GREEN stopper in rack</i>
BSP128	Provetta tipo pediatrico in PS	<i>Pediatric tube in PS</i>
BSS555	Provetta 13x75mm con K3 EDTA x 3ml di sangue tappo verde in gomma perforabile AUTOMED®	<i>AUTOMED® Test tube 13x75mm with K3EDTA x 3ml of blood pierceable violet rubber stopper</i>
BSS555/VI	Provetta 13x75mm con K3 EDTA x 3ml di sangue tappo viola in gomma perforabile AUTOMED®	<i>AUTOMED® Test tube 13x75mm with K3EDTA x 3ml of blood pierceable violet rubber stopper</i>
BST106	Provetta 12x86mm con 0,25 ml di sodio citrato tappo rosa in rack da 50 pezzi	<i>Tube 12x86mm with 0.25 ml sodium citrate cap pink in rack of 50 pieces</i>



Biosigma s.r.l.
a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY DECLARATION DE CONFORMITÉ/ KONFORMITÄTSERKLÄRUNG

BSTA18	Provetta 10 ml in PP con granuli e acceleratore tappo BIANCO in rack	<i>Tube 10 ml in PP with granulars and coagulation accelerator WHITE stopper in rack</i>
BSTA18/AZ	Provetta 10 ml in PP con granuli e acceleratore tappo AZZURRO in rack	<i>Tube 10 ml in PP with granulars and coagulation accelerator LIGHT BLUE stopper in rack</i>
BSTA18/BLU	Provetta 10 ml in PP con granuli e acceleratore tappo BLU in rack	<i>Tube 10 ml in PP with granulars and coagulation accelerator BLUE stopper in rack</i>
BSTA18/MAR	Provetta 10 ml in PP con granuli e acceleratore tappo MARRONE in rack	<i>Tube 10 ml in PP with granulars and coagulation accelerator BROWN stopper in rack</i>
BSTA18/NE	Provetta 10 ml in PP con granuli e acceleratore tappo NERO in rack	<i>Tube 10 ml in PP with granulars and coagulation accelerator BLACK stopper in rack</i>
BSTA18/RA	Provetta 10 ml in PP con granuli e acceleratore tappo ROSA in rack	<i>Tube 10 ml in PP with granulars and coagulation accelerator PINK stopper in rack</i>
BSTA18/RO	Provetta 10 ml in PP con granuli e acceleratore tappo ROSSO in rack	<i>Tube 10 ml in PP with granulars and coagulation accelerator RED stopper in rack</i>
BSTA18/VC	Provetta 10 ml in PP con granuli e acceleratore tappo VERDE CHIARO in rack	<i>Tube 10 ml in PP with granulars and coagulation accelerator LIGHT GREEN stopper in rack</i>
BSTA18/VI	Provetta 10 ml in PP con granuli e acceleratore tappo VIOLA in rack	<i>Tube 10 ml in PP with granulars and coagulation accelerator VIOLET stopper in rack</i>
BSTA18/VS	Provetta 10 ml in PP con granuli e acceleratore tappo VERDE SCURO in rack	<i>Tube 10 ml in PP with granulars and coagulation accelerator DARK GREEN stopper in rack</i>
BSTA19	Provetta 10 ml in PP con granuli e acceleratore tappo BIANCO etichettata in rack	<i>Tube 10 ml in PP labelled with granulars and coagulation accelerator WHITE stopper in rack</i>
BSTA19/AZ	Provetta 10 ml in PP con granuli e acceleratore tappo AZZURRO etichettata in rack	<i>Tube 10 ml in PP labelled with granulars and coagulation accelerator LIGHT BLUE stopper in rack</i>
BSTA19/BLU	Provetta 10 ml in PP con granuli e acceleratore tappo BLU etichettata in rack	<i>Tube 10 ml in PP labelled with granulars and coagulation accelerator BLUE stopper in rack</i>
BSTA19/RO	Provetta 10 ml in PP con granuli e acceleratore tappo ROSSO etichettata in rack	<i>Tube 10 ml in PP labelled with granulars and coagulation accelerator RED stopper in rack</i>
BSTA20/RO	Provetta in PP 12X75mm con granuli e acceleratore di coagulazione / tappo ROSSO	<i>Tube in PP 12x75mm with granulars and coagulation accelerator/ RED stopper</i>
BSTA21	Provetta 10 ml in PMMA con granuli e acceleratore tappo BIANCO in rack	<i>Tube 10 ml in PMMA with granulars and coagulation accelerator WHITE stopper in rack</i>
BSTA21/AR	Provetta 10 ml in PMMA con granuli e acceleratore tappo ARANCIO in rack	<i>Tube 10 ml in PMMA with granulars and coagulation accelerator ORANGE stopper in rack</i>
BSTA21/AZ	Provetta 10 ml in PMMA con granuli e acceleratore tappo AZZURRO in rack	<i>Tube 10 ml in PMMA with granulars and coagulation accelerator LIGHT BLUE stopper in rack</i>
BSTA21/BLU	Provetta 10 ml in PMMA con granuli e acceleratore tappo BLU in rack	<i>Tube 10 ml in PMMA with granulars and coagulation accelerator BLUE stopper in rack</i>
BSTA21/G	Provetta 10 ml in PMMA con granuli e acceleratore tappo GIALLO in rack	<i>Tube 10 ml in PMMA with granulars and coagulation accelerator YELLOW stopper in rack</i>
BSTA21/MAR	Provetta 10 ml in PMMA con granuli e acceleratore tappo MARRONE in rack	<i>Tube 10 ml in PMMA with granulars and coagulation accelerator BROWN stopper in rack</i>



Biosigma s.r.l.

a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY DECLARATION DE CONFORMITÉ/ KONFORMITÄTSERKLÄRUNG

		<i>rack</i>
BSTA21/NE	Provetta 10 ml in PMMA con granuli e acceleratore tappo NERO in rack	<i>Tube 10 ml in PMMA with granulars and coagulation accelerator BLACK stopper in rack</i>
BSTA21/RA	Provetta 10 ml in PMMA con granuli e acceleratore tappo ROSA in rack	<i>Tube 10 ml in PMMA with granulars and coagulation accelerator PINK stopper in rack</i>
BSTA21/RO	Provetta 10 ml in PMMA con granuli e acceleratore tappo ROSSO in rack	<i>Tube 10 ml in PMMA with granulars and coagulation accelerator RED stopper in rack</i>
BSTA21/VC	Provetta 10 ml in PMMA con granuli e acceleratore tappo VERDE CHIARO in rack	<i>Tube 10 ml in PMMA with granulars and coagulation accelerator LIGHT GREEN stopper in rack</i>
BSTA21/VI	Provetta 10 ml in PMMA con granuli e acceleratore tappo VIOLA in rack	<i>Tube 10 ml in PMMA with granulars and coagulation accelerator VIOLET stopper in rack</i>
BSTA21/VS	Provetta 10 ml in PMMA con granuli e acceleratore tappo VERDE SCURO in rack	<i>Tube 10 ml in PMMA with granulars and coagulation accelerator DARK GREEN stopper in rack</i>
BSTA22	Provetta 10 ml in PMMA con granuli e acceleratore tappo BIANCO etichettata in rack	<i>Tube 10 ml in PMMA labelled with granulars and coagulation accelerator WHITE stopper in rack</i>
BSTA22/AZ	Provetta 10 ml in PMMA con granuli e acceleratore tappo AZZURRO etichettata in rack	<i>Tube 10 ml in PMMA labelled with granulars and coagulation accelerator LIGHT BLUE stopper in rack</i>
BSTA22/BLU	Provetta 10 ml in PMMA con granuli e acceleratore tappo BLU etichettata in rack	<i>Tube 10 ml in PMMA labelled with granulars and coagulation accelerator BLUE stopper in rack</i>
BSTA22/G	Provetta 10 ml in PMMA con granuli e acceleratore tappo GIALLO etichettata in rack	<i>Tube 10 ml in PMMA labelled with granulars and coagulation accelerator YELLOW stopper in rack</i>
BSTA22/MAR	Provetta 10 ml in PMMA con granuli e acceleratore tappo MARRONE etichettata in rack	<i>Tube 10 ml in PMMA labelled with granulars and coagulation accelerator BROWN stopper in rack</i>
BSTA22/RO	Provetta 10 ml in PMMA con granuli e acceleratore tappo ROSSO etichettata in rack	<i>Tube 10 ml in PMMA labelled with granulars and coagulation accelerator RED stopper in rack</i>
BSTA22/VI	Provetta 10 ml in PMMA con granuli e acceleratore tappo VIOLA etichettata in rack	<i>Tube 10 ml in PMMA labelled with granulars and coagulation accelerator VIOLET stopper in rack</i>
BSTA23	Provetta 12x86mm in metacrilato per sierare con granuli separatori ed acceleratore tappo bianco	<i>Test tube 12x86mm in PMMA with separating granules + coagulation accelerator white stopper</i>
BSV006	SEDI-TUBE provetta in vetro con Sodio Citrato 4NC	<i>Glass SEDI-TUBE with Sodium Citrate 4NC</i>
U9005	Conical tubes 10 ml with cap	<i>Conical tubes 10 ml with cap</i>
U9030	Conical tubes 6.5 ml	<i>Conical tubes 6.5 ml</i>
4204	URIN TEK- provette per urina	<i>URIN TEK – urine tubes</i>
4204/U	Provetta conica 10 ml in PS tappo BLU etichettata in rack	<i>Conical tube 10ml labelled in PS BLUE stopper in rack</i>
080053	Tube fond rond 5ml PS 13x75mm 500/sachet	<i>Round base tube 5ml PS 13x75mm 500/bag</i>
080071	Tube fond rond 5ml PP 13x75mm 500/sachet	<i>Round base tube 5ml PP 13x75mm 500/bag</i>
BSP0607	Provetta conica in PP da 10 ml etichettata con tappo alettato	<i>Conical tube 10 ml labelled in PP, stopper with wings</i>



Biosigma s.r.l.
a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY DECLARATION DE CONFORMITÉ/ KONFORMITÄTSERKLÄRUNG

CL482	Provetta 15 ml conica, tappo a vite in conf. da 25 pz - Sterile - Clearline	Conical tube 15 ml, screw cap, in bag of 25 pcs. - Sterile - Clearline
CL490	Provetta 15 ml conica, tappo a vite in rack da 50 pz - Sterile - Clearline	Conical tube 15 ml, screw cap, in rack of 50 pcs. - Sterile - Clearline
CL477	Provetta 50 ml con base, tappo a vite in conf. da 25 pz - Sterile - Clearline	Self-standing tube 50 ml, screw cap in bag of 25 pcs - Sterile - Clearline
CL474	Provetta 50 ml conica, tappo a vite in conf. da 25 pz - Sterile - Clearline	Conical tube 50 ml, screw cap, in bag of 25 pcs. - Sterile - Clearline
CL475	Provetta 50 ml conica, tappo a vite in rack da 25 pz - Sterile - Clearline	Conical tube 50 ml, screw cap, in rack of 25 pcs. - Sterile - Clearline
REF	TAPPI	STOPPERS
BSA065	TAPPO alettato diametro 13mm in PE	Stopper with wings Ø 13mm in PE
BSO025	TAPPO diametro 16 mm in PE ROSSO	Stopper Ø 16mm in PE RED
BSO026	TAPPO diametro 16 mm in PE BIANCO	Stopper Ø 16mm in PE WHITE
BSO027	TAPPO diametro 16 mm in PE GIALLO	Stopper Ø 16mm in PE YELLOW
BSO027/100	TAPPO diametro 16 mm in PE GIALLO	Stopper Ø 16mm in PE YELLOW
BSO027/200	TAPPO diametro 16 mm in PE GIALLO	Stopper Ø 16mm in PE YELLOW
BSO028	TAPPO diametro 16 mm in PE VERDE SCURO	Stopper Ø 16mm in PE DARK GREEN
BSO029	TAPPO diametro 16 mm in PE AZZURRO	Stopper Ø 16mm in PE LIGHT BLUE
BSO030	TAPPO diametro 16 mm in PE BLU	Stopper Ø 16mm in PE BLUE
BSO031	TAPPO diametro 16 mm in PE VERDE CHIARO	Stopper Ø 16mm in PE LIGHT BLUE
BSO032	TAPPO diametro 16 mm in PE ROSA	Stopper Ø 16mm in PE PINK
BSO033	TAPPO diametro 16 mm in PE VIOLA	Stopper Ø 16mm in PE VIOLET
BSO034	TAPPO diametro 16 mm in PE ARANCIO	Stopper Ø 16mm in PE ORANGE
BSO035	TAPPO diametro 16 mm in PE NERO	Stopper Ø 16mm in PE BLACK
BSO049	TAPPO diametro 16 mm in PE MARRONE	Stopper Ø 16mm in PE BROWN
BSO100	TAPPO con tenuta di sicurezza per provette sorwall in PE bianco	Stopper with security holder for sorwall tubes in PE white
BSO100/BLU	TAPPO con tenuta di sicurezza per provette sorwall in PE blu	Stopper with security holder for sorwall tubes in PE blue
BSO100/GR	TAPPO con tenuta di sicurezza per provette sorwall in PE grigio	Stopper with security holder for sorwall tubes in PE Grey
BSO100RO	TAPPO con tenuta di sicurezza per provette sorwall in PE rosso	Stopper with security holder for sorwall tubes in PE red
BSO100/V	TAPPO con tenuta di sicurezza per provette sorwall in PE verde	Stopper with security holder for sorwall tubes in PE green
BSO113	Tappo Vacustopper Ø 13mm per richiusura provette vuoto predeterminato	Vacustopper universal in PE Ø13mm for predetermined vacuum tubes re-closing
BSP277	Scatola di stoccaggio per provette sottovuoto Ø 13x100 mm	Storage box for Ø 13x100 mm vacuum tubes
BSP278	Scatola di stoccaggio per provette sottovuoto Ø 13x75 mm	Storage box for Ø 13x75 mm vacuum tubes
REF	MICROPROVETTE E ACCESSORI	MICROTUBES AND ACCESSORIES
BSA021	Conetto tipo Olly 1,5 ml in PP NEUTRO	Microtube Olly type 1,5ml in PP neutral
BSA022	Conetto tipo Eppendorf 1,5 ml PP NEUTRO	Microtube Eppendorf type 1,5ml in PP neutral



Biosigma s.r.l.

a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY DECLARATION DE CONFORMITÉ/ KONFORMITÄTSEKTLÄRUNG

BSA022V	Conetto tipo Eppendorf 1,5 ml PP NEUTRO	<i>Microtube Eppendorf type 1,5ml in PP neutral</i>
BSA022/XX	Conetto tipo Eppendorf 1,5 ml PP colorato	<i>Microtube Eppendorf type 1,5ml in PP coloured</i>
BSA024	Conetto tipo Cobas PE NEUTRO	<i>Microtube Cobas type in PE neutral</i>
BSA024V	Conetto tipo Cobas PE NEUTRO	<i>Microtube Cobas type in PE neutral</i>
BSA024/XX	Conetto tipo Cobas PE colorato	<i>Microtube Cobas type in PE coloured</i>
BSA0270SL	Microtubo per centrifuga 5 ml in PP con tappo di sicurezza	<i>Micro-centrifuge tube 5 ml in PP with safety lock cap</i>
BSA0269	Conetto 2 ml con tappo in PP	<i>Neutral microtube 2 ml with cap in PP</i>
366380B	/	<i>Screw cap neutral with o-ring</i>
366381B	/	<i>Screw cap with o-ring WHITE</i>
366383B	/	<i>Screw cap with o-ring RED</i>
366384B	/	<i>Screw cap with o-ring LIGHT BLUE</i>
366385B	/	<i>Screw cap with o-ring GREEN</i>
366386B	/	<i>Screw cap with o-ring YELLOW</i>
716201B	/	<i>Microtube 1,5 ml in PP conical, engraved graduation</i>
716261B	/	<i>Microtube 1,5 ml conical in PP, engraved graduation, with white screw cap + o-ring, sterile</i>
716261B/500	/	<i>Microtube 1,5 ml conical in PP, engraved graduation, with white screw cap + o-ring, sterile (14 bags of 500pcs)</i>
717201B	/	<i>Microtube 1,5 ml in PP self-standing, engraved graduated</i>
717261B	/	<i>Microtube 1,5 ml self-standing in PP, engraved graduation, with white screw cap + o-ring, sterile</i>
717261B/500	/	<i>Microtube 1,5 ml self-standing in PP, engraved graduation, with white screw cap + o-ring, sterile (14 bags of 500pcs)</i>
722201B	/	<i>Microtube 2 ml in PP self-standing, engraved graduation</i>
722261B	/	<i>Microtube 2 ml self-standing in PP, engraved graduation, with screw cap + o-ring, sterile</i>
722261B/500	/	<i>Microtube 2 ml self-standing in PP, engraved graduation, with screw cap + o-ring, sterile(14 bags of 500pcs)</i>
BSM910/XX	CAP DISK colorato per microprovette	<i>Coloured CAP DISK for microtubes</i>
CL022	Conetto con tappo tipo Eppendorf 1,5ml in PP - CLEARLINE	<i>Microtube with cap Eppendorf type 1.5ml in PP - CLEARLINE</i>
CL022B	Conetto con tappo tipo Eppendorf 1,5ml in PP- CLEARLINE – in bulk	<i>Microtube with cap Eppendorf type 1.5ml in PP - CLEARLINE – in bulk</i>
CL022BB	Conetto con tappo tipo Eppendorf 1,5ml in PP- CLEARLINE – in bulk (20 sacchetti da 500 pz)	<i>Microtube with cap Eppendorf type 1.5ml in PP - CLEARLINE – in bulk (20 bags of 500pcs)</i>
CL0269	Conetto 2ml con tappo in PP - CLEARLINE	<i>Microtube 2ml with cap in PP – CLEARLINE</i>
CL0269B	Conetto 2ml con tappo in PP - CLEARLINE – in bulk	<i>Microtube 2ml with cap in PP - CLEARLINE – in bulk</i>



Biosigma s.r.l.
a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY DECLARATION DE CONFORMITÉ/ KONFORMITÄTSERKLÄRUNG

CL0269BB	Conetto 2ml con tappo in PP - CLEARLINE – in bulk (16 sacchetti da 500 pz)	<i>Microtube 2ml with cap in PP - CLEARLINE – in bulk (16 bags of 500pcs)</i>
CL600	Microprovetta da 0,5 ml conica, CLEARLINE	<i>Microtube 0,5 ml conical, CLEARLINE</i>
CL601	Microprovetta da 0,5 ml conica, tappo a vite, CLEARLINE	<i>Microtube 0,5 ml conical, screw cap, CLEARLINE</i>
CL601/S	Microprovetta da 0,5 ml conica, tappo a vite, sterile, CLEARLINE	<i>Microtube 0,5 ml conical, screw cap, sterile, CLEARLINE</i>
CL602	Microprovetta da 0,5 ml conica, tappo a vite per CAP DISK, CLEARLINE	<i>Microtube 0,5 ml conical, screw cap suitable suitable for CAP DISK , CLEARLINE</i>
CL602/S	Microprovetta da 0,5 ml conica, tappo a vite per CAP DISK, sterile, CLEARLINE	<i>Microtube 0,5 ml conical, screw cap suitable suitable for CAP DISK , sterile, CLEARLINE</i>
CL603	Microprovetta da 0,5 ml conica, tappo a vite per CAP DISK con collare, CLEARLINE	<i>Microtube 0,5 ml conical, screw cap suitable suitable for CAP DISK with collar, CLEARLINE</i>
CL603/S	Microprovetta da 0,5 ml conica, tappo a vite per CAP DISK con collare, sterile, CLEARLINE	<i>Microtube 0,5 ml conical, screw cap suitable suitable for CAP DISK with collar, sterile, CLEARLINE</i>
CL604	Microprovetta da 0,5 ml conica, tappo a vite con o-ring, CLEARLINE	<i>Microtube 0,5 ml conical, screw cap with o-ring, CLEARLINE</i>
CL604/S	Microprovetta da 0,5 ml conica, tappo a vite con o-ring, sterile, CLEARLINE	<i>Microtube 0,5 ml conical, screw cap with o-ring, sterile, CLEARLINE</i>
CL605	Microprovetta da 0,5 ml conica, tappo a vite per CAP DISK con o-ring, CLEARLINE	<i>Microtube 0,5 ml conical, screw cap suitable suitable for CAP DISK with o-ring, CLEARLINE</i>
CL605/S	Microprovetta da 0,5 ml conica, tappo a vite per CAP DISK con o-ring, sterile, CLEARLINE	<i>Microtube 0,5 ml conical, screw cap suitable suitable for CAP DISK with o-ring, sterile, CLEARLINE</i>
CL606	Microprovetta da 0,5 ml conica, tappo a vite per CAP DISK con collare e o-ring, CLEARLINE	<i>Microtube 0,5 ml conical, screw cap suitable suitable for CAP DISK with collar and o-ring, CLEARLINE</i>
CL606/S	Microprovetta da 0,5 ml conica, tappo a vite per CAP DISK con collare e o-ring, sterile, CLEARLINE	<i>Microtube 0,5 ml conical, screw cap suitable suitable for CAP DISK with collar and o-ring, sterile, CLEARLINE</i>
CL610	Microprovetta da 0,5 ml con base, CLEARLINE	<i>Microtube 0,5 ml self-standing, CLEARLINE</i>
CL611	Microprovetta da 0,5 ml con base, tappo a vite, CLEARLINE	<i>Microtube 0,5 ml self-standing , screw cap, CLEARLINE</i>
CL611/S	Microprovetta da 0,5 ml con base, tappo a vite, sterile, CLEARLINE	<i>Microtube 0,5 ml self-standing, screw cap, sterile, CLEARLINE</i>
CL612	Microprovetta da 0,5 ml con base, tappo a vite per CAP DISK, CLEARLINE	<i>Microtube 0,5 ml self-standing, screw cap suitable suitable for CAP DISK , CLEARLINE</i>
CL612/S	Microprovetta da 0,5 ml con base, tappo a vite per CAP DISK, sterile, CLEARLINE	<i>Microtube 0,5 ml self-standing, screw cap suitable suitable for CAP DISK , sterile, CLEARLINE</i>
CL613	Microprovetta da 0,5 ml con base, tappo a vite per CAP DISK con collare, CLEARLINE	<i>Microtube 0,5 ml self-standing, screw cap suitable suitable for CAP DISK with collar, CLEARLINE</i>
CL613/S	Microprovetta da 0,5 ml con base, tappo a vite per CAP DISK con collare, sterile, CLEARLINE	<i>Microtube 0,5 ml self-standing, screw cap suitable suitable for CAP DISK with collar, sterile, CLEARLINE</i>
CL614	Microprovetta da 0,5 ml con base, tappo a vite con o-ring, CLEARLINE	<i>Microtube 0,5 ml self-standing, screw cap with o-ring, CLEARLINE</i>



Biosigma s.r.l.

a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY DECLARATION DE CONFORMITÉ/ KONFORMITÄTSERKLÄRUNG

CL614/S	Microprovetta da 0,5 ml con base, tappo a vite con o-ring, sterile, CLEARLINE	<i>Microtube 0,5 ml self-standing, screw cap with o-ring, sterile, CLEARLINE</i>
CL615	Microprovetta da 0,5 ml con base, tappo a vite per CAP DISK con o-ring, CLEARLINE	<i>Microtube 0,5 ml self-standing, screw cap suitable suitable for CAP DISK with o-ring, CLEARLINE</i>
CL615/S	Microprovetta da 0,5 ml con base, tappo a vite per CAP DISK con o-ring, sterile, CLEARLINE	<i>Microtube 0,5 ml self-standing, screw cap suitable suitable for CAP DISK with o-ring, sterile, CLEARLINE</i>
CL616	Microprovetta da 0,5 ml con base, tappo a vite per CAP DISK con collare e o-ring, CLEARLINE	<i>Microtube 0,5 ml self-standing, screw cap suitable suitable for CAP DISK with collar and o-ring, CLEARLINE</i>
CL616/S	Microprovetta da 0,5 ml con base, tappo a vite per CAP DISK con collare e o-ring, sterile, CLEARLINE	<i>Microtube 0,5 ml self-standing, screw cap suitable suitable for CAP DISK with collar and o-ring, sterile, CLEARLINE</i>
CL720	Microprovetta da 1,5 ml conica, graduata, CLEARLINE	<i>Microtube 1,5 ml conical, graduated, CLEARLINE</i>
CL721	Microprovetta da 1,5 ml conica, tappo a vite, graduata, CLEARLINE	<i>Microtube 1,5 ml conical, screw cap, graduated, CLEARLINE</i>
CL721/S	Microprovetta da 1,5 ml conica, tappo a vite, graduata, sterile, CLEARLINE	<i>Microtube 1,5 ml conical, screw cap, graduated, sterile, CLEARLINE</i>
CL722	Microprovetta da 1,5 ml conica, tappo a vite per CAP DISK, graduata, CLEARLINE	<i>Microtube 1,5 ml conical, screw cap suitable suitable for CAP DISK, graduated, CLEARLINE</i>
CL722/S	Microprovetta da 1,5 ml conica, tappo a vite per CAP DISK, graduata, sterile, CLEARLINE	<i>Microtube 1,5 ml conical, screw cap suitable suitable for CAP DISK, graduated, sterile, CLEARLINE</i>
CL723	Microprovetta da 1,5 ml conica, tappo a vite per CAP DISK con collare, graduata, CLEARLINE	<i>Microtube 1,5 ml conical, screw cap suitable suitable for CAP DISK with collar, graduated, CLEARLINE</i>
CL723/S	Microprovetta da 1,5 ml conica, tappo a vite per CAP DISK con collare, graduata, sterile, CLEARLINE	<i>Microtube 1,5 ml conical, screw cap suitable suitable for CAP DISK with collar, graduated, sterile, CLEARLINE</i>
CL724	Microprovetta da 1,5 ml conica, tappo a vite con o-ring, graduata, CLEARLINE	<i>Microtube 1,5 ml conical, screw cap with o-ring, graduated, CLEARLINE</i>
CL724/S	Microprovetta da 1,5 ml conica, tappo a vite con o-ring, graduata, sterile, CLEARLINE	<i>Microtube 1,5 ml conical, screw cap with o-ring, graduated, sterile, CLEARLINE</i>
CL725	Microprovetta da 1,5 ml conica, tappo a vite per CAP DISK con o-ring, graduata, CLEARLINE	<i>Microtube 1,5 ml conical, screw cap suitable suitable for CAP DISK with o-ring, graduated, CLEARLINE</i>
CL725/S	Microprovetta da 1,5 ml conica, tappo a vite per CAP DISK con o-ring, graduata, sterile, CLEARLINE	<i>Microtube 1,5 ml conical, screw cap suitable suitable for CAP DISK with o-ring, graduated, sterile, CLEARLINE</i>
CL726	Microprovetta da 1,5 ml conica, tappo a vite per CAP DISK con collare e o-ring, graduata, CLEARLINE	<i>Microtube 1,5 ml conical, screw cap suitable suitable for CAP DISK with collar and o-ring, graduated, CLEARLINE</i>
CL726/S	Microprovetta da 1,5 ml conica, tappo a vite per CAP DISK con collare e o-ring, graduata, sterile, CLEARLINE	<i>Microtube 1,5 ml conical, screw cap suitable suitable for CAP DISK with collar and o-ring, graduated, sterile, CLEARLINE</i>
CL730	Microprovetta da 1,5 ml con base, graduata, CLEARLINE	<i>Microtube 1,5 ml self-standing, graduated, CLEARLINE</i>
CL731	Microprovetta da 1,5 ml con base, tappo a vite, graduata, CLEARLINE	<i>Microtube 1,5 ml self-standing, screw cap, graduated, CLEARLINE</i>



Biosigma s.r.l.
a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY DECLARATION DE CONFORMITÉ/ KONFORMITÄTSERKLÄRUNG

CL731/S	Microprovetta da 1,5 ml con base, tappo a vite, graduata, sterile, CLEARLINE	<i>Microtube 1,5 ml self-standing, screw cap, graduated, sterile, CLEARLINE</i>
CL732	Microprovetta da 1,5 ml con base, tappo a vite per CAP DISK, graduata, CLEARLINE	<i>Microtube 1,5 ml self-standing, screw cap suitable suitable for CAP DISK , graduated, CLEARLINE</i>
CL732/S	Microprovetta da 1,5 ml con base, tappo a vite per CAP DISK, graduata, sterile, CLEARLINE	<i>Microtube 1,5 ml self-standing, screw cap suitable suitable for CAP DISK , graduated, sterile, CLEARLINE</i>
CL733	Microprovetta da 1,5 ml con base, tappo a vite per CAP DISK con collare, graduata, CLEARLINE	<i>Microtube 1,5 ml self-standing, screw cap suitable suitable for CAP DISK with collar, graduated, CLEARLINE</i>
CL733/S	Microprovetta da 1,5 ml con base, tappo a vite per CAP DISK con collare, graduata, sterile, CLEARLINE	<i>Microtube 1,5 ml self-standing, screw cap suitable suitable for CAP DISK with collar, graduated, sterile, CLEARLINE</i>
CL734	Microprovetta da 1,5 ml con base, tappo a vite con o-ring, graduata, CLEARLINE	<i>Microtube 1,5 ml self-standing, screw cap with o-ring, graduated, CLEARLINE</i>
CL734/S	Microprovetta da 1,5 ml con base, tappo a vite con o-ring, graduata, sterile, CLEARLINE	<i>Microtube 1,5 ml self-standing, screw cap with o-ring, graduated, sterile, CLEARLINE</i>
CL735	Microprovetta da 1,5 ml con base, tappo a vite per CAP DISK con o-ring, graduata, CLEARLINE	<i>Microtube 1,5 ml self-standing, screw cap suitable suitable for CAP DISK with o-ring, graduated, CLEARLINE</i>
CL735/S	Microprovetta da 1,5 ml con base, tappo a vite per CAP DISK con o-ring, graduata, sterile, CLEARLINE	<i>Microtube 1,5 ml self-standing, screw cap suitable suitable for CAP DISK with o-ring, graduated, sterile, CLEARLINE</i>
CL736	Microprovetta da 1,5 ml con base, tappo a vite per CAP DISK con collare e o-ring, graduata, CLEARLINE	<i>Microtube 1,5 ml self-standing, screw cap suitable suitable for CAP DISK with collar and o-ring, graduated, CLEARLINE</i>
CL736/S	Microprovetta da 1,5 ml con base, tappo a vite per CAP DISK con collare e o-ring, graduata, sterile, CLEARLINE	<i>Microtube 1,5 ml self-standing, screw cap suitable suitable for CAP DISK with collar and o-ring, graduated, sterile, CLEARLINE</i>
CL740	Microprovetta da 2 ml conica, graduata, CLEARLINE	<i>Microtube 2 ml conical, graduated, CLEARLINE</i>
CL741	Microprovetta da 2 ml conica, tappo a vite graduata	<i>Microtube 2 ml conical, screw cap, graduated, CLEARLINE</i>
CL741/S	Microprovetta da 2 ml conica, tappo a vite graduata, sterile, CLEARLINE	<i>Microtube 2 ml conical, screw cap, graduated, sterile, CLEARLINE</i>
CL742	Microprovetta da 2 ml conica, tappo a vite per CAP DISK, graduata, CLEARLINE	<i>Microtube 2 ml conical, screw cap suitable suitable for CAP DISK , graduated, CLEARLINE</i>
CL742/S	Microprovetta da 2 ml conica, tappo a vite per CAP DISK, graduata, sterile, CLEARLINE	<i>Microtube 2 ml conical, screw cap suitable suitable for CAP DISK , graduated, sterile, CLEARLINE</i>
CL743	Microprovetta da 2 ml conica, tappo a vite per CAP DISK con collare, graduata, CLEARLINE	<i>Microtube 2 ml conical, screw cap suitable suitable for CAP DISK with collar, graduated, CLEARLINE</i>
CL743/S	Microprovetta da 2 ml conica, tappo a vite per CAP DISK con collare, graduata, sterile, CLEARLINE	<i>Microtube 2 ml conical, screw cap suitable suitable for CAP DISK with collar, graduated, sterile, CLEARLINE</i>
CL744	Microprovetta da 2 ml conica, tappo a vite con o-ring, graduata, CLEARLINE	<i>Microtube 2 ml conical, screw cap with o-ring, graduated, CLEARLINE</i>
CL744/S	Microprovetta da 2 ml conica, tappo a vite con o-ring, graduata, sterile, CLEARLINE	<i>Microtube 2 ml conical, screw cap with o-ring, graduated, sterile, CLEARLINE</i>



Biosigma s.r.l.
a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY DECLARATION DE CONFORMITÉ/ KONFORMITÄTSERKLÄRUNG

CL745	Microprovetta da 2 ml conica, tappo a vite per CAP DISK con o-ring, graduata, CLEARLINE	<i>Microtube 2 ml conical, screw cap suitable suitable for CAP DISK with o-ring, graduated, CLEARLINE</i>
CL745/S	Microprovetta da 2 ml conica, tappo a vite per CAP DISK con o-ring, graduata, sterile, CLEARLINE	<i>Microtube 2 ml conical, screw cap suitable suitable for CAP DISK with o-ring, graduated, sterile, CLEARLINE</i>
CL746	Microprovetta da 2 ml conica, tappo a vite per CAP DISK con collare e o-ring, graduata, CLEARLINE	<i>Microtube 2 ml conical, screw cap suitable suitable for CAP DISK with collar and o-ring, graduated, CLEARLINE</i>
CL746/S	Microprovetta da 2 ml conica, tappo a vite per CAP DISK con collare e o-ring, graduata, sterile, CLEARLINE	<i>Microtube 2 ml conical, screw cap suitable suitable for CAP DISK with collar and o-ring, graduated, sterile, CLEARLINE</i>
CL750	Microprovetta 2 ml con base, graduata, CLEARLINE	<i>Microtube 2 ml self-standing, graduated, CLEARLINE</i>
CL751	Microprovetta 2 ml con base, con tappo a vite, graduata, CLEARLINE	<i>Microtube 2 ml self-standing, screw cap, Graduated, CLEARLINE</i>
CL751/S	Microprovetta 2 ml con base, con tappo a vite, graduata, sterile, CLEARLINE	<i>Microtube 2 ml self-standing, screw cap, Graduated, sterile, CLEARLINE</i>
CL752	Microprovetta 2 ml con base, tappo a vite per CAP DISK, graduata, CLEARLINE	<i>Microtube 2 ml self-standing, screw cap suitable suitable for CAP DISK , graduated, CLEARLINE</i>
CL752/S	Microprovetta 2 ml con base, tappo a vite per CAP DISK, graduata, sterile, CLEARLINE	<i>Microtube 2 ml self-standing, screw cap suitable suitable for CAP DISK , graduated, sterile, CLEARLINE</i>
CL753	Microprovetta 2 ml con base, tappo a vite per CAP DISK con collare, graduata, CLEARLINE	<i>Microtube 2 ml self-standing, screw cap suitable suitable for CAP DISK with collar, graduated, CLEARLINE</i>
CL753/S	Microprovetta 2 ml con base, tappo a vite per CAP DISK con collare, graduata, sterile, CLEARLINE	<i>Microtube 2 ml self-standing, screw cap suitable suitable for CAP DISK with collar, graduated, sterile, CLEARLINE</i>
CL754	Microprovetta 2 ml con base, tappo a vite con o-ring, graduata, CLEARLINE	<i>Microtube 2 ml self-standing, screw cap with o-ring, Graduated, CLEARLINE</i>
CL754/S	Microprovetta 2 ml con base, tappo a vite con o-ring, graduata, sterile, CLEARLINE	<i>Microtube 2 ml self-standing, screw cap with o-ring, Graduated, sterile, CLEARLINE</i>
CL755	Microprovetta 2 ml con base, tappo a vite per CAP DISK con o-ring, graduata, CLEARLINE	<i>Microtube 2 ml self-standing, screw cap suitable suitable for CAP DISK with o-ring, graduated, CLEARLINE</i>
CL755/S	Microprovetta 2 ml con base, tappo a vite per CAP DISK con o-ring, graduata, sterile, CLEARLINE	<i>Microtube 2 ml self-standing, screw cap suitable suitable for CAP DISK with o-ring, graduated, sterile, CLEARLINE</i>
CL756	Microprovetta 2 ml con base, tappo a vite per CAP DISK con collare e o-ring, graduata, CLEARLINE	<i>Microtube 2 ml self-standing, screw cap suitable suitable for CAP DISK with collar and o-ring, graduated, CLEARLINE</i>
CL756/S	Microprovetta 2 ml con base, tappo a vite per CAP DISK con collare e o-ring, graduata, sterile, CLEARLINE	<i>Microtube 2 ml self-standing, screw cap suitable suitable for CAP DISK with collar and o-ring, graduated, sterile, CLEARLINE</i>
CL800	Tappo a vite per microprovette, CLEARLINE	<i>Screw cap for microtubes, CLEARLINE</i>
CL800/XX	Tappo a vite per microprovette, CLEARLINE/colorato	<i>Coloured screw cap for microtubes, CLEARLINE</i>
CL801	Tappo a vite per CAP DISK per microprovette, CLEARLINE	<i>Screw cap suitable suitable for CAP DISK for microtubes, CLEARLINE</i>



Biosigma s.r.l.
a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



**DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY
DECLARATION DE CONFORMITE' / KONFORMITÄTSERKLÄRUNG**

CL801/XX	Tappo a vite per CAP DISK per microprovette, CLEARLINE/colorato	<i>Coloured Screw cap suitable for CAP DISK for microtubes, CLEARLINE</i>
CL802	Tappo a vite per CAP DISK con collare per microprovette, CLEARLINE	<i>Screw cap suitable suitable for CAP DISK with collar for microtubes, CLEARLINE</i>
CL802/XX	Tappo a vite per CAP DISK con collare per microprovette, CLEARLINE/colorato	<i>Coloured Screw cap suitable for CAP DISK with collar for microtubes, CLEARLINE</i>
CL810	Tappo a vite con o-ring per microprovette, CLEARLINE	<i>Screw cap with o-ring for microtubes, CLEARLINE</i>
CL810/XX	Tappo a vite con o-ring per microprovette, CLEARLINE/colorato	<i>Coloured Screw cap with o-ring for microtubes, CLEARLINE</i>
CL811	Tappo a vite con o-ring per CAP DISK per microprovette, CLEARLINE	<i>Screw cap with o-ring suitable for CAP DISK for microtubes, CLEARLINE</i>
CL811/XX	Tappo a vite con o-ring per CAP DISK per microprovette, CLEARLINE/colorato	<i>Coloured Screw cap with o-ring suitable for CAP DISK for microtubes, CLEARLINE</i>
CL812	Tappo a vite con o-ring per CAP DISK con collare per microprovette, CLEARLINE	<i>Screw cap with o-ring suitable for CAP DISK with collar for microtubes, CLEARLINE</i>
CL812/XX	Tappo a vite per CAP DISK con collare e o-ring per microprovette, CLEARLINE/colorato	<i>Coloured Screw cap with o-ring suitable for CAP DISK with collar for microtubes, CLEARLINE</i>
BSM811	Tappo a vite per CAP DISK con o-ring per microprovette	<i>Screw cap suitable for insert with o-ring for microtubes</i>
CL05.001.0500	MICROTUBI CLEAR-LOCK 0,5 ml CLEARLINE	<i>MICROTUBES CLEAR-LOCK 0,5 ml CLEARLINE</i>
CL15.002.0500	MICROTUBI CLEAR-LOCK 1,5 ml CLEARLINE	<i>MICROTUBES CLEAR-LOCK 1,5 ml CLEARLINE</i>
CL20.002.0500	MICROTUBI CLEAR-LOCK 2 ml CLEARLINE	<i>MICROTUBES CLEAR-LOCK 2 ml CLEARLINE</i>
CL50.002.0100	MICROTUBI CLEAR-LOCK 5 ml CLEARLINE	<i>MICROTUBES CLEAR-LOCK 5 ml CLEARLINE</i>
CL05.010.1000B	MICROTUBI CLEAR-LOCK 0,5 ml CLEARLINE - IN BULK	<i>MICROTUBES CLEAR-LOCK 0,5 ml CLEARLINE - IN BULK</i>
CL05.020.0500B	MICROTUBI CLEAR-LOCK 0,5 ml CLEARLINE - IN BULK	<i>MICROTUBES CLEAR-LOCK 0,5 ml CLEARLINE - IN BULK</i>
CL15.010.1000B	MICROTUBI CLEAR-LOCK 1,5 ml CLEARLINE - IN BULK	<i>MICROTUBES CLEAR-LOCK 1,5 ml CLEARLINE - IN BULK</i>
CL15.020.0500B	MICROTUBI CLEAR-LOCK 1,5 ml CLEARLINE - IN BULK	<i>MICROTUBES CLEAR-LOCK 1,5 ml CLEARLINE - IN BULK</i>
CL20.008.1000B	MICROTUBI CLEAR-LOCK 2 ml CLEARLINE - IN BULK	<i>MICROTUBES CLEAR-LOCK 2 ml CLEARLINE - IN BULK</i>
CL20.016.0500B	MICROTUBI CLEAR-LOCK 2 ml CLEARLINE - IN BULK	<i>MICROTUBES CLEAR-LOCK 2 ml CLEARLINE - IN BULK</i>
CL50.006.0500B	MICROTUBI CLEAR-LOCK 5 ml CLEARLINE - IN BULK	<i>MICROTUBES CLEAR-LOCK 5 ml CLEARLINE - IN BULK</i>
CL50.030.0100B	MICROTUBI CLEAR-LOCK 5 ml CLEARLINE - IN BULK	<i>MICROTUBES CLEAR-LOCK 5 ml CLEARLINE - IN BULK</i>
CL05.050.0001S	MICROTUBI CLEAR-LOCK 0,5 ml CLEARLINE STERILI	<i>MICROTUBES CLEAR-LOCK 0,5 ml CLEARLINE STERILI</i>
CL05.005.0100S	MICROTUBI CLEAR-LOCK 0,5 ml CLEARLINE STERILI	<i>MICROTUBES CLEAR-LOCK 0,5 ml CLEARLINE STERILI</i>



Biosigma s.r.l.
a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



**DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY
DECLARATION DE CONFORMITÉ/ KONFORMITÄTSEKTLÄRUNG**

CL15.100.0001S	MICROTUBI CLEAR-LOCK 1,5 ml CLEARLINE STERILI	<i>MICROTUBES CLEAR-LOCK 1,5 ml CLEARLINE STERILI</i>
CL15.004.0050S	MICROTUBI CLEAR-LOCK 1,5 ml CLEARLINE STERILI	<i>MICROTUBES CLEAR-LOCK 1,5 ml CLEARLINE STERILI</i>
CL20.100.0001S	MICROTUBI CLEAR-LOCK 2 ml CLEARLINE STERILI	<i>MICROTUBES CLEAR-LOCK 2 ml CLEARLINE STERILI</i>
CL20.004.0050S	MICROTUBI CLEAR-LOCK 2 ml CLEARLINE STERILI	<i>MICROTUBES CLEAR-LOCK 2 ml CLEARLINE STERILI</i>
CL50.050.0001S	MICROTUBI CLEAR-LOCK 5 ml CLEARLINE STERILI	<i>MICROTUBES CLEAR-LOCK 5 ml CLEARLINE STERILI</i>
CL50.010.0020S	MICROTUBI CLEAR-LOCK 5 ml CLEARLINE STERILI	<i>MICROTUBES CLEAR-LOCK 5 ml CLEARLINE STERILI</i>
BSA05.010.1000	MICROTUBI 0,5 ML CON CHIUSURA DI SICUREZZA	<i>MICROTUBES 0.5 ML WITH LOCKING SYSTEM</i>
BSA05.020.0500	MICROTUBI 0,5 ML CON CHIUSURA DI SICUREZZA	<i>MICROTUBES 0.5 ML WITH LOCKING SYSTEM</i>
BSA15.010.1000	MICROTUBI 1,5 ML CON CHIUSURA DI SICUREZZA	<i>MICROTUBES 1.5 ML WITH LOCKING SYSTEM</i>
BSA15.020.0500	MICROTUBI 1,5 ML CON CHIUSURA DI SICUREZZA	<i>MICROTUBES 1.5 ML WITH LOCKING SYSTEM</i>
BSA20.008.1000	MICROTUBI 2 ML CON CHIUSURA DI SICUREZZA	<i>MICROTUBES 2 ML WITH LOCKING SYSTEM</i>
BSA20.016.0500	MICROTUBI 2 ML CON CHIUSURA DI SICUREZZA	<i>MICROTUBES 2 ML WITH LOCKING SYSTEM</i>
BSA50.006.0500	MICROTUBI 5 ML CON CHIUSURA DI SICUREZZA	<i>MICROTUBES 5 ML WITH LOCKING SYSTEM</i>
BSA50.030.100	MICROTUBI 5 ML CON CHIUSURA DI SICUREZZA	<i>MICROTUBES 5 ML WITH LOCKING SYSTEM</i>
REF	PUNTALI	TIPS
BSR001	PUNTALE GIALLO TIPO EPPENDORF- BRAND-SOCOREX DA 0 A 200 µl	<i>YELLOW TIP EPPENDORF-BRAND- SOCOREX TYPE FROM 0 TO 200 µl</i>
BSR002	PUNTALE AZZURRO TIPO EPPENDORF- BRAND-SOCOREX DA 201 A 1000 µl	<i>LIGHT BLUE TIP EPPENDORF-BRAND- SOCOREX TYPE FROM 201 TO 1000 µl</i>
BSR066	PUNTALE GIALLO TIPO GILSON DA 0 A 200 ml	<i>YELLOW TIP GILSON TYPE FROM 0 TO 20 0µl</i>
REF	STACCACOAGULO	BREAKCLOT
BSA066	Bacchetta staccacoagulo in PS	<i>Breakclot in PS</i>
REF	STERNHEIMER & MALBIN	STERNHEIMER & MALBIN
BVS135	FAST-READ STAIN 35ml	<i>FAST-READ STAIN 35ml</i>
4252	COLORANTE PER SEDIMENTO URINARIO	<i>SEDIMENT COLOURING</i>
REF	TAZZINE PORTACAMPIONI	SAMPLE CUPS
BSA027	Tazzina 1,5 ml in PS per automazione tipo Technicon	<i>Sample cup 1,5 ml in PS for automation Technicon type</i>
BSA030	Tazzina 0,5 ml in PS per automazione tipo Genseac	<i>Sample cup 0.5 ml in PS for automation Genseac type</i>



Biosigma s.r.l.
a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



**DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY
DECLARATION DE CONFORMITÉ/ KONFORMITÄTSEKTLÄRUNG**

BSA031S	Tazzina tipo Hitachi 3ml in PS	Sample cup Hitachi type 3 ml in PS
BSA031V	Tazzina tipo Hitachi 3ml in PS	Sample cup Hitachi type 3 ml in PS
BSA031VS	Tazzina tipo Hitachi 3ml in PS	Sample cup Hitachi type 3 ml in PS
BSA046	Tazzina tipo Axsym in PS	Sample cup Axsym type in PS
BSA102	Coppetta per coagulazione tipo Coag-Mat	Coagulation Cup Coag-Mat type
BSA5258	Tazzina portacampioni 10x 30 mm in PS bordo ottagonale	Sample cup 10x30 mm in PS octagonal edge
BSA636S	Riduttore 1ml in PS per provette Ø 13 mm	Sample Nest cup 1 ml in PS for tube Ø13/16mm
BSA636V	Riduttore 1ml in PS per provette Ø 13 mm	Sample Nest cup 1 ml in PS for tube Ø13/16mm
BSA636VS	Riduttore 1ml in PS per provette Ø 13 mm	Sample Nest cup 1 ml in PS for tube Ø13/16mm
BSA646S	Riduttore 2ml in PS per provette Ø 16 mm	Sample Nest cup 2 ml in PS for tube Ø16 mm
BSA651	Tazzina portacampione utilizzabile su Olympus AU400/AU640	Sample Cup utilizable on Olympus AU400/AU640
BSA658	Tazzina CEX tipo Express 550 e Plus	CEX sample cup Express 550 and Plus type
BSA659	Tazzina CEX pediatrica	CEX pediatric cup Express 550 and Plus type
BSA69	Coppetta portacampione in PP tipo Biotechnica Apparatus	Sample Cup in PP Biotechnica Apparatus type
BSA70	Piastra di reazione 30 fori 900 S tipo Plus	Reaction Plate 30 holes 900 S type
BSA73	Coppetta per coagulazione in PS tipo Optin, Benk, Turbotimer	Coagulation Cup in PS Option, Benk, Turbotimer type
BSA731	Coppetta per coagulazione in PS tipo Optin, Benk, Turbotimer con sfere metalliche	Coagulation Cup in PS Option, Benk, Turbotimer type with metallic ball mixer
BSA74	Coppetta per coagulazione in PS con mixer metallico tipo Accu Coa Data	Coagulation Cup in PS with metallic steel stick, Accu Coa Data type
BSA77	Tazzina per coagulazione tipo Seac Clot S (2)	Coagulation cup Seac Clot S (2) type
BSA86	Coppetta pediatrica per Hitachi tipo Elecsys, Cobas Integra, Core	Pediatric Sample Cup for Hitachi, Elecsys, Cobas Integra and Core type
R6044	Sample cups 1,5 ml	Sample cups 1,5 ml
U5505	Sample cups 0.5 ml	Sample cups 0.5 ml
4162315001	Coppetta tipo Hitachi 3ml in PS	Sample cup Hitachi type 3 ml in PS
39019	Coppette portacampioni 3 ml	Sample cups 3 ml
BSA88	Coppetta portacampione in PS tipo A.A.	Sample cup A.A. type in PS
REF	PIPETTE PER V.E.S.	E.S.R PIPETTES
BSV012	Pipette per V.E.S.	E.S.R. Pipettes
REF	PORTAOGGETTO Dim. ~ 26 x 76 mm	SLIDES Dim. ~ 26 x 76 mm
NORMALI		
VBS653/A	TAGLIATI	CUT EDGES
VBS654/A	TAGLIATI con DOPPIA AREA SMERIGLIATA	CUT EDGES WITH DOUBLE FROSTED AREA
VBS655/A	MOLATI	GROUND EDGES
SILVER QUALITY		
VBS653	TAGLIATI – SILVER QUALITY	CUT EDGES - SILVER QUALITY
VBS653/50	TAGLIATI – SILVER QUALITY - conf. 50 pz.	CUT EDGES - SILVER QUALITY – conf. 50 pcs.



Biosigma s.r.l.

a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY DECLARATION DE CONFORMITÉ/ KONFORMITÄTSEKKLÄRUNG

VBS654	TAGLIATI con DOPPIA AREA SMERIGLIATA – SILVER QUALITY	CUT EDGES WITH DOUBLE FROSTED AREA - SILVER QUALITY
VBS654/50	TAGLIATI con DOPPIA AREA SMERIGLIATA – SILVER QUALITY- conf. 50 pz.	CUT EDGES WITH DOUBLE FROSTED AREA - SILVER QUALITY- conf. 50 pcs.
VBS655	MOLATI – SILVER QUALITY	GROUND EDGES - SILVER QUALITY
VBS655/50	MOLATI – SILVER QUALITY- conf. 50 pz.	GROUND EDGES - SILVER QUALITY- conf. 50 pcs.
VBS656	MOLATI con DOPPIA AREA SMERIGLIATA- SILVER QUALITY	GROUND EDGES WITH DOUBLE FROSTED AREA - SILVER QUALITY
VBS656/50	MOLATI con DOPPIA AREA SMERIGLIATA - SILVER QUALITY- conf. 50 pz.	GROUND EDGES WITH DOUBLE FROSTED AREA - SILVER QUALITY- conf. 50 pcs.
VBS656/AR	MOLATI con AREA COLORATA ARANCIONE SILVER QUALITY	GROUND EDGES WITH ORANGE COULORED AREA - SILVER QUALITY
VBS656/AR50	MOLATI con AREA COLORATA ARANCIONE - SILVER QUALITY- conf. 50 pz.	GROUND EDGES WITH ORANGE COULORED AREA - SILVER QUALITY- conf. 50 pcs.
VBS656/AZ	MOLATI con AREA COLORATA AZZURRA SILVER QUALITY	GROUND EDGES WITH LIGHT BLUE COULORED AREA - SILVER QUALITY
VBS656/AZ50	MOLATI con AREA COLORATA AZZURRA - SILVER QUALITY- conf. 50 pz.	GROUND EDGES WITH LIGHT BLUE COULORED AREA - SILVER QUALITY- conf. 50 pcs.
VBS656/B	MOLATI con AREA COLORATA BIANCA SILVER QUALITY	GROUND EDGES WITH WHITE COULORED AREA - SILVER QUALITY
VBS656/B50	MOLATI con AREA COLORATA BIANCA - SILVER QUALITY- conf. 50 pz	GROUND EDGES WITH WHITE COULORED AREA - SILVER QUALITY- conf. 50 pcs.
VBS656/G	MOLATI con AREA COLORATA GIALLA SILVER QUALITY	GROUND EDGES WITH YELLOW COULORED AREA - SILVER QUALITY
VBS656/G50	MOLATI con AREA COLORATA GIALLA - SILVER QUALITY- conf. 50 pz.	GROUND EDGES WITH YELLOW COULORED AREA - SILVER QUALITY- conf. 50 pcs.
VBS656/RA	MOLATI con AREA COLORATA ROSA SILVER QUALITY	GROUND EDGES WITH PINK COULORED AREA - SILVER QUALITY
VBS656/RA50	MOLATI con AREA COLORATA ROSA – SILVER QUALITY- conf. 50 pz.	GROUND EDGES WITH PINK COULORED AREA - SILVER QUALITY- conf. 50 pcs.
VBS656/V	MOLATI con AREA COLORATA VERDE SILVER QUALITY	GROUND EDGES WITH GREEN COULORED AREA - SILVER QUALITY
VBS656/V50	MOLATI con AREA COLORATA VERDE - SILVER QUALITY- conf. 50 pz.	GROUND EDGES WITH GREEN COULORED AREA - SILVER QUALITY- conf. 50 pcs.
REF	COPRIOGGETTO	COVER SLIPS
VBS633	18 x 18 mm	18 x 18 mm
VBS634	20 x 20 mm	20 x 20 mm
VBS636	22 x 22 mm	22 x 22 mm
VBS637	24 x 24 mm	24 x 24 mm
VBS638	24 x 32 mm	24 x 32 mm
VBS649	24 x 40 mm	24 x 40 mm
VBS650	24 x 50 mm	24 x 50 mm
VBS651	24 x 60 mm	24 x 60 mm
REF	VETRINI PER SEDIMENTO URINARIO	URINARY SEDIMENT SLIDES
BVS100	Fast-Read 102	Fast-Read 102
BVS100H	Fast-Read 102	Fast-Read 102
U9027	URIGLASS	URIGLASS



Biosigma s.r.l.
a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY DECLARATION DE CONFORMITÉ/ KONFORMITÄTSERKLÄRUNG

4163419001	Fast-Read 102	Fast-Read 102
390497	Fast-Read 102	Fast-Read 102
4242	Fast-Read 102	Fast-Read 102
REF	PIPETTE SIEROLOGICHE	SEROLOGICAL PIPETTES
CL100	Pipetta sierologica 1 ml in 1/100 in confezione singola - Sterile - CLEARLINE	Serological pipette 1 ml in 1/100 individually packed - Sterile - CLEARLINE
CL101	Pipetta sierologica 2 ml in 1/100 in confezione singola - Sterile - CLEARLINE	Serological pipette 2 ml in 1/100 individually packed - Sterile - CLEARLINE
CL102	Pipetta sierologica 5 ml in 1/10 in confezione singola - Sterile - CLEARLINE	Serological pipette 5 ml in 1/10 individually packed - Sterile - CLEARLINE
CL103	Pipetta sierologica 10 ml in 1/10 in confezione singola - Sterile - CLEARLINE	Serological pipette 10 ml in 1/10 individually packed - Sterile - CLEARLINE
CL111	Pipetta sierologica 25 ml in 2/10 in confezione singola - Sterile - CLEARLINE	Pipetta sierologica 25 ml in 2/10 individually packed - Sterile - CLEARLINE
CL112	Pipetta sierologica 50 ml in 1 in confezione singola - Sterile - CLEARLINE	Serological pipette 50 ml in 1 Individually packed - Sterile - CLEARLINE
CL113	Pipetta sierologica 100 ml in 1 in confezione singola - Sterile - CLEARLINE	Serological pipette 100 ml in 1 individually packed - Sterile - CLEARLINE
CL105	Pipetta sierologica 1 ml in 1/100 in sacchetti da 25 - Sterile - CLEARLINE	Serological pipette 1 ml in 1/100 in sacchetti da 25 - Sterile - CLEARLINE
CL106	Pipetta sierologica 2 ml in 1/100 in sacchetti da 25 - Sterile - CLEARLINE	Serological pipette 2 ml in 1/100 in sacchetti da 25 - Sterile - CLEARLINE
CL107	Pipetta sierologica 5 ml in 1/10 in sacchetti da 25 - Sterile - CLEARLINE	Serological pipette 5 ml in 1/10 in sacchetti da 25 - Sterile - CLEARLINE
CL108	Pipetta sierologica 10 ml in 1/10 in sacchetti da 25 - Sterile - CLEARLINE	Serological pipette 10 ml in 1/10 in sacchetti da 25 - Sterile - CLEARLINE
CL109	Pipetta sierologica 25 ml in 2/10 in sacchetti da 25 - Sterile - CLEARLINE	Serological pipette 25 ml in 2/10 in sacchetti da 25 - Sterile - CLEARLINE
REF	SACCHETTI BIOHAZARD	BIOHAZARD BAGS
BSA2401	Sacchetto a doppia tasca per trasporto campioni con simbolo biohazard - trasparente	Bag with double pocket for specimens transport with biohazard symbol - transparent
BSA2101	Sacchetto a doppia tasca per trasporto campioni con simbolo biohazard - giallo	Bag with double pocket for specimens transport with biohazard symbol - yellow
BSA2201	Sacchetto a doppia tasca per trasporto campioni con simbolo biohazard - verde	Bag with double pocket for specimens transport with biohazard symbol - green
BSA2301	Sacchetto a doppia tasca per trasporto campioni con simbolo biohazard - azzurro	Bag with double pocket for specimens transport with biohazard symbol - light blue
BSA2461	Sacchetto a doppia tasca per trasporto campioni con simbolo biohazard - rosso	Bag with double pocket for specimens transport with biohazard symbol - red
BSA2601	Sacchetto a doppia tasca per trasporto campioni con simbolo biohazard - arancione	Bag with double pocket for specimens transport with biohazard symbol - orange
BSA2701	Sacchetto a doppia tasca per trasporto campioni con simbolo biohazard - rosa	Bag with double pocket for specimens transport with biohazard symbol - pink
BSA2801	Sacchetto a doppia tasca per trasporto campioni con simbolo biohazard - viola	Bag with double pocket for specimens transport with biohazard symbol - violet



Biosigma s.r.l.
a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626



DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY DECLARATION DE CONFORMITÈ/ KONFORMITÄTSERKLÄRUNG

Il presente Elenco si riferisce ai dispositivi medico-diagnostici in vitro registrati da Biosigma S.r.l. presso il Ministero della Sanità alla data riportata nella Dichiarazione di Conformità.

This List refers to in vitro diagnostic medical devices registered by Biosigma S.r.l. to Ministry for Health at the date reported on the Declaration of Conformity.



Biosigma s.r.l.
a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626