


FICHE TECHNIQUE

Référence	390494
Désignation	Portoir pour cryotubes



Description :

- ✓ Polycarbonate workstation
- ✓ For 40 Clearline cryotubes
- ✓ Autoclavable at 121°C (steam) for 20 min
- ✓ Locking of Clearline tubes inside the workstation
- ✓ Non-slip rubber feet
- ✓ Dimensions w x d x h : 215 x 100 x 25 mm
- ✓ Hole diameter : 12.8 mm
- ✓ Color : blue

Unité de vente :

1 Portoir
