

**Microtube Clearline pp 2 ml rnase free**

- ClearLine® Microtubes RNase-DNase free
- In medical grade virgin polypropylene
- Non-sterile
- Natural PP, optimal transparency
- Centrifugeable up to 25 000 x g
- Autoclavable (open) at +121°C for 20 min
- Ideal for PCR and molecular biology applications
- Graduated, with writing surface on the cap
- Packaged in zip-resealable bags of 500 or 1000
- Certified free from human DNA, RNase-DNase, ATP, pyrogens and PCR inhibitors
- Volume : 2 ml
- $\varnothing$  x h (mm) : 10.8 x 42
- Pack details : 500 x2

**Units : 1000**