



DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY DECLARATION DE CONFORMITÈ/ KONFORMITÄTSEKKLÄRUNG

Nome e indirizzo della ditta
Name and address of the firm
Name und Adresse der Firma
Nom et adresse de l'entreprise

BIOSIGMA S.r.l.
Via Valletta, 6 – Zona Ind. Cantarana
30010 - Cona (VE) - ITALY
Tel. +39.0426.302224 - Fax +39.0426.302228
E-mail: info@biosigmaeu.com
<http://www.biosigma.com>

Dichiariamo sotto nostra responsabilità che il dispositivo medico-diagnostico in vitro:

We declare under our sole responsibility that the in vitro diagnostic medical device:

Wir erklären in alleiniger Verantwortung, dass das Medizinprodukt für die In-vitro-Diagnostik:

Nous déclarons sous notre propre responsabilité que le dispositif médical de diagnostique in vitro:

**Famiglia: PRODOTTI PER CRIOGENIA
CRYOGENICS' PRODUCTS
PRODUITS POUR CRYOGÉNIE**

Vedi elenco allegato / *See the enclosed list*
Produkt der Liste beiliegend / Voir liste ci-joint /

Codice CND W0503010201 *Provette in materiale plastico per analisi*
EDMA code: 51091001 *Other containers for sample of human origin*
GMDN code: 0 *No Equivalent GMDN Term*

Classificazione secondo la direttiva relativa ai dispositivi medico-diagnostici in vitro 98/79/CE.

Classification according to the directive on in vitro diagnostic medical devices 98/79/EC.

Klassifizierung nach der Richtlinie über In-vitro-Diagnostika 98/79/EC.

Classification selon la directive relative aux dispositifs médicaux de diagnostic in vitro 98/79/CE.

Dispositivo dell'elenco A, allegato II / Produkt der Liste A, Anhang II / Dispositif de la liste B, annexe II / Device of List A, Annex II

Dispositivo dell'elenco B, allegato II / Device of List B, Annex II / Produkt der Liste B, Anhang II / Dispositif de la liste B, annexe II

Dispositivo per testautodiagnostico non elencato nell'allegato II / Device for self-testing not listed in Annex II / Produkt zur Eigenanwendung, das nicht in Anhang II genannt ist / Dispositif destiné à l'autodiagnostic non listé dans l'annexe II

Altro dispositivo / Other device / Sonstiges Produkt / Autre dispositif

soddisfa tutte le disposizioni della direttiva relativa ai dispositivi medico-diagnostici in vitro 98/79/CE che lo riguardano.

meets all the provisions of the directive on in vitro diagnostic medical devices 98/79/EC which apply to it.

allen Anforderungen der Richtlinie über In-vitro-Diagnostika 98/79/EC entspricht, die anwendbar sind.

remplit toutes les exigences de la directive relative aux dispositifs médicaux de diagnostic in vitro 98/79/CE qui le concernent.

Procedimento di valutazione della conformità:

Conformity assessment procedure:

Konformitätsbewertungsverfahren:

Procédure d'évaluation de la conformité:

Allegato III (escluso punto 6)

Annex III (excluding point 6)

Anhang III (mit Ausnahme von Punkt 6)

Annexe III (à l'exclusion, point 6)

Organo incaricato della valutazione della conformità (se consultato):

Notified Body (if consulted):

Konformitätsbewertungsstelle (falls beigezogen):

Organe resp. de l'évaluat. de la conform.(si consult):

N.A.



Biosigma s.r.l.
Via Valletta, 6
30010, Cantarana di Cona (VE) Italy
Part. IVA 03328440270
C.F. 02173800281
Tel. 0426 302224
Fax 0426 302228

CONA 07/06/2017

MARCOLIN MARTINO – Consigliere Delegato

Luogo, data / Place, date / Ort, Datum / Lieu, date

Nome e funzione / Name and function / Name und Funktion / Nom et fonction



Biosigma s.r.l.

a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) | Fax ++39 0426 748030

| Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |

P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT0903000005843 | R.N.P.P.A. IT11080P00002626

**Allegato alla Dichiarazione di Conformità CE - Famiglia: PRODOTTI PER CRIOGENIA**
Attached to Declaration of Conformity CE – Family: CRYOGENICS' PRODUCTS

Ultimo aggiornamento 08/05/2017

10001TS	Tappo esterno MARIA STELLA	MARIA STELLA external cap
10011TS	Tappo interno MARIO STELLA	MARIO STELLA internal cap
1BETS	Provetta CryoGen da 1,2 ml con base, per tappo Stella esterno	CryoGen Tube 1.2 ml self standing, for external Stella cap
2BETS	Provetta CryoGen da 2 ml con base, per tappo Stella esterno	CryoGen Tube 2 ml self standing, for external Stella cap
2BI	Provetta CryoGen 2ml con base, filetto interno	CryoGen Tube 2ml self standing, internal thread
2ETS	Provetta CryoGen da 2 ml, per tappo Stella esterno	CryoGen Tube 2 ml, for external Stella cap
2I	Provetta CryoGen 2ml filetto interno	CryoGen Tube 2ml internal thread
3BETS	Provetta CryoGen da 3 ml con base, per tappo Stella esterno	CryoGen Tube 3 ml self standing, for external Stella cap
4BETS	Provetta CryoGen da 4 ml con base, per tappo Stella esterno	CryoGen Tube 4 ml self standing, for external Stella cap
4BI	Provetta CryoGen da 4 ml con base, filetto interno	CryoGen Tube 4ml self standing, internal thread
4I	Provetta CryoGen da 4ml, filetto interno	CryoGen Tube 4ml, internal thread
5BETS	Provetta CryoGen da 5 ml con base, per tappo Stella esterno	CryoGen Tube 5 ml self standing, for external Stella cap
5I	Provetta CryoGen da 5ml filetto interno	CryoGen Tube 5ml internal thread
BSM58025/X	CRYOGEN® BOX25 per provette da 1,2 e 2 ml, griglia vari colori	CRYOGEN® BOX25 for 1.2 and 2 ml tubes, coloured grid
BSM58025B/8X	CRYOGEN® BOX25 per provette da 1,2 e 2 ml, fondo bianco, griglia vari colori - conf.8 pz	CRYOGEN® BOX25 for 1.2 and 2 ml tubes, white bottom, coloured grid - conf.8 pcs
BSM58025B/X	CRYOGEN® BOX25 per provette da 1,2 e 2 ml, fondo bianco, griglia vari colori	CRYOGEN® BOX25 for 1.2 and 2 ml tubes, white bottom, coloured grid
BSM58025B1D/8X	CRYOGEN® BOX25 per provette da 1,2 e 2 ml, fondo bianco con barcode, griglia vari colori - conf.8 pz	CRYOGEN® BOX25 for 1.2 and 2 ml tubes, white bottom with barcode, coloured grid - conf.8 pcs
BSM58025B1D/X	CRYOGEN® BOX25 per provette da 1,2 e 2 ml, fondo bianco con barcode, griglia vari colori	CRYOGEN® BOX25 for 1.2 and 2 ml tubes, white bottom with barcode, coloured grid
BSM58061/X	CRYOGEN® BOX81 per provette da 1,2 e 2 ml, griglia vari colori	CRYOGEN® BOX81 for 1,2 and 2 ml tubes, coloured grid
BSM580612D/5X	CRYOGEN® BOX81-2D per provette da 2 ml 2D, griglia vari colori - conf.5 pz	CRYOGEN® BOX81-2D for 2 ml tubes 2D, coloured grid - conf.5 pcs
BSM580612D/X	CRYOGEN® BOX81-2D per provette da 2 ml 2D, griglia vari colori	CRYOGEN® BOX81-2D for 2 ml tubes 2D, coloured grid
BSM58061B/5X	CRYOGEN® BOX81 per provette da 1,2 e 2 ml, fondo bianco, griglia vari colori - conf.5 pz	CRYOGEN® BOX81 for 1,2 and 2 ml tubes, white bottom, coloured grid - conf.5 pcs
BSM58061B/X	CRYOGEN® BOX81 per provette da 1,2 e 2 ml, fondo bianco, griglia vari colori	CRYOGEN® BOX81 for 1,2 and 2 ml tubes, white bottom, coloured grid
BSM58061B1D/5X	CRYOGEN® BOX81 per provette da 1,2 e 2 ml, fondo bianco con barcode, griglia vari colori - conf.5 pz	CRYOGEN® BOX81 for 1,2 and 2 ml tubes, white bottom with barcode, coloured grid - conf.5 pcs
BSM58061B1D/X	CRYOGEN® BOX81 per provette da 1,2 e 2 ml, fondo bianco con barcode, griglia vari colori	CRYOGEN® BOX81 for 1,2 and 2 ml tubes, white bottom with barcode, coloured grid
BSM58071/X	CRYOGEN® BOX81 per provette da 3, 4 e 5 ml, griglia vari colori	CRYOGEN® BOX81 for 3, 4 and 5 ml tubes, coloured grid
BSM580712D/4X	CRYOGEN® BOX81-2D per provette da 5 ml 2D, griglia vari colori - conf.4 pz	CRYOGEN® BOX81-2D for 5 ml tubes 2D, coloured grid - conf.4 pcs
BSM580712D/X	CRYOGEN® BOX81-2D per provette da 5 ml 2D, griglia vari colori	CRYOGEN® BOX81-2D for 5 ml tubes 2D, coloured grid
BSM58071B/4X	CRYOGEN® BOX81 per provette da 3, 4 e 5 ml, fondo bianco, griglia vari colori - conf.4 pz	CRYOGEN® BOX81 for 3, 4 and 5 ml tubes, white bottom, coloured grid - conf.4 pcs
BSM58071B/X	CRYOGEN® BOX81 per provette da 3, 4 e 5 ml, fondo bianco, griglia vari colori	CRYOGEN® BOX81 for 3, 4 and 5 ml tubes, white bottom, coloured grid
BSM58071B1D/4X	CRYOGEN® BOX81 per provette da 3, 4 e 5 ml, fondo bianco con barcode, griglia vari colori - conf.4 pz	CRYOGEN® BOX81 for 3, 4 and 5 ml tubes, white bottom with barcode, coloured grid - conf.4 pcs
BSM58071B1D/X	CRYOGEN® BOX81 per provette da 3, 4 e 5 ml, fondo bianco con barcode, griglia vari colori	CRYOGEN® BOX81 for 3, 4 and 5 ml tubes, white bottom with barcode, coloured grid
BSM58081/X	CRYOGEN® BOX81 per provette da 1,2 e 2 ml, griglia vari colori, coperchio inciso	CRYOGEN® BOX81 for 1,2 and 2 ml tubes, coloured grid, engraved cover
BSM58081B/X	CRYOGEN® BOX81 per provette da 1,2 e 2 ml, griglia vari colori, coperchio inciso	CRYOGEN® BOX81 for 1,2 and 2 ml tubes, coloured grid, engraved cover
BSM58091/X	CRYOGEN® BOX81 per provette da 3, 4 e 5 ml, griglia vari colori, coperchio inciso	CRYOGEN® BOX81 for 3, 4 and 5 ml tubes, coloured grid, engraved cover
BSM58091B/X	CRYOGEN® BOX81 per provette da 3, 4 e 5 ml, griglia vari colori, coperchio inciso	CRYOGEN® BOX81 for 3, 4 and 5 ml tubes, coloured grid, engraved cover
BSM5810/X	WORKSTATION per provette CryoGen a 40 posti, vari colori	WORKSTATION for CryoGen tubes 40 holes, coloured

BSM58100/5X	CRYOGEN® BOX100 per provette da 1,2 e 2 ml, fondo vari colori - conf.5 pz	CRYOGEN® BOX100 for 1,2 and 2 ml tubes, coloured bottom conf.5 pcs
BSM58100/X	CRYOGEN® BOX100 per provette da 1,2 e 2 ml, fondo vari colori	CRYOGEN® BOX100 for 1,2 and 2 ml tubes, coloured bottom
BSM581001D/B	CRYOGEN® BOX100 per provette da 1,2 e 2 ml, fondo vari colori con barcode	CRYOGEN® BOX100 for 1,2 and 2 ml tubes, coloured bottom with barcode
BSM581001D/5B	CRYOGEN® BOX100 per provette da 1,2 e 2 ml, fondo vari colori con barcode- conf.5 pz	CRYOGEN® BOX100 for 1,2 and 2 ml tubes, coloured bottom with barcode- conf.5 pcs
BSM581002D/5X	CRYOGEN® BOX100-2D per provette da 2 ml 2D, fondo vari colori - conf.5 pz	CRYOGEN® BOX100-2D for 2 ml tubes 2D, coloured bottom - conf.5 pcs
BSM581002D/X	CRYOGEN® BOX100-2D per provette da 2 ml 2D, fondo vari colori	CRYOGEN® BOX100-2D for 2 ml tubes 2D, coloured bottom
BSM581052D/4X	CRYOGEN® BOX100-2D per provette da 5 ml 2D, fondo vari colori - conf.4 pz	CRYOGEN® BOX100-2D for 5 ml tubes 2D, coloured bottom - conf.4 pcs
BSM581052D/X	CRYOGEN® BOX100-2D per provette da 5 ml 2D, fondo vari colori	CRYOGEN® BOX100-2D for 5 ml tubes 2D, coloured bottom
CL1ARBEPSTS	Provetta CryoGen® da 1,2 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	CryoGen® Tube 1,2 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®
CL1ARBEPSTS/B	Provetta CryoGen® da 1,2 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	CryoGen® Tube 1,2 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK
CL1ARBEPSTS/BX	Provetta CryoGen® da 1,2 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	CryoGen® Tube 1,2 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK
CL1ARBEPSTS/X	Provetta CryoGen® da 1,2 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	CryoGen® Tube 1,2 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®
CL1ARBIPSTS	Provetta CryoGen® da 1,2 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	CryoGen® Tube 1,2 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®
CL1ARBIPSTS/B	Provetta CryoGen® da 1,2 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	CryoGen® Tube 1,2 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK
CL1ARBIPSTS/BX	Provetta CryoGen® da 1,2 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	CryoGen® Tube 1,2 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK
CL1ARBIPSTS/X	Provetta CryoGen® da 1,2 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	CryoGen® Tube 1,2 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®
CL2ARBEP2D	Provetta CryoGen® da 2 ml con base, tappo Stella esterno, graduata, bar-code + 2D, sterile - CLEARLINE®	CryoGen® Tube 2 ml self standing, external Stella cap, graduated, with bar-code + 2D, sterile - CLEARLINE®
CL2ARBEP2D/B	Provetta CryoGen® da 2 ml con base, tappo Stella esterno, graduata, bar-code + 2D, sterile - CLEARLINE® - IN BULK	CryoGen® Tube 2 ml self standing, external Stella cap, graduated, with bar-code + 2D, sterile - CLEARLINE® - IN BULK
CL2ARBEPSTS	Provetta CryoGen® da 2 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	CryoGen® Tube 2 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®
CL2ARBEPSTS/B	Provetta CryoGen® da 2 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	CryoGen® Tube 2 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK
CL2ARBEPSTS/BX	Provetta CryoGen® da 2 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	CryoGen® Tube 2 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK
CL2ARBEPSTS/X	Provetta CryoGen® da 2 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	CryoGen® Tube 2 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®
CL2ARBIPS2D	Provetta CryoGen® da 2 ml con base, tappo Stella interno, graduata, bar-code + 2D, sterile - CLEARLINE®	CryoGen® Tube 2 ml self standing, internal Stella cap, graduated, with bar-code + 2D, sterile - CLEARLINE®
CL2ARBIPS2D/B	Provetta CryoGen® da 2 ml con base, tappo Stella interno, graduata, bar-code + 2D, sterile - CLEARLINE® - IN BULK	CryoGen® Tube 2 ml self standing, internal Stella cap, graduated, with bar-code + 2D, sterile - CLEARLINE® - IN BULK
CL2ARBIPSTS	Provetta CryoGen® da 2 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	CryoGen® Tube 2 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®
CL2ARBIPSTS/B	Provetta CryoGen® da 2 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	CryoGen® Tube 2 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK
CL2ARBIPSTS/BX	Provetta CryoGen® da 2 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	CryoGen® Tube 2 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK
CL2ARBIPSTS/X	Provetta CryoGen® da 2 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	CryoGen® Tube 2 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®

CL2AREPSTS	Provetta CryoGen® da 2 ml , tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2AREPSTS/B	Provetta CryoGen® da 2 ml , tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2AREPSTS/BX	Provetta CryoGen® da 2 ml , tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2AREPSTS/X	Provetta CryoGen® da 2 ml , tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARIPSTS	Provetta CryoGen® da 2 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL2ARIPSTS/B	Provetta CryoGen® da 2 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2ARIPSTS/BX	Provetta CryoGen® da 2 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 2 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL2ARIPSTS/X	Provetta CryoGen® da 2 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 2 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL3ARBEPSSN	Provetta CryoGen® da 3 ml con base, tappo a vite esterno bi-iniettato, graduata, con area di scrittura, sterile - CLEARLINE® con bar-code a sequenza personalizzata	<i>CryoGen® Tube 3ml self standing, external bi-injected screw cap, graduated, with writing area, sterile - CLEARLINE® - with bar-code custom sequence</i>
CL3ARBEPSTS	Provetta CryoGen® da 3 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 3 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL3ARBEPSTS/B	Provetta CryoGen® da 3 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 3 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL3ARBEPSTS/BX	Provetta CryoGen® da 3 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 3 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL3ARBEPSTS/X	Provetta CryoGen® da 3 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 3 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARBEPSTS	Provetta CryoGen® da 4 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARBEPSTS/B	Provetta CryoGen® da 4 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARBEPSTS/BX	Provetta CryoGen® da 4 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARBEPSTS/X	Provetta CryoGen® da 4 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARBIPSTS	Provetta CryoGen® da 4 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARBIPSTS/B	Provetta CryoGen® da 4 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARBIPSTS/BX	Provetta CryoGen® da 4 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARBIPSTS/X	Provetta CryoGen® da 4 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARIPSTS	Provetta CryoGen® da 4 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL4ARIPSTS/B	Provetta CryoGen® da 4 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARIPSTS/BX	Provetta CryoGen® da 4 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 4 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL4ARIPSTS/X	Provetta CryoGen® da 4 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 4 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARBEP2D	Provetta CryoGen® da 5 ml con base, tappo Stella esterno, graduata, bar-code + 2D, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, external Stella cap, graduated, with bar-code + 2D, sterile - CLEARLINE®</i>
CL5ARBEP2D/B	Provetta CryoGen® da 5 ml con base, tappo Stella esterno, graduata, bar-code + 2D, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, external Stella cap, graduated, with bar-code + 2D, sterile - CLEARLINE® - IN BULK</i>

CL5ARBEPSTS	Provetta CryoGen® da 5 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARBEPSTS/B	Provetta CryoGen® da 5 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARBEPSTS/BX	Provetta CryoGen® da 5 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARBEPSTS/X	Provetta CryoGen® da 5 ml con base, tappo Stella esterno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, colored external Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARBIPS2D	Provetta CryoGen® da 5 ml con base, tappo Stella interno, graduata, bar-code + 2D, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, internal Stella cap, graduated, with bar-code + 2D, sterile - CLEARLINE®</i>
CL5ARBIPS2D/B	Provetta CryoGen® da 5 ml con base, tappo Stella interno, graduata, bar-code + 2D, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, internal Stella cap, graduated, with bar-code + 2D, sterile - CLEARLINE® - IN BULK</i>
CL5ARBIPSSN	Provetta CryoGen® da 5ml con base, tappo a vite interno bi-iniettato, graduata, con area di scrittura, sterile - CLEARLINE® con bar-code a sequenza personalizzata	<i>CryoGen® Tube 5ml self standing, internal bi-injected screw cap, graduated, with writing area, sterile - CLEARLINE® - with bar-code custom sequence</i>
CL5ARBIPSTS	Provetta CryoGen® da 5 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARBIPSTS/B	Provetta CryoGen® da 5 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARBIPSTS/BX	Provetta CryoGen® da 5 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARBIPSTS/X	Provetta CryoGen® da 5 ml con base, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml self standing, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARIPSTS	Provetta CryoGen® da 5 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CL5ARIPSTS/B	Provetta CryoGen® da 5 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARIPSTS/BX	Provetta CryoGen® da 5 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE® - IN BULK	<i>CryoGen® Tube 5 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE® - IN BULK</i>
CL5ARIPSTS/X	Provetta CryoGen® da 5 ml, tappo Stella interno colorato, graduata, bar-code, sterile - CLEARLINE®	<i>CryoGen® Tube 5 ml, colored internal Stella cap, graduated, with bar-code, sterile - CLEARLINE®</i>
CRY910TS/X	CAP DISK per tappo STELLA, vari colori	<i>CAPDISK for STELLA cap / colour</i>
BSM58025	GRIGLIA IN PC PER FONDO CRYOGEN BOX25	<i>GRID IN PC FOR BOTTOM CRYOGEN BOX25</i>
BSM58081	GRIGLIA IN PC PER FONDO CRYOGEN BOX81	<i>GRID IN PC FOR BOTTOM CRYOGEN BOX81</i>
BSM580100	GRIGLIA IN PC PER FONDO CRYOGEN BOX100	<i>GRID IN PC FOR BOTTOM CRYOGEN BOX100</i>
I2D	IgenomixXera Tube 2,0 ml 2D	<i>IgenomixXera Tube 2,0 ml 2D</i>