



# Transport Systems

ADR COMPLIANCE

THE TRANSPORT SYSTEMS PRODUCTS ALLOW FOR THE SECURE TRANSPORT OF BIOLOGICAL SUBSTANCE (UN 3373) CATEGORY B - FROM COLLECTION OF THE SPECIMEN RIGHT UP TO ANALYSIS IN THE LABORATORY.

## SECONDARY CONTAINERS VACUETTE® Transport Container (VTC) REF BS800105 - BS800110



REF	<b>BS800105</b> <b>VACUETTE® Transport Container (VTC) with foam insert for up to 12 tubes, with shipping box.</b>
	<b>BS800110</b> <b>VACUETTE® Transport Container (VTC) with foam insert for up to 12 tubes, without shipping box</b>
FEATURES	VACUETTE® Transport Containers ensure the safe transport of sample material from the collection point right up to analysis in the laboratory. The VACUETTE® Transport Containers can store twelve tubes and are available with or without shipping box. Dimension: 113 x 180 mm For transport of up to 12 tubes With screw cap and absorbent pad Foam insert suitable for tubes diameters from 9 to 16 mm
APPLICATION	secondary container for biological samples transport.
MATERIAL	PP polypropylene
DESCRIPTION	For tubes 13x75 – 13x100 – 16x100 with foam insert Not Sterile
MANUFACTURER	GREINER Bio-One GmbH
CE MARKING	Yes (Dir. 98/79/CE)