

Transport Systems

TERTIARY CONTAINERS

Bags for transport biological samples in accordance with ADR

- **External** in Polyester 600/900 PVC coated.
- **Inner** in Nylon 210/210 PVC coated. *
- **Semi-rigid structure** insulating with 8 mm padding in Elastolen. *
- **Zipper closure.**
- **Inside pocket** for housing the cooler element. *
- **External pocket** for documents. *
- **Symbols:** UN3373 - Biological Material Category "B"
- **Disinfectable** with aqueous solutions of the most common disinfectants
- **Can be sterilized** with aqueous solutions of cold sterilants
- **Washable** at 40 ° C
- **Blue color**

| | REF | Capacity | Dimensions cm | Internal pocket for |
|---|------------------|---|------------------|---------------------|
|  | BS720205E | 1 container BS800110 12 tubes | Ø12 x 15 | NA |
|  | BS720201E | 1 container BS472040 40 tubes | 28 x 18 x 21 | 1 cooler element |
|  | BS720202E | 2 containers BS472040 80 tubes | 32 x 30 x 22 | 2 cooler elements |
|  | BS720203E | 3 containers BS472040 120 tubes | 43 x 30 x 20 | 3 cooler elements |
|  | BS720204E | 4 containers BS472041 312 tubes | 35 x 37 x 39 | 4 cooler elements |
|  | BS720206E | 6 containers BS472040 240 tubes With shoulder strap | 42 x 27 x 39 | 4 cooler elements |