



Transport Systems

TERTIARY CONTAINERS

Bags for transport biological samples in accordance with ADR

- **External** in Polyester 600/900 PVC coated.
- **Inner** in Nylon 210/210 PVC coated. *
- **Semi-rigid structure** insulating with 8 mm padding in Elastolen. *
- **Zipper closure.**
- **Inside pocket** for housing the cooler element. *
- **External pocket** for documents. *
- **Symbols:** UN3373 - Biological Material Category "B"
- **Disinfectable** with aqueous solutions of the most common disinfectants
- **Can be sterilized** with aqueous solutions of cold sterilants
- **Washable** at 40 ° C
- **Blue color**

	REF	Capacity	Dimensions cm	Internal pocket for
	BS720205E	1 container BS800110 12 tubes	Ø12 x 15	NA
	BS720201E	1 container BS472040 40 tubes	28 x 18 x 21	1 cooler element
	BS720202E	2 containers BS472040 80 tubes	32 x 30 x 22	2 cooler elements
	BS720203E	3 containers BS472040 120 tubes	43 x 30 x 20	3 cooler elements
	BS720206E	6 containers BS472040 240 tubes With shoulder strap	42 x 27 x 39	4 cooler elements