

VACUETTE® Transport Box (VTB)

REF BS472040



REF	BS472040 VACUETTE® Transport Box (VTB) with foam insert for up to 40 tubes, without shipping box	
APPLICATION	Secondary container for biological samples transport. VACUETTE® Transport Box ensure the safe transport of sample material from the collection point right up to analysis in the laboratory. The VACUETTE® Transport Box can store 40 tubes	
DIMENSIONS:	270 mm × 125 mm x 180 mm	
VOLUME:	for up to 40 tubes	
RAW MATERIAL:	Transparent cover:	Semitransparent polypropylene
	Box:	Red polypropylene
	Seal:	Silicone
	Foam insert:	Polyurethane sponge for 40 tubes
	Dimensions foam insert:	220 mm x 95 mm x 70 mm
	Absorbent pad:	Polypropylene impregnated with absorbent polymer
	Dimensions absorbent pad:	190 mm x 85 mm x 3 mm
MANUFACTURER	GREINER Bio-One GmbH	
CE MARKING	Yes (Dir. 98/79/CE)	
STORAGE:	Temperatures from -18°C to +50°C Avoid direct sunlight	
RECOMMENDATIONS: DISINFECTION:	To clean, soak the VTB in a disinfectant solution recommended for medical instruments, according to the guidelines of your structure.	
AUTOCLAVABLE	No	
REFRIGERABLE	Yes - 20°C	