

TECHNICAL SPECIFICATIONS

Trade name	Microscope slide boxes
Use	Easy handling of 76x26 mm slides
Material	Polystyrene PS

DIMENSIONAL SPECIFICATIONS

Dimensions	mm. 96,70 x 229,60 x 34,00
General dimensional tolerance	mm. +/- 0,1
Capacity	50 places

GENERAL FEATURES

CE Mark	No
Autoclavability	No
Sterilization	To be effected by radiations
Foodstuff suitability	Yes

CHEMICAL RESISTANCES

Categories of substances	Reference	Concen. %	T. 20° C	T. 40° C	T.° 60 C
Inorganic Acid	Sulphuric acid	98	C	-	-
Organic Acid	Benzoic acid	100	B	B	B
Alcohol	Ethanol	100	B	C	-
Aldehydes	Acetaldeyde	100	C	-	-
Inorganic Alkali	Isodium hydroxide	50	A	A	A
Organic Alkali	Aniline	100	C	-	-
Ketone	Acetone	100	-	-	-
Esther	Ethyl acetate	100	C	-	-
Hydrocarbons, halogenated	Dichloroethylene	100	C	-	-
Hydrocarbons,aromatic	Benzene	100	C	-	-
Hydrocarbons, linear	Hexane	100	C	-	-
Supersolvent	Tetrahydrofuran	100	C	-	-

Legenda (not applicable if neglected):

- **A: fair resistance; exposure (30 days) to the chemical does not cause any damage.**
- **C: poor resistance; exposure to chemical is not allowed, or causes immediate permanent damages.**
- **B: sufficient resistance; exposure causes damage of poor importance, which sometimes is only temporary.**

Results of testing with reference materials are to be considered as indication: in case of specific use it is recommended to carry out preliminary testing.

PACKAGING

Bag	Pcs. 1
Box	Pcs. 40

Quality Assurance

Massimiliano Capitanio

Technical Direction

Guido Borona

DOMINIQUE DUTSHER SAS