

Noviglio, 01/01/2014

Reference: PRODUCTS "CE" MARK

**Description: IN VITRO DIAGNOSTIC MEDICAL DEVICES
VESSELS
DISPOSABLES
TIPS**

Products List: see attached list 1

We undersigned, Kartell S.P.A., Labware Division, declare under our exclusive responsibility that our above mentioned product belonging to our production range, is conforming to product specifications specified in 98/79/CE Directive.

Assicurazione Qualità

M. Capitanio

LIST 1

CODE	DESCRIPTION	CODE	DESCRIPTION
279-XX	Graduated micro test tubes w/cap	2819	Serological Pipette 10 ML Sterile
297	Micro test tubes type Eppendorf®	2820	Serological Pipette 10 ML Sterile
298-XX	Micro test tubes w/cap type Eppendorf®	2821	Serological Pipette 10 ML
300	Conical centrifuge tubes ø 16 10 ml	2822	Serological Pipette 25 ML Sterile
301	Conical centrifuge tubes, ø 18 15 ml	2823	Serological Pipette 25 ML Sterile
302	Centrifuge tubes, cylindrical ø 12 7 ml	2824	Serological Pipette 25 ML
303	Centrifuge tubes, cylindrical ø 18 16 ml	2825	Serological Pipette 50 ML Sterile
305	Centrifuge tubes, cylindrical ø 24 31 ml	2830	Bacteriological Pipette 1,1 ML Sterile
306	Centrifuge tubes, cylindrical ø 30 48 ml	2831	Bacteriological Pipette 1,1 ML Sterile
307	Centrifuge tubes, cylindrical ø 35 70 ml	2832	Bacteriological Pipette 2,2 ML Sterile
308	Centrifuge tubes, cylindrical ø 40 110 ml	2833	Bacteriological Pipette 2,2 ML Sterile
315	Disposable dropping pipettes 1,5 ml	2840	Open end pipette 1 ML Sterile
316	Disposable dropping pipettes 5 ml	2842	Open end pipette 2 ML Sterile
329	PASTEUR Pipette STER. Capillary	2844	Open end pipette 5 ML Sterile
330	PASTEUR Pipette STER. 1 ml	2846	Open end pipette 10 ML Sterile
334	PASTEUR Pipette STER. 3 ml	2899-03	Disposable syringes for dispenser / diluter
358	PETRI dishes with triple vents ø 90	2900-03	Disposable syringes for dispenser / diluter
359	PETRI dishes with triple vents ø 100	2901-03	Disposable syringes for dispenser / diluter
360	PETRI dishes without triple vents ø 60	2902-03	Disposable syringes for dispenser / diluter
361	PETRI dishes without triple vents ø 120	2903-03	Disposable syringes for dispenser / diluter
362	PETRI dishes ø 90 2 sections	2904-03	Disposable syringes for dispenser / diluter
363	PETRI dishes ø 90 3 sections	2905-03	Disposable syringes for dispenser / diluter
364	PETRI dishes ø 90 4 sections	2906-03	Disposable syringes for dispenser / diluter
365	PETRI dishes with triple vents. Ø 150	2907-03	Disposable syringes for dispenser / diluter
366	PETRI dishes with triple vents ø 90	2908-03	Disposable syringes for dispenser / diluter
367	PETRI dishes type Contact ø 65	2909-03	Disposable syringes for dispenser / diluter
368	PETRI dishes with triple vents ø 60	2910-03	Disposable syringes for dispenser / diluter
369	PETRI dishes square 120 x 120	2911-03	Disposable syringes for dispenser / diluter
380	PETRI dishes type Contact ø 65	2914-03	Disposable syringes for dispenser / diluter
434	Disconnecter – tapered dia. 4/5/6	2915-03	Disposable syringes for dispenser / diluter
435	Disconnecter – tapered dia. 6/7/8	2916	Adapter for 25 and 50 ml dispenser tips
436	Disconnecter – tapered dia. 8/9/10	2919-XX	Embedding rings
459	“T” untapered connector dia. 4	2921-XX	Tissue embedding cassettes
460	“T” untapered connector dia. 6	2922-XX	Tissue embedding sponge
461	“T” untapered connector dia. 8	2923	Metal trays for istology
462	“T” untapered connector dia. 10	2924	Metal trays for biopsy
463	“T” untapered connector dia. 12	2925-XX	Tissue embedding cassettes mega type
464	“Y” untapered connector dia. 4	2940-XX	Micro test tubes w/cap type Cobas Bio®
465	“Y” untapered connector dia. 6	3360	Petri dish without triple vents
466	“Y” untapered connector dia. 8	5550	Screw cap container 30 ml

Assicurazione Qualità

M. Capitanio

467	"Y" untapered connector dia. 10	5551	Screw cap container 30 ml Sterile
468	"Y" untapered connector dia. 12	5560	Screw cap container ml. 60
481	Container for urine	5561	Screw cap container ml.60 - sterile
482	Sampling bottle for urine 24h	5570	Urine container 60 ml
510	Straight connector tapered dia. 4/5/6	5571	Urine container 60 ml Sterile
511	Straight connector tapered dia. 6/7/8	5600	Test tubes w/screw caps 12 ml
512	Straight connector tapered dia. 8/9/10	5601	Test tubes w/screw caps 12 ml
513	Straight connector tapered dia. 10/11/12	5620	Test tubes w/screw caps 30 ml
514	Straight connector tapered dia. 12/13/14	5621	Screw cap container 30 ml
515	Straight connector tapered dia. 14/15/16	5622	Test tubes w/screw caps 30 ml
520	"T" untapered connector dia. 14	5623	Screw cap container 30 ml
521	"T" untapered connector dia. 16	5630	Container for urine and biological specimens 60 ml
525	"Y" untapered connector dia. 14	5631	Containers for urine and biological specimens 60 ml
526	"Y" untapered connector dia. 16	5632	Screw cap container 60 ml
527	120° "Y" tapered connector 4/5/6	5633	Faeces container 60 ml - sterile
528	120° "Y" tapered connector 6/7/8	5640	Urine Container 150 ml
529	120° "Y" tapered connector 8/9/10	5641	Urine container 150 ml Sterile
530	120° "Y" tapered connector 10/11/12	5642	Urine container 150 ml Aseptic
531	120° "Y" tapered connector 12/13/14	5643	Urine container 200 ml
532	120° "Y" tapered connector 14/15/16	5644	Urine container 200 ml Sterile
535	Disconnecter – tapered dia. 10/11/12	9585	Disposable tips for micropipettes rack (PP) type Gilson® yellow
536	Disconnecter – tapered dia. 12/13/14	9605	Disposable tips for micropipettes rack (PP) yellow UNIVERSAL
537	Disconnecter – tapered dia. 14/15/16	9615	Disposable tips for micropipettes rack (PP) blue UNIVERSAL
569	Safety Box	9624	Disposable tips for micropipettes rack (PP) type CRYSTAL
589	Douche cans or irrigators 1000 ml	9629	Sterile disposable filter tips for micropipettes rack 1:170µl tipo Univers.
590	Douche cans or irrigators 2000 ml	9630	Sterile disposable filter tips for micropipettes rack 1:10µl t.Eppendorf®
591	Male bed bottle with cap	9631	Sterile disposable filter tips for micropipettes rack 100:1000µl t. Epp.®
597	Sputum collection container ml 55	9632	Sterile disposable filter tips for micropipettes rack 5:100µl t. Eppend.®
600	Faeces containers 15 ml	9633	Sterile disposable filter tips for micropipettes rack 100:1000µl t.Gilson®
650	Test tube caps ø 10,5	9634	Sterile disposable filter tips for micropipettes rack 2:30µl t.Gilson® grad.
651	Test tube caps ø 14	9635	Sterile disposable filter tips for micropipettes rack 2:100µl t.Gilson®grad.
652	Test tube caps ø 14,5	9636	Sterile disposable filter tips for micropipettes rack 2:150µl t.Gilson®grad.
654	Test tube caps ø 16	9637	Sterile disposable filter tips for micropipettes rack 20:300µl Universale
655	Test tube caps ø 22	9638	Sterile disposable filter tips for micropipettes rack1:10µl t.Gilson®
		9639	Sterile disposable filter tips for micropipettes rack new Eppend.® 2:20µl

Assicurazione Qualità

M. Capitanio

656	Test tube caps ø 28	9640	Disposable tips for micropipettes rack t. Eppend.® blu 100:1000µl
657	Test tube caps ø 33	9641	Disposable tips for micropipettes rack t. Eppend.® yellow 5:200µl
658	Test tube caps ø 38	9642	Disposable tips for micropipettes rack type Gilson® giallo 5:200µl
659	Test tube caps ø 15,5	9643	Disposable tips for micropipettes rack type Gilson® grad. 2:200µl
877	Unequal straight connector tapered dia. 4/6/8/10/12	9644	Disposable tips for micropipettes rack type Biohit® 2:300µl
878	Unequal straight connector tapered dia. 4/6/8/12/14/16	9645	Disposable tips for micropipettes rack type Gilson® blu 100:1000µl
879	Unequal straight connector tapered dia. 8/10/12/14/16	9646	Disposable tips for micropipettes rack type Gilson® 1:10µl
933	Scintillation vials macro	9647	Disposable tips for micropipettes rack t. Eppend.® grad. 10:300µl
934	Scintillation vials micro	9648	Disposable tips for micropipettes rack MACRO 1:5 ml
935	Disposable tips for micropipettes Universal grad 10:300 µl	9649	Disposable tips for micropipettes rack MACRO 1:5 ml
936	Beakers for cell counters type COULTER®	9650	Disposable tips for micropipettes rack type Eppendorf® 1:10µl
939	Disposable stirrer for cuvettes	9657	Disposable tips for micropipettes rack type Eppendorf® 2:20µl
940	Disposable tips for micropipettes type Gilson® grad. bevelled 1:200 µl	9745	Refill 96 Disposable tips for micropipettes t. Eppendorf® Orig. 2:20µl
941	Disposable tips for micropipettes type Gilson® Orig. blu 100:1000 µl	9746	Refill 96 Disposable tips for micropipettes t. Eppendorf® Cristal Orig.
942	Disposable tips for micropipettes t. Finnpiette® Orig. 1000:5.000 µl	9747	Refill 96 Disposable tips for micropipettes type Gilson® Orig. 1:10µl
943	Disposable tips for micropipettes MACRO 2000:10.000 µl	9748	Refill 96 Disposable tips for micropipettes type Eppendorf® Blu Orig.
944	Disposable fine tips 200 µl	9749	Refill 96 Disposable tips for micropipettes type Eppendorf® Giallo
956	Disposable tips for micropipettes type Gilson® MICRO	9750	Refill 96 Disposable tips for micropipettes t. Gilson® yellow Orig. 5:200µl
957	Disposable tips for micropipettes blu (per 9615)	9751	Refill 96 Disposable tips for micropipettes t. Gilson® grad. 2:200µl
958	Disposable tips for micropipettes type Gilson® 5-200 µl	9752	Refill 96 Disposable tips for micropipettes type Biohit® Orig. 5:300µl
959	Disposable tips for micropipettes MICRO	9753	Refill 96 Disposable tips for micropipettes t. Gilson® blu Orig. 101:1000µl
960	Disposable tips for micropipettes yellow universal	9754	Refill 96 Disposable tips for micropipettes Univ. grad. neutral 10:300 µl
961	Disposable tips for micropipettes BLU UNIVERSAL	9847	Disposable tips for micropipettes rack MACRO 1:5 ml
963	Disposable tips for micropipettes OXFORD® green old model	9848	Disposable tips for micropipettes rack MACRO 2:10 ml
964	Disposable tips for micropipettes SUPER LARGE	12960-00	Pluripet II 2
965	Disposable tips for micropipettes type BECKMAN® 1000 µl	13000-00	Pluripet II 20
966	Disposable tips for micropipettes type OXFORD® slim line	13010-00	Pluripet II 200
967	Disposable tips for micropipettes type OXFORD® blu slim line	13190-00	Pluripet II 10
968	Disposable tips for micropipettes type MLA® 200 µl	13210-00	Pluripet II 100
969	Disposable tips for micropipettes type MLA® 1000 µl	13220-00	Pluripet II 1000

Assicurazione Qualità

M. Capitanio

971	Disposable tips for micropipettes 2-300 ul neutral	13250-00	Pluripet II 5000
979	Disposable tips for micropipettes type Eppendorf® Originali 2:20µl	13290-00	Pluripet II 10000
980	Disposable tips for micropipettes yellow graduated	70434	Disconnecter – tapered dia. 4/5/6
984	Disposable tips for micropipettes 1000-5000µl	70435	Disconnecter – tapered dia. 6/7/8
1250	“L” untapered connector dia. 4	70436	Disconnecter – tapered dia. 8/9/10
1251	“L” untapered connector dia. 6	70459	“T” untapered connector dia. 4
1252	“L” untapered connector dia. 8	70460	“T” untapered connector dia. 6
1253	“L” untapered connector dia. 10	70461	“T” untapered connector dia. 8
1254	“L” untapered connector dia. 12	70462	“T” untapered connector dia. 10
1255	“L” untapered connector dia. 14	70463	“T” untapered connector dia. 12
1256	“L” untapered connector dia. 16	70464	“Y” untapered connector dia. 4
1285-00	Minifix 5ul	70465	“Y” untapered connector dia. 6
1286-00	Minifix 10ul	70466	“Y” untapered connector dia. 8
1287-00	Minifix 20ul	70467	“Y” untapered connector dia. 10
1288-00	Minifix 25ul	70468	“Y” untapered connector dia. 12
1289-00	Minifix 50ul	70510	Straight connector tapered dia. 4/5/6
1290-00	Minifix 100ul	70511	Straight connector tapered dia. 6/7/8
1291-00	Minifix 200ul	70512	Straight connector tapered dia. 8/9/10
1292-00	Minifix 250ul	70513	Straight connector tapered dia. 10/11/12
1293-00	Minifix 500ul	70514	Straight connector tapered dia. 12/13/14
1294-00	Minifix 1000ul	70515	Straight connector tapered dia. 14/15/16
1295-10	Seropet	70520	“T” untapered connector dia. 14
1296-10	Pluripet PL 2	70521	“T” untapered connector dia. 16
1297-10	Labofill	70525	“Y” untapered connector dia. 14
1298	Micro test tubes w/cap type Vitraton Akes®	70526	“Y” untapered connector dia. 16
1299	Micro test tubes w/cap type Beckman®	70527	120° “Y” tapered connector 4/5/6
1300-10	Pluripet PL 20	70528	120° “Y” tapered connector 6/7/8
1301-10	Pluripet PL 200	70529	120° “Y” tapered connector 8/9/10
1303-10	Fixopet FX 5	70530	120° “Y” tapered connector 10/11/12
1304-10	Fixopet FX 10	70531	120° “Y” tapered connector 12/13/14
1305-10	Fixopet FX 20	70532	120° “Y” tapered connector 14/15/16
1306-10	Fixopet FX 25	70535	Disconnecter – tapered dia. 10/11/12
1307-10	Fixopet FX 50	70536	Disconnecter – tapered dia. 12/13/14
1308-10	Fixopet FX 100	70537	Disconnecter – tapered dia. 14/15/16
1309-10	Fixopet FX 200	70877	Unequal straight connector tapered dia. 4/6/8/10/12
1310-10	Fixopet FX 250	70878	Unequal straight connector tapered dia. 4/6/8/12/14/16
1311-10	Fixopet FX 500	70879	Unequal straight connector tapered dia. 8/10/12/14/16
1312-10	Fixopet FX 750	71250	“L” untapered connector dia. 4
1313-10	Fixopet FX 1000	71251	“L” untapered connector dia. 6
1319-10	Pluripet PL 10	71252	“L” untapered connector dia. 8

Assicurazione Qualità

M. Capitanio

1321-10	Pluripet PL 100	71253	"L" untapered connector dia. 10
1322-10	Pluripet PL 1001	71254	"L" untapered connector dia. 12
1325-10	Pluripet PL 5000	71255	"L" untapered connector dia. 14
1326-10	Multipet 8 channel	71256	"L" untapered connector dia. 16
1328-10	Multipet 12 channel	71410	4 way connector dia. 3,5
1329-10	Pluripet PL 10000	71411	4 way connector dia. 6
1369-00	Labofix 5ul	71412	4 way connector dia. 8
1370-00	Labofix 10ul	71413	4 way connector dia. 10
1371-00	Labofix 20ul	71414	4 way connector dia. 12
1372-00	Labofix 25ul	71415	4 way connector dia. 13,5
1373-00	Labofix 50ul	71416	4 way connector dia. 15
1374-00	Labofix 100ul	81937	Visible range cuvettes 4,5 ml
1375-00	Labofix 200ul	81938	Visible range cuvettes 1,5 ml SEMI-MICRO
1376-00	Labofix 250ul	81939	Cuvette 4,5 ml UV grade
1377-00	Labofix 500ul	81940	Visible range cuvettes 2,5 ml SEMI-MICRO
1378-00	Labofix 1000ul	81941	Cuvette 2,5 ml SEMI-MICRO UV grade
1387	Graduated conical centrifuge tubes TPX® ø16 10 ml	81948	Cuvette 1,5 ml SEMI-MICRO UV grade
1388	Graduated conical centrifuge tubes TPX® ø18 15 ml	81960	Visible range cuvettes 4,5 ml 4 o.f.
1410	4 way connector dia. 3,5	81961	Cuvette 4,5 ml 4 o.f. UV grade
1411	4 way connector dia. 6	82631	Cuvette for AMELUNG®
1412	4 way connector dia. 8	83002	Dosipet 2
1413	4 way connector dia. 10	83005	Dosipet 5
1414	4 way connector dia. 12	83010	Dosipet 10
1415	4 way connector dia. 13,5	83025	Dosipet 25
1416	4 way connector dia. 15	83050	Dosipet 50
1417-10	Multipet 8 channel	83060	Dosipet 60
1418-10	Multipet 12 channel	83070	Dosipet S 10
1500	Cobas® Mira segments in rack	83071	Dosipet S 25
1501	Cobas® Mira segments in rack	83072	Dosipet S 50
1563	Cylindrical jar with ribbed cap – ml. 70	83080	Dosipet HF 10
1564	Cylindrical jar with ribbed cap – ml. 120	83081	Dosipet HF 25
1565	Cylindrical jar with ribbed cap – ml. 250	83082	Dosilab ml. 2,5
1566	Cylindrical jar with ribbed cap – ml. 500	83083	Dosilab ml. 5
1567	Cylindrical jar with ribbed cap – ml. 1000	83084	Dosilab ml. 10
1568	Cylindrical jar with ribbed cap – ml. 2000	83085	Dosilab ml. 30
1569	Cylindrical jar with ribbed cap – ml. 5000	83086	Dosilab ml. 60
1937	Visible range cuvettes 4,5 ml	83087	Dosilab ml. 100
1938	Visible range cuvettes 1,5 ml SEMI-MICRO	83250	Digital Burettes
1939	Cuvette 4,5 ml UV grade	84000	Graduated conic. test tubes w/screw cap15ml
1940	Visible range cuvettes 2,5 ml SEMI-MICRO	84001	Graduated conical test tubes w/screw cap ster. 15 ml
1941	Cuvette 2,5 ml SEMI-MICRO UV grade	84002	Graduated conical test tubes w/screw cap50 ml

Assicurazione Qualità

M. Capitanio

1948	Cuvette 1,5 ml SEMI-MICRO UV grade	84003	Graduated conical test tubes w/screw cap ster. 50 ml
1960	Visible range cuvettes 4,5 ml 4 o.f.	84004	Conical grad test tubes with self-standing bottom w/screw cap 50 ml
1961	Cuvette 4,5 ml 4 o.f. UV grade	84005	Conical grad test tubes with self-standing bottom w/screw cap ster. 50 ml
1962	Cap for cuvettes 10x10	84006	Graduated conical centrifuge tube 15 ml. multipack
2300	Graduated conical centrifuge tubes ø1610 ml	84007	Graduated conical centrifuge tube 15 ml. multipack - sterile
2302	Graduated conical centrifuge tubes ø18 15 ml	84008	Graduated conical centrifuge tube 15 ml sterile
2501	Sample cups type Technicon® 1,5 ml	84009	Graduated conical centrifuge tube 15 ml. multipack on PS tray
2502	Sample cups type Technicon® 2 ml	84010	Graduated conical centrifuge tube 15 ml. multipack on PS tray - sterile
2503	Sample cups type Technicon RA 1000®	84011	Graduated conical centrifuge tube 50 ml. multipack
2508	Sample cups type Technicon® 4 ml	84012	Graduated conical centrifuge tube 50 ml. multipack - sterile
2510	Sample cups type CentrifChem® 0,25 ml	84013	Graduated conical centrifuge tube 50 ml sterile
2511	Sample cups type Gemsac® 0,5 ml	84014	Graduated conical centrifuge tube 50 ml. multipack on PS tray
2512	Anti-evaporation cap special cross cut cap for 2510/2511	84015	Graduated conical centrifuge tube 50 ml. multipack on PS tray - sterile
2514	Anti-evaporation cap push-on for 2501/2502/2503/2510/2512	84016	Test tube single pack ml. 15 - sterile
2580	Snap-on lid container – ml. 30	84017	Test tube single pack ml. 50 - sterile
2585	Snap-on lid container – ml. 35	84018	Test tube with base ml. 50 - sterile
2595	Feaces containers 60 ml	88108	Test tubes cylindrical 55x11,50
2596	Containers for urine and biological specimens 200 ml	88200	Test tube screw cap – ml. 15
2598	Specimen container ml 32	88201	Test tube screw cap – ml. 15 - sterile
2600	sample cups for HITACHI®	88202	Test tube screw cap – ml. 10 - sterile
2602	sample cups for Boehringer ES300®/ES600®	88205	Test tube screw cap – ml. 15
2618	Cuvetta type OLLI®	88207	Test tube screw cap – ml. 15 - sterile
2620	Microtitration plates wells	88210	Test tube screw cap – ml. 20
2620/S	Microtitration plates wells, Sterile	88212	Test tube screw cap – ml. 15 - sterile
2621	Microtitration plates wells U	88301	Test tubes conical 16x105
2621/S	Microtitration plates wells U, Sterile	88302	Test tubes cylindrical 16x100
2622	Microtitration plates wells V	88303-XX	Test tube with srew caps
2622/S	Microtitration plates wells V, Sterile	88304-XX	Test tube screw caps – sterile
2623	Lids for microtitratyion plates	88305-XX	Stopper for disposable test tubes ø15 17 mm
2623/S	Lids for microtitratyion plates, Sterile	88306	Cylindrical test tube 13x75
2630	Fibrintimer cuvettes	88307	Test tubes cylindrical 75x12
2631	Amelung cuvettes	88308	Test tubes cylindrical 150x24
2647	L-shape spreader - sterile	88309-XX	Stopper for disposable test tubes ø 11-13mm
2648	L-shape spreader - sterile	88310	Test tubes R.I.A. 11x70
2650	Inoculating loop and needle – ml. 1	88317	Test tubes cylindrical 87x12
2651	Inoculating loop and nedle – ml. 10	88318	Stopper for disposable test tubes for R.I.A.
2652	Rod for clot detachment and extraction	88319	Test tubes cylindrical 86x12

Assicurazione Qualità

M. Capitanio

2695	Feaces Containers 60 ml. Sterile	88320	Test tubes cylindrical 150x16
2696	Containers for urine and biological specimens 200 ml	88321	Stopper for disposable test tubes for 88320
2810	Serological Pipette 1 ML Sterile	88322	Test tubes cylindrical 99,30x16
2811	Serological Pipette 1 ML Sterile	88323	Test tubes cylindrical 100x16
2812	Serological Pipette 1 ML	88324	Test tubes conical 10 ml grad.
2813	Serological Pipette 2 ML Sterile	88325	Test tubes conical 10 ml grad.
2814	Serological Pipette 2 ML Sterile	88335	Pasteur pipette capillar
2815	Serological Pipette 2 ML	88336	Pasteur pipette 1 ml
2816	Serological Pipette 5 ML Sterile	88337	Pasteur pipette 3 ml
2817	Serological Pipette 5 ML Sterile	88435	PP test tube 12x86 mm with pink cap
2818	Serological Pipette 5 ML	88436	10 places numbered rack for SEDI RATE system
		88437	Pipette in PS

Assicurazione Qualità

M. Capitanio