

**TECHNICAL SPECIFICATIONS**

Trade name	Draining rack
Use	For laboratory draining
Material	Polystyrene PS

**DIMENSIONAL SPECIFICATIONS**

Dimension	mm. 450 x 630 x 110
Sides	72
General dimensional tolerance	mm. +/- 0,30

**GENERAL FEATURES**

CE mark	No
Autoclavability	No
Sterile	No
Food contact suitability of the raw material	Yes

**CHEMICAL RESISTANCE**

Categories of substances	Reference	Concen. %	T. 20°	T. 40°	T.° 60
Inorganic acid	Sulphuric acid	98	C	-	-
Organic acid	Benzoic acid	100	B	B	B
Alcohol	Ethanol	100	B	C	-
Aldehydes	Acetaldehyde	100	C	-	-
Inorganic alkali	Sodium hydroxide	50	A	A	A
Organic alkali	Aniline	100	C	-	-
Ketone	Acetone	100	-	-	-
Esther	Ethyl acetate	100	C	-	-
Hydrocarbons, halogenated	Dichloroethylene	100	C	-	-
Hydrocarbons, aromatic	Benzene	100	C	-	-
Hydrocarbons, linear	Hexane	100	C	-	-
Supersolvent	Tetrahydrofuran	100	C	-	-

**Legenda** (not applicable if neglected):

- **A: fair resistance**; exposure (30 days) to the chemical does not cause any damage.
- **B: sufficient resistance**; exposure causes damage of poor importance, which sometimes is only temporary.
- **C: poor resistance**; exposure to chemical is not allowed, or causes immediate permanent damage.

Results of testing with reference materials are to be considered as indication: in case of specific use it is recommended to carry out preliminary testing.

**PACKAGING**

Case	Pz. 1
------	-------

**Quality Assurance**

Massimiliano Capitanio

This paper is an electronic print-out: signature is not requested