

NEW SLIDE STORING SYSTEM

After many years of honored service, we decided to dismiss the well-known cod.930 and to present the new Slide Storing System (shown in the pictures below).

As the previous one, the new model is designed for the classification and storage of microscope slides. It is composed of 4 sliding drawers, a closing cover and four fixed support feet or, as alternative, four swivel castors to guarantee mobility.



This is a completely new storing system, made of impact-resistant PMMA and steel.

It is a real design object that will also be appreciated by users who are attentive to the aesthetic aspect.

The semi-transparent ice color gives an attractive appearance while the metal structure makes it stable, while ensuring the value of lightness, functionality, strength and durability.

Each single drawer can hold up to 13 supports for slides (code 921) and transparent cover plates (code 920), to be ordered separately.

Each slides support can contain 100 microscope slides, easily extractable. Each cavity is numbered. The complete Storing System can therefore contain up to 5,200 slides.

The transparent lids guarantee an easy check of the contents, protecting the slides from dust and other foreign bodies.

The Slide Storing System has the following ordering codes:

946L1 - version with four fixed support feet

The overall dimensions, including feet, are: 490x475x700h mm

945L1 - version with four swivel castors

The overall dimensions, including swivel castors, are: 490x475x780h mm

A Kit of Four swivel casters is available to be retrofitted on code 946L1. Ordering code **947L1** for any future need.

A Kit of Four fixed feet is available to be retrofitted on code 945L1. Ordering code **948L1** for any future need.

The sooner we will update the images and the presentation text of the product on the website.

Feel free to contact us for any further information you might need.

Best Regards,
Kartell Labware

DOMINIQUE DUTSCHER SAS