

## Stacker Microtube Rack

(#61048)

### 1. Description

SPL Stacker Microtube Racks are uniquely designed for handling and storage of Microcentrifuge Tubes.

### 2. Features

- ❖ Stackable feature for space saving
- ❖ 5 colors: Pink, Yellow, Green, Orange, Blue
- ❖ Non-autoclavable


### 3. Product Range

Cat. No.	Material	Description	Packaging
61048	ABS	6 x 8 ( 48 Holes)	10

### 4. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

### 5. Product Photographs




### 6. Specifications & Technical Drawing

- Specifications

Dimensions	Material (Resin)	Packaging
Length: 123.00 mm Width: 183.00 mm Height: 26.00 mm	ABS	▪ Pieces/Case: 10

- Technical Drawing



[www.spllifesciences.com](http://www.spllifesciences.com)

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: [spl@ispl.co.kr](mailto:spl@ispl.co.kr)

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: [business@ispl.co.kr](mailto:business@ispl.co.kr)

SPL-TDS-STTURA

*For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.*

© SPL Life Sciences Co., Ltd., Version 2, July 20<sup>th</sup>, 2016