

Omni Box (#80010)

1. Description

Omni Boxes of SPL provides a simple and convenient solution for storing small experiment-related equipments in a single container.

2. Features

- ❖ For the storage of magnetic bar, tubing, etc.

3. Product Range

Cat. No.	Material	Dimension l x w x h (mm)	Packaging
80100	PS	183.00 x 97.00 x 33.00	1

4. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

5. Product Photographs




6. Specifications & Technical Drawing

- Specifications

Dimensions	Material (Resin)	Packaging
Length: 86.40 mm Width: 175.60 mm Height: 30.00 mm	Polystyrene (PS)	▪ Pieces/Case: 1

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-OMNIBO

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

© SPL Life Sciences Co., Ltd., Version21, July 20th, 2016