

Microlab 625 Advanced Dual Syringe Diluter

Specification Sheet (Part/REF # ML625-DIL)



Quickly perform analytical dilutions of standards and samples. Minimize reagent and labor costs.

Product Specifications

Power Rating	24 VDC, 2.5A
Power Requirements	100-240V 1.5A max 50-60 Hz
Program Memory	2 GB
Syringe Drive Mechanism	Stepper motor driven high precision lead screw with encoder
Controller	Advanced
Number of Syringe Drives	Two
Probe	Concorde CT Hand Probe
Syringes Included	Includes two user selected syringes (default 250 μ L and 2.5 mL)

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/laboratory-products/microlab-600/ml625-dil>