

# Microlab 600 Universal Valve Cross Tube Assembly

Specification Sheet (Part/REF # 61498-01)



The Microlab Universal Valves are used on the Single Syringe Dispenser, Dual Syringe Diluter, Dual Syringe Dispenser and Continuous Dispenser. Adjust the Cross Tube and Valve Plugs to achieve the desired application.

## Product Specifications

<b>Threads</b>	1/4"-28 UNF
----------------	-------------

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/valves/61498-01>

**Spec. Version B**

DOMINIQUE DUTSCHER SAS