






# Méthanol

Synonyms:

Alcool méthylique

CH <sub>3</sub> OH	Classification	Danger
Molecular weight	32 ONU	1230  
CAS	67-56-1 Transport hazard class:	3 
EEC-N:	200-659-6 Packing	II

Méthanol

> RS - Pour HPLC Isocratique - ACS -  
Reag.Ph.Eur R2- Reag.USP

Grade Spécifique

Description	Clear colourless liquid	-
Identification (I.R.)	Positive	-
Density at 20°C	0.791 - 0.793	-
Refractive index at 20°C	1.3278 - 1.3298	-
Boiling point	64 - 65	?°C
Residue after evaporation	<= 5	ppm
Water (K.F.)	<= 0.05	%



Acidity	<= 0.0005	meq/g
Alcalinity	<= 0.0002	meq/g
Subst. darkened by H2SO4	Conform	-
Subs. reducing KMnO4	Conform	-
Water solubility	Conform	-
Carbonyl compounds (as CO)	<= 0.001	%
UV Absorbance ACS	Conform	-
HPLC Gradient (ACS)	Conform	-
Assay (GC)	>= 99.9	%
UV Transmittance	-	-
at 210 nm	>= 20	%
at 235 nm	>= 80	%
at 260 nm	>= 98	%
UV Absorbance	-	-
at 225 nm	<= 0.17	AU
Description	Clear colourless liquid	-
Identification (I.R.)	Positive	-
Density at 20°C	0.791 - 0.793	-
Refractive index at 20°C	1.3278 - 1.3298	-
Boiling point	64 - 65	°C



Residue after evaporation	<= 5	ppm
Water (K.F.)	<= 0.05	%
Acidity	<= 0.0005	meq/g
Alcalinity	<= 0.0002	meq/g
Subst. darkened by H2SO4	Conform	-
Subs. reducing KMnO4	Conform	-
Water solubility	Conform	-
Carbonyl compounds (as CO)	<= 0.001	%
UV Absorbance ACS	Conform	-
HPLC Gradient (ACS)	Conform	-
Assay (GC)	>= 99.9	%
UV Transmittance	-	-
at 210 nm	>= 20	%
at 235 nm	>= 80	%
at 260 nm	>= 98	%
UV Absorbance	-	-
at 225 nm	<= 0.17	AU

**Méthanol**

> RS - Pour HPLC Isocratique - ACS -  
Reag.Ph.Eur R2- Reag.USP

**Grade Spécifique**

ID	Size	Packaging
000000000000412532	2.5 l	Flacon verre
000000000000412535	2.5 l	Flacon verre sécurisé
000000000000412531	1 l	Flacon verre sécurisé
000000000000412533	1 l	Flacon verre



**CARLO ERBA Reagents S.A.S.**  
Chaussée du Vexin – 27106 Val de  
Reuil cedex  
N° TVA: FR 63391048824  
Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



*CARLO ERBA Reagents  
operates with a Certified  
Quality Management  
System*