

# Standard Analysis Certificate

**PRODUCT** :SILVER NITRATE 1 N (1 mol/L) RPE**CODE** :424030000

TEST	U.M.	SPECIFICATION
Description Assay (potentiometry)	- N	Clear colourless liquid 0.999 - 1.001

This volumetric solution is checked by sodium chloride volumetric standard, traceable to the SRM 399b of NIST (National Institute of Standards and Technology)

**Date**

:11/01/2020

**QUALITY CONTROL RESPONSIBLE****B. COULANGE (VDR)**