

Standard Analysis Certificate

PRODUCT :CHLOROFORM STAB.AMYLENE RS SPECTROSOL
CODE :438590

TEST	U.M.	SPECIFICATION
Clear, colourless liq.appearance	-	Conform
Identification	-	Conform
Colour	Apha	<= 10
Refractive index at 20°C	-	1.444 - 1.448
Water content (K.F.)	mg/Kg	<= 100
Free acid (as HCl)	mg/Kg	<= 5
Non volatile residue	mg/Kg	<= 5
Stabilizer (Amylene)	mg/Kg	5 - 50
Carbon tetrachloride	mg/Kg	<= 100
Dichloromethane	mg/Kg	<= 100
UV transmittance at 250 nm	%	>= 50
UV transmittance at 260 nm	%	>= 85
UV transmittance at 280 nm	%	>= 95
Assay (GC)	%	>= 99.9

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE
B. COULANGE (VDR)