

# Standard Analysis Certificate

**PRODUCT** :HEPTANE Mixture of isomers RS PESTIPUR - For pesticide analysis  
**CODE** :446840  
**METHOD** :LAB118

TEST	U.M.	SPECIFICATION
Refractive index at 20°C	-	1.387 - 1.407
Water content (K.F.)	mg/Kg	<= 50
Non volatile residue	mg/Kg	<= 5
Colour	Hazen	<= 10
Aromatic compounds	mg/Kg	<= 100
GC-ECD.Individual peak (Lindane)	ng/l	<= 3
Retention time trichlorobenzene to mirex	-	-

For chlorinated compounds analysis; extraction of pesticides.

Date

:11/01/2020

**QUALITY CONTROL RESPONSIBLE**

**B. COULANGE (VDR)**