

# Standard Analysis Certificate

**PRODUCT** :n-HEXANE 99 % RS SPECTROSOL  
**CODE** :447050  
**METHOD** :447050

TEST	U.M.	SPECIFICATION
Clear,colourless liq.appearance	-	Conform
Identification	-	Conform
Colour	Apha	<= 10
Refractive index at 20°C	-	1.373 - 1.377
Water content (K.F.)	mg/Kg	<= 150
Total sulphur (S)	ppm	<= 5
Aromatic compounds	mg/Kg	<= 5
Non volatile residue	mg/Kg	<= 5
Assay (GC)	%	>= 99
UV transmittance at 210 nm	%	>= 60
UV transmittance at 215 nm	%	>= 70
UV transmittance at 220 nm	%	>= 80
UV transmittance at 230 nm	%	>= 94
UV transmittance at 245 nm	%	>= 98

DOMINIQUE DUTSCHER SAS

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

B. COULANGE (VDR)