

Standard Analysis Certificate

PRODUCT :ETHYLE ACETATE RS SPECTROSOL
CODE :448270
METHOD :8509

TEST	U.M.	SPECIFICATION
Description	-	Clear liquid
Colour (APHA)	-	<=10
Identification	-	Positive
Density at 20° C	-	0.898 ÷ 0.902
Refractive index at 20°C	-	1.3699 ÷ 1.3739
Boiling point	° C	76.9 ÷ 77.4
Water (K.F.)	ppm	<=100
Residue on evaporation	ppm	<=10
Acidity	meq/g	<=0.0005
Alcalinity	meq/g	<=0.0002
Ethyl alcohol	ppm	<= 100
Methyl alcohol	ppm	<=100
Assay (GLC)	%	>=99.8
U.V. Transmittance	-	-
at 255 nm	%	>=15
at 260 nm	%	>=75
at 270 nm	%	>=95
at 280 nm	%	>=98

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

B. COULANGE (VDR)