

EC Declaration of Conformity

Conformity to DIRECTIVE 98/79/EC OF THE EUROPEAN PARLIAMENT AND
OF THE COUNCIL of 27 October 1998 on in vitro diagnostic medical devices

Manufacturer: Greiner Bio-One GmbH
Bad Haller Str. 32
4550 Kremsmünster
Austria

Production Location: Greiner Bio-One GmbH
Bad Haller Str. 32
4550 Kremsmünster
Austria

Greiner Bio-One
North America Inc.
4238 Capital Drive
Monroe
NC 28110
United States of America

Greiner Bio-One Brasil
Produtos Medicos
Hospitalares Ltda.
Av. Affonso Pansan no. 1.967
13473-620
Vila Bertini Americana
Sao Paulo
Brasil

Greiner Bio-One (Thailand) Ltd.
700/172 Moo. 1
Amata Nakorn Industrial Estate
Tambon Bankao
Amphur Phanthong
Chonburi 20160
Thailand

Product/ Product Group: VACUETTE® TUBES
(for details please refer to page 2-33)
Classification: Other device (all devices except Annex II and except self-testing devices)
GMDN Codes (for details please refer to page 2-33)

We herewith declare under our sole responsibility that the above mentioned products meet the provisions of the following EC Council Directives and applicable Standards. All supporting documentations are retained under the premises of the manufacturer.

Conformity Assessment procedure acc. to: Annex III of the Directive 98/79/EC of the European Parliament and of the Council on in vitro diagnostic medical devices

Standards:
Refer to the List of applicable (harmonized) standards in the Technical Documentation

Kremsmünster, 16.08.2018



Signature:

Georg Sambs
Reg. Affairs Manager

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
4NC	ESR sodium citrate 3.2%	VACUETTE® TUBE 2 ml 4NC ESR sodium citrate 3,2% 13x75 black cap-black ring, PET, PREMIUM	454073	42585
4NC	ESR sodium citrate 3.2%	VACUETTE® TUBE 2 ml 4NC ESR sodium citrate 3,2% 13x75 black cap-black ring, PET, special label, PREMIUM	454073CN	42585
4NC	ESR sodium citrate 3.2%	VACUETTE® TUBE 2 ml 4NC ESR sodium citrate 3,2% 13x75 black cap-black ring, PET, special label, PREMIUM, SGH	454073SGH	42585
4NC	ESR sodium citrate 3.2%	VACUETTE® TUBE 1.5 ml 4NC ESR sodium citrate 3.2% 13x75 black cap-black ring, PREMIUM	454284	42585
4NC	ESR sodium citrate 3.2%	VACUETTE® TUBE 2 ml 4NC ESR sodium citrate 3.2% 13x75 black cap-black ring, barcode label, PREMIUM	484516	42585
4NC	ESR sodium citrate 3.2%	VACUETTE® TUBE 2.75 ml 4NC ESR sodium citrate 3.2% 9x120 black cap, PP	729070	42585
4NC	ESR sodium citrate 3.2%	VACUETTE® TUBE 1.5 ml 4NC ESR sodium citrate 3.2% 9x120 black cap, PP	729073	42585
4NC	ESR sodium citrate 3.2%	VACUETTE® TUBE 1.5 ml 4NC ESR sodium citrate 3.2% 9x120 black cap, PP, High Altitude	729073H	42585
4NC	ESR sodium citrate 3.2%	VACUETTE® TUBE 2.9 ml 4NC ESR sodium citrate 3.2% 9x120 black cap, glass	729090	42585
4NC	ESR sodium citrate 3.2%	VACUETTE® TUBE 1.6 ml 4NC ESR sodium citrate 3.2% 9x120 black cap, glass	729093	42585
4NC	ESR sodium citrate 3.2%	VACUETTE® TUBE 1.6 ml 4NC ESR sodium citrate 3.2% 9x120 black cap, barcode label, glass	769093	42585
4NC	ESR sodium citrate 3.2%	VACUETTE® TUBE 1.5 ml 4NC ESR sodium citrate 3.2% 9x120 black cap, PP, Barcode Label	769094	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 1 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, PREMIUM	454320	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, PREMIUM	454321	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, PREMIUM, China	454321CN	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, special label, PREMIUM, SGH	454321SGH	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, non-ridged	454322	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, non-ridged, High Altitude	454322H	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, non-ridged, China	454322CN	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, non-ridged, High Altitude	454323	42585

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-black ring, sandwich tube, PREMIUM	454325	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 4.5 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-black ring, non-ridged	454326	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-black ring, sandwich tube, PREMIUM	454327	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 4.5 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-black ring, PREMIUM	454329	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3.2% 13x75 light blue cap-black ring, sandwich tube, non-ridged, High Altitude	454331	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-black ring, sandwich tube, non-ridged	454332	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, special label, non-ridged, CSS	454332CSS	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3 ml 9NC Coagulation sodium citrate 3,2% 13x75 blue cap-black ring, sandwich tube, transparent label, non-ridged	454333	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3 ml 9NC Coagulation sodium citrate 3,2% 13x75 blue cap-black ring, sandwich tube, transparent special label, non-ridged	454333CN	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3 ml 9NC Coagulation sodium citrate 3,2% 13x75 blue cap-black ring, sandwich tube, non-ridged	454334	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3ml 9NC Coagulation sodium citrate 3.2% 13x75 light blue cap-black Ring, sandwich tube, non-ridged, Canada	454334CA	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3 ml 9NC Coagulation sodium citrate 3,2% 13x75 blue cap-black ring, sandwich tube, non-ridged, China	454334CN	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-black ring, sandwich tube, non-ridged, High- Altitude	454337H	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3 ml 9NC Coagulation sodium citrate 3.2% 13x75 light blue cap-black ring, sandwich tube, non-ridged, High Altitude	454351	42585

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, non-ridged	454352	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 1.8 ml 9NC Coagulation sodium citrate 3.2% 13x75 black cap-white ring, sandwich tube, transparent special label, with circle fill line, non-ridged, Japan	454354JP	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, transparent label, non-ridged	454360	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 1.8 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, transparent special label, non-ridged, Japan	454360JP	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-black ring, sandwich tube, transparent label, non-ridged	454471	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 9 ml 9NC Coagulation sodium citrate 3,2% 16x100 blue cap-black ring, non-ridged	455322	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, transparent label, PREMIUM	474321	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 1.8 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, transparent special label, PREMIUM, Japan	474321JP	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 1.8 ml 9NC Coagulation sodium citrate 3.2% 13x75 black cap-white ring, sandwich tube, transparent special label, PREMIUM, Japan	474322JP	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-black ring, sandwich tube, transparent label, PREMIUM	474327	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3 ml 9NC Coagulation sodium citrate 3,2% 13x75 blue cap-black ring, sandwich tube, barcode label, non- ridged	484334	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3 ml 9NC Coagulation sodium citrate 3,2% 13x75 blue cap-black ring, sandwich tube, barcode label, non- ridged	484508	42585

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-white ring, sandwich tube, barcode label, G- Code, PREMIUM	484517	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3.2% 13x75 blue cap-black ring, sandwich tube, barcode label, G- Code, PREMIUM	484518	42585
9NC	Coagulation sodium citrate 3.2%	VACUETTE® TUBE 5 ml 9NC Coagulation sodium citrate 3.2% 13x100 blue cap-black ring, barcode labe, G-Code, PREMIUM	486505	42585
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3,8% 13x75 blue cap-white ring, sandwich tube, PREMIUM	454381	42585
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3,8% 13x75 blue cap-white ring, sandwich tube, PREMIUM, China	454381CN	42585
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3,8% 13x75 blue cap-white ring, sandwich tube, non-ridged	454382	42585
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 2 ml 9NC Coagulation sodium citrate 3,8% 13x75 blue cap-white ring, sandwich tube, special label, PREMIUM, High Altitude	454383CN	42585
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 3 ml 9NC Coagulation sodium citrate 3,8% 13x75 blue cap-black ring, sandwich tube, PREMIUM	454385	42585
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 3 ml 9NC Coagulation sodium citrate 3,8% 13x75 blue cap-black ring, sandwich tube, PREMIUM, China	454385CN	42585
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 4.5 ml 9NC Coagulation sodium citrate 3,8% 13x75 blue cap-black ring, non-ridged	454386	42585
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3,8% 13x75 blue cap-black ring, sandwich tube, PREMIUM	454387	42585
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 4.5 ml 9NC Coagulation sodium citrate 3,8% 13x75 blue cap-black ring, PREMIUM	454389	42585
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3,8% 13x75 blue cap-black ring, sandwich tube, non-ridged	454392	42585
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 3 ml 9NC Coagulation sodium citrate 3,8% 13x75 blue cap-black ring, sandwich tube, non-ridged	454393	42585
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 9 ml 9NC Coagulation sodium citrate 3,8% 16x100 blue cap-black ring, non-ridged	455382	42585

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
9NC	Coagulation sodium citrate 3.8%	VACUETTE® TUBE 3.5 ml 9NC Coagulation sodium citrate 3.8% 13x75 blue cap-black ring, sandwich tube, barcode label, non- ridged	484524	42585
AH	Ammonium Heparin	VACUETTE® TUBE 9 ml AH Ammonium Heparin 16x100 green cap-red ring, non-ridged	455031	47586
CAT	Crossmatch Serum Clot Activator	VACUETTE® TUBE 9 ml CAT Crossmatch Serum Clot Activator 16x100 pink cap-black ring, special label, non-ridged	455093	42386
CAT	Crossmatch Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Crossmatch Serum Clot Activator 13x100 pink cap-black ring, special label non-ridged	456097	42386
CAT	Crossmatch Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Crossmatch Serum Clot Activator 13x100 pink cap-black ring, special label, non-ridged	456297	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4.5 ml CAT Serum Clot Activator 13x75 red cap-black ring, PREMIUM	454027	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 red cap-black ring, PREMIUM	454092	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 red cap-black ring, PREMIUM, China	454092CN	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 3 ml CAT Serum Clot Activator 13x75 red cap-black ring, PREMIUM	454095	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 3 ml CAT Serum Clot Activator 13x75 red cap-black ring, PREMIUM, China	454095CN	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 2 ml CAT Serum Clot Activator 13x75 red cap-white ring, PREMIUM	454096	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 1 ml CAT Serum Clot Activator 13x75 red cap-white ring, PREMIUM	454098	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4.5 ml CAT Serum Clot Activator 13x75 red Snap Cap, without vacuum, non sterile, PREMIUM	454102	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 white cap-black ring, PREMIUM	454200	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 yellow cap-black ring, PREMIUM	454201	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 pink cap-black ring, PREMIUM	454202	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 red cap-black ring, non-ridged	454204	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 red cap-black ring, non-ridged, High Altitude	454204H	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 red cap-black ring, non-ridged, China	454204CN	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 royal blue cap-black ring, non-ridged	454206	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 2 ml CAT Serum Clot Activator 13x75 red cap-white ring, non-ridged	454236	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 pink cap-black ring, non-ridged	454262	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 orange cap-black ring, non-ridged	454272	42386

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 orange cap-black ring, PREMIUM	454278	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4,5 ml CAT Serum Clot Activator 13x75 red cap-black ring, non-ridged	454298	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 gold cap-royalblue ring, non-ridged	454313	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 gold cap-black ring, non-ridged	454314	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 red cap-black ring, transparent label, non-ridged	454433	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 red cap-black ring, special label, non-ridged	454451	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 2 ml CAT Serum Clot Activator 13x75 red cap-white ring, transparent label, non-ridged	454461	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Clot Activator 16x100 red cap-black ring, transparent label, non-ridged	455009	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Clot Activator 16x100 red cap-black ring, non-ridged	455092	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Clot Activator 16x100 red cap-black ring, non-ridged, High Altitude	455092H	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Clot Activator 16x100 red cap-black ring, non-ridged, China	455092CN	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Clot Activator 16x100 red cap-black ring, special label, non-ridged, Japan	455092JP	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Clot Activator 16x100 red cap-black ring, special label, non-ridged	455092CSS	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Clot Activator 16x100 red cap-black ring, special label, non-ridged, SGH	455092SGH	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Clot Activator 16x100 gold cap-royal blue ring, non-ridged	455202	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Clot Activator 16x100 gold cap-gold ring, non-ridged	455204	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Clot Activator 16x100 green cap-black ring, non-ridged	455205	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Clot Activator 16x100 red cap-black ring, special label, non-ridged	455228	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 pink cap-black ring, pink label, PREMIUM	456014	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 gold cap-black ring, black label, PREMIUM	456015	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, transparent label, non-ridged	456048	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, transparent label, non-ridged, Korea	456048KR	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-white ring, non-ridged	456066	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, non-ridged	456089	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, non-ridged, High Altitude	456089H	42386

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, non-ridged, China	456089CN	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, non-ridged, Korea	456089KR	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 orange cap-black ring, non-ridged	456091	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, PREMIUM	456092	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 white cap-black ring, non-ridged	456245	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 pink cap-black ring, non-ridged	456248	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-red ring, non-ridged	456249	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, special label, non-ridged	456304	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4.5 ml CAT Serum Clot Activator 13x75 red cap-black ring, transparent label, PREMIUM	474027	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 red cap-black ring, transparent label, PREMIUM	474092	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, transparent label, PREMIUM	476092	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, transparent special label, PREMIUM, Japan	476092JP	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 gold cap-gold ring, transparent label, PREMIUM	476093	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 2 ml CAT Serum Clot Activator 13x75 red cap-white ring, barcode label, PREMIUM	484515	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Clot Activator 13x75 red cap-black ring, PREMIUM, G-Barcode	484525	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Clot Activator 16x100 red cap-black ring, barcode label, G-Code, non-ridged	485504	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, barcode label, non-ridged	486008	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, barcode label, non-ridged	486501	42386
CAT	Serum Clot Activator	VACUETTE® TUBE 6 ml CAT Serum Clot Activator 13x100 red cap-black ring, barcode label, G-Code, PREMIUM	486504	42386
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 white cap-yellow ring, black label, non-ridged	454004	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 2.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, PREMIUM	454028	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 2.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, PREMIUM, China	454028CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-red ring, non-ridged	454066	41128

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, non-ridged	454067	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, non-ridged, USA	454067P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, non-ridged	454067R	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, non-ridged, China	454067CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, non-ridged, High Altitude, USA	454067PH	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, PREMIUM	454071	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, PREMIUM, China	454071CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, PREMIUM, Hong Kong	454071HK	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, special label, PREMIUM, OMB	454071OMB	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, non-ridged	454078	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, special label, non-ridged	454078CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 2.5 ml CAT Serum Separator Clot Activator 13x75 yellow cap-white ring, black label, non-ridged	454205	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 pink cap-yellow ring, pink label, PREMIUM	454212	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 gold cap-gold ring, PREMIUM	454214	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 2.5 ml CAT Serum Separator Clot Activator 13x75 red cap-white ring, transparent label, non-ridged	454219	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 2.5 ml CAT Serum Separator Clot Activator 13x75 red cap-white ring, transparent special label, non-ridged	454219CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, non-ridged	454225	41128

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, non-ridged	454225P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 gold cap-gold ring, non-ridged	454226	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 gold cap-gold ring, black label, non-ridged	454228	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 gold cap-gold ring, black label, non-ridged, USA	454228P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 gold cap-gold ring, black label, non-ridged, China	454228CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 2.5 ml CAT Serum Separator Clot Activator 13x75 red cap-white ring, non-ridged	454243	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 2.5 ml CAT Serum Separator Clot Activator 13x75 red cap-white ring, non-ridged, USA	454243P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 2.5 ml CAT Serum Separator Clot Activator 13x75 red cap-white ring, non-ridged, China	454243CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 gold cap-white ring, black label, non-ridged	454253	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 white cap-yellow ring, PREMIUM	454260	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 royal blue cap-yellow ring, non-ridged	454263	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 royal blue cap-yellow ring, non-ridged	454263P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 yellow cap-yellow ring, non-ridged	454264	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 gold cap-black ring, non-ridged	454319	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 gold cap-gold ring, transparent label, non-ridged	454420	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 gold cap-gold ring, transparent label, non-ridged, China	454420CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 2.5 ml CAT Serum Separator Clot Activator 13x75 gold cap-gold ring, transparent label, non-ridged	454473	41128

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, non-ridged	455010	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, non-ridged, USA	455010P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, non-ridged, China	455010CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Separator Clot Activator 16x100 gold cap-gold ring, non-ridged	455034	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Separator Clot Activator 16x100 gold cap-gold ring, non-ridged	455034P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8.5 ml CAT Serum Separator Clot Activator 16x100 gold cap-gold ring, non-ridged	455035	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, non-ridged	455057	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, transparent label, non-ridged	455059	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, transparent special label, non-ridged, Japan	455059JP	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, non-ridged	455071	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, non-ridged, Gamma	455071G	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, non-ridged	455071P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, non-ridged	455071R	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, special label, non-ridged, Korea	455071KR	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 yellow cap-black ring, non-ridged	455077	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 white cap-black ring, black label, non-ridged	455078	41128

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 7 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, non-ridged, double gel	455099	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 7 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, non-ridged, double gel, USA	455099P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, transparent label, non-ridged	455211	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, transparent special label, non-	455211KR	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Separator Clot Activator 16x100 gold cap-yellow ring, transparent label, non-ridged	455215	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Separator Clot Activator 16x100 white cap-yellow ring, transparent label, non-ridged	455231	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, non-ridged	455232	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, non-ridged, USA	455232P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8.5 ml CAT Serum Separator Clot Activator 16x100 pink cap-yellow ring, transparent label, non-ridged	455335	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 white cap-yellow ring, black label, non-ridged	456004	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 yellow cap-yellow ring, black label, non-ridged	456005	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 yellow cap-yellow ring, black special label, non-ridged	456005CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 royal blue cap-yellow ring, blue label, non-ridged	456006	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, PREMIUM	456010	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, black label, non-ridged	456017	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, black label, non-ridged	456018	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, black label, non-ridged	456018P	41128

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, black special label, non-ridged	456018CA	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-yellow ring, transparent special label, with circle fill line, PREMIUM, Japan	456018JP	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, black label, non-ridged, High Altitude, USA	456018PH	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, black special label, non-ridged	456018PCA	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 pink cap-yellow ring, non-ridged	456019	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 black cap-yellow ring, non-ridged	456044	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 black cap-yellow ring, non-ridged, USA	456044P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 orange cap-black ring, black label, non-ridged	456045P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 yellow cap-yellow ring, PREMIUM	456049	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, transparent label, non-ridged	456069	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, transparent special label, non- ridged	456069CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, transparent special label, non- ridged	456069JP	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, transparent label, non-ridged, Korea	456069KR	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, PREMIUM	456071	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, PREMIUM	456071P	41128

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, PREMIUM, China	456071CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, non-ridged	456073	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, non-ridged	456073P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, non-ridged	456073R	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, non-ridged, China	456073CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, non-ridged, Korea	456073KR	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, non-ridged, High Altitude, USA	456073PH	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 grey cap-yellow ring, non-ridged	456076P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 yellow cap-black ring, black label, PREMIUM	456077	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 white cap-black ring, black label, PREMIUM	456078	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 white cap-black ring, black special label, PREMIUM	456078CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 grey cap-yellow ring, PREMIUM	456208	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 royal blue cap-yellow ring, transparent label, non-ridged	456210	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, transparent label, non-ridged	456234	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, transparent label, non-ridged, China	456234CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-red ring, non-ridged	456239	41128

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-yellow ring, non-ridged, USA	456254	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 white cap-white ring, black label, non-ridged	456278	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 pink cap-black ring, special label, non-ridged	456288	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 yellow cap-black ring, special label, non-ridged	456289	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-black ring, special label, non-ridged	456290	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, non-ridged	456292	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 4 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, non-ridged, USA	456292P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 black cap-yellow ring, non-ridged	456293	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 black cap-yellow ring, non-ridged, USA	456293P	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-white ring, non-ridged	456295	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, transparent label, PREMIUM	474071	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, transparent label, PREMIUM, China	474071CN	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 gold cap-gold ring, transparent label, PREMIUM	474214	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, transparent label, PREMIUM	476010	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-yellow ring, transparent special label, with circle fill line, PREMIUM, Japan	476010JP	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, transparent label, PREMIUM	476071	41128

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, transparent special label, PREMIUM, Japan	476071JP	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 2.5 ml CAT Serum Separator Clot Activator 13x75 red cap-white ring, barcode label, non-ridged	484250	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 2.5 ml CAT Serum Separator Clot Activator 13x75 yellow cap-white ring, black barcode label, non-ridged	484485	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, barcode label, non-ridged	484502	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 gold cap-gold ring, barcode label, non-ridged	484503	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 2.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, barcode label, G-Code, PREMIUM	484513	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 3.5 ml CAT Serum Separator Clot Activator 13x75 red cap-yellow ring, barcode label, G-Code, PREMIUM	484514	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, barcode tube, non-ridged	485020	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator 16x100 red cap-yellow ring, barcode label, G-Code, non-ridged	485502	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 9 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, barcode label, non-ridged	485509	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 white cap-yellow ring, barcode label, non-ridged	486004	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, barcode label, non-ridged	486005	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 royal blue cap-yellow ring, barcode label, non-ridged	486006	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 yellow cap-yellow ring, barcode label, non-ridged	486007	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, black barcode label, non-ridged	486018	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 pink cap-yellow ring, barcode label, non-ridged	486019	41128

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 gold cap-gold ring, black barcode label, non-ridged	486500	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, barcode label, G-Code, PREMIUM	486503	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, bi-barcode label, non-ridged	495071	41128
CAT	Serum Sep Clot Activator	VACUETTE® TUBE 5 ml CAT Serum Separator Clot Activator 13x100 red cap-yellow ring, bi-barcode label, PREMIUM	496503	41128
CAT	Serum Sep Clot Activator FAST	VACUETTE® TUBE 3.5 ml CAT Serum Fast Separator 13x75 orange cap-yellow ring, non-ridged	454591	41128
CAT	Serum Sep Clot Activator FAST	VACUETTE® TUBE 5 ml CAT Serum Fast Separator 13x100 orange cap-yellow ring, non-ridged	456310	41128
FE	Glucose Test only	VACUETTE® TUBE 2 ml FE Glucose Test only 13x75 grey cap-white ring, PREMIUM	454211	58090
FE	Glucose Test only	VACUETTE® TUBE 2 ml FE Glucose Test only 13x75 grey cap-white ring, non-ridged	454252	58090
FE	Sodium Fluoride/ K3EDTA	VACUETTE® TUBE 4 ml FE Sodium Fluoride / K3EDTA 13x75 grey cap-black ring, non-ridged	454033	58090
FE	Sodium Fluoride/ K3EDTA	VACUETTE® TUBE 2 ml FE Sodium Fluoride / K3EDTA 13x75 grey cap-white ring, PREMIUM	454085	58090
FE	Sodium Fluoride/ K3EDTA	VACUETTE® TUBE 4 ml FE Sodium Fluoride / K3EDTA 13x75 grey cap-black ring, PREMIUM	454091	58090
FE	Sodium Fluoride/ K3EDTA	VACUETTE® TUBE 2 ml FE Sodium Fluoride / K3EDTA 13x75 grey cap-white ring, non-ridged	454221	58090
FE	Sodium Fluoride/ K3EDTA	VACUETTE® TUBE 2 ml FE Sodium Fluoride / K3EDTA 13x75 grey cap-white ring, transparent label, PREMIUM	454344	58090
FE	Sodium Fluoride/ K3EDTA	VACUETTE® TUBE 4 ml FE Sodium Fluoride / K3EDTA 13x75 grey cap-black ring, barcode label, non-ridged	484033	58090
FE	Sodium Fluoride/ K3EDTA	VACUETTE® TUBE 2 ml FE Sodium Fluoride / K3EDTA 13x75 grey cap-white ring, barcode label G-Code, PREMIUM	484520	58090
FE	Sodium Fluoride/ K3EDTA	VACUETTE® TUBE 4 ml FE Sodium Fluoride / K3EDTA 13x75 grey cap-black ring, PREMIUM, G-Barcode	484527	58090
FH	Sodium Fluoride/ Sodium Heparin	VACUETTE® TUBE 2 ml FH Sodium Fluoride / Sodium Heparin 13x75 grey cap-white ring, PREMIUM	454218	42836

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
FH	Sodium Fluoride/ Sodium Heparin	VACUETTE® TUBE 2 ml FH Sodium Fluoride / Sodium Heparin 13x75 grey cap-white ring, PREMIUM, China	454218CN	42836
FH	Sodium Fluoride/ Sodium Heparin	VACUETTE® TUBE 2 ml FH Sodium Fluoride / Sodium Heparin 13x75 grey cap-white ring, non-ridged	454430	42836
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 2 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-white ring, PREMIUM	454061	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 2 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-white ring, PREMIUM, China	454061CN	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 2 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-white ring, PREMIUM, Japan	454061JP	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 4 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-black ring, PREMIUM	454062	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 4 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-black ring, PREMIUM, China	454062CN	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 3 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-black ring, PREMIUM	454094	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 2 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-white ring, transparent label, non-ridged	454231	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 2 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-white ring, transparent special label, non- ridged	454231CN	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 2 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-white ring, non-ridged	454238	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 2 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-white ring, special label, non-ridged, Japan	454238JP	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 4 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-black ring, non-ridged	454297	47591

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 4 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-black ring, non-ridged, Canada	454297CA	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 2 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-yellow ring, non-ridged	454427	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 4 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-black ring, non-ridged, additionally french text on label	454437	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 4 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-black ring, transparent label, non-ridged	454438	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 2 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-white ring, special label, non-ridged	454453	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 5,5 ml FX Sodium Fluoride / Potassium Oxalate 13x100 grey cap-black ring, non-ridged	456020	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 5,5 ml FX Sodium Fluoride / Potassium Oxalate 13x100 grey cap-black ring, french label, non-ridged	456020CA	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 6 ml FX Sodium Fluoride / Potassium Oxalate 13x100 grey cap-black ring, PREMIUM	456062	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 4 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-black ring, barcode label, non-ridged	484505	47591
FX	Sodium Fluoride/ Potassium Oxalate	VACUETTE® TUBE 2 ml FX Sodium Fluoride / Potassium Oxalate 13x75 grey cap-white ring, barcode label, PREMIUM	484519	47591
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-black ring, non-ridged	454009	43865
K2E	K2EDTA	VACUETTE® TUBE 3 ml K2E K2EDTA 13x75 lavender cap-black ring, PREMIUM	454020	43865
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-black ring, PREMIUM	454023	43865
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-black ring, special label, PREMIUM, Japan	454023JP	43865
K2E	K2EDTA	VACUETTE® TUBE 2 ml K2E K2EDTA 13x75 lavender cap-white ring, PREMIUM	454024	43865

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
K2E	K2EDTA	VACUETTE® TUBE 2 ml K2E K2EDTA 13x75 lavender cap-white ring, PREMIUM, China	454024CN	43865
K2E	K2EDTA	VACUETTE® TUBE 2 ml K2E K2EDTA 13x75 lavender cap-white ring, special label, PREMIUM, Japan	454024JP	43865
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-yellow ring, PREMIUM	454032	43865
K2E	K2EDTA	VACUETTE® TUBE 2 ml K2E K2EDTA 13x75 lavender cap-white ring, non-ridged	454047	43865
K2E	K2EDTA	VACUETTE® TUBE 1 ml K2E K2EDTA 13x75 lavender cap-white ring, non-ridged	454052	43865
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 pink cap-black ring, non-ridged	454208	43865
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-black ring, non-ridged	454209	43865
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-black ring, non-ridged, High Altitude	454209H	43865
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-black ring, non-ridged, Canada	454209CA	43865
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-black ring, non-ridged, China	454209CN	43865
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-white ring, PREMIUM	454210	43865
K2E	K2EDTA	VACUETTE® TUBE 3 ml K2E K2EDTA 13x75 lavender cap-red ring, non-ridged	454234	43865
K2E	K2EDTA	VACUETTE® TUBE 3 ml K2E K2EDTA 13x75 lavender cap-black ring, non-ridged	454246	43865
K2E	K2EDTA	VACUETTE® TUBE 3 ml K2E K2EDTA 13x75 lavender cap-black ring, special label, non-ridged	454246CA	43865
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-white ring, non-ridged	454286	43865
K2E	K2EDTA	VACUETTE® TUBE 3 ml K2E K2EDTA 13x75 purple cap-black ring, non-ridged	454346	43865
K2E	K2EDTA	VACUETTE® TUBE 2 ml K2E K2EDTA 13x75 lavender cap-white ring, transparent label, non-ridged	454409	43865
K2E	K2EDTA	VACUETTE® TUBE 2 ml K2E K2EDTA 13x75 lavender cap-white ring, transparent special label, non-ridged, Japan	454409JP	43865
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-black ring, transparent label, non-ridged	454410	43865
K2E	K2EDTA	VACUETTE® TUBE 3 ml K2E K2EDTA 13x75 lavender cap-black ring, transparent label, non-ridged	454411	43865
K2E	K2EDTA	VACUETTE® TUBE 2 ml K2E K2EDTA 13x75 lavender cap-white ring, non-ridged, USA	454428	43865
K2E	K2EDTA	VACUETTE® TUBE 3 ml K2E K2EDTA 13x75 gold cap-black ring, non-ridged	454503	43865
K2E	K2EDTA	VACUETTE® TUBE 3 ml K2E K2EDTA 13x75 yellow cap-black ring, non-ridged	454504	43865
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 pink cap-black ring, transparent label, non-ridged	454505	43865

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
K2E	K2EDTA	VACUETTE® TUBE 2 ml K2E K2EDTA 13x75 yellow cap-white ring, non-ridged	454506	43865
K2E	K2EDTA	VACUETTE® TUBE 9 ml K2E K2EDTA 16x100 lavender cap-black ring, non-ridged	455045	43865
K2E	K2EDTA	VACUETTE® TUBE 10 ml K2E K2EDTA 16x100 lavender cap-black ring, non-ridged	455046	43865
K2E	K2EDTA	VACUETTE® TUBE 6 ml K2E K2EDTA 13x100 lavender cap-black ring, non-ridged, USA	456002	43865
K2E	K2EDTA	VACUETTE® TUBE 6 ml K2E K2EDTA 13x100 lavender cap-black ring, special label, non-ridged	456002CA	43865
K2E	K2EDTA	VACUETTE® TUBE 6 ml K2E K2EDTA 13x100 lavender cap-black ring, PREMIUM	456023	43865
K2E	K2EDTA	VACUETTE® TUBE 6 ml K2E K2EDTA 13x100 lavender cap-black ring, non-ridged	456043	43865
K2E	K2EDTA	VACUETTE® TUBE 6 ml K2E K2EDTA 13x100 pink cap-black ring, transparent special label, non-ridged	456068	43865
K2E	K2EDTA	VACUETTE® TUBE 6 ml K2E K2EDTA 13x100 lavender cap-black ring, transparent label, non-ridged	456243	43865
K2E	K2EDTA	VACUETTE® TUBE 6 ml K2E K2EDTA 13x100 red cap-black ring, PREMIUM	456244	43865
K2E	K2EDTA	VACUETTE® TUBE 6 ml K2E K2EDTA 13x100 pink cap-pink ring, non-ridged	456279	43865
K2E	K2EDTA	VACUETTE® TUBE 6 ml K2E K2EDTA 13x100 pink cap-pink ring, transparent label, non-ridged	456280	43865
K2E	K2EDTA	VACUETTE® TUBE 3 ml K2E K2EDTA 13x75 lavender cap-black ring, barcode label, non-ridged	484252	43865
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 pink cap-black ring, barcode label, non-ridged	484506	43865
K2E	K2EDTA	VACUETTE® TUBE 4 ml K2E K2EDTA 13x75 lavender cap-black ring, barcode label, non-ridged	484507	43865
K2E	K2EDTA	VACUETTE® TUBE 3 ml K2E K2EDTA 13x75 lavender cap-black ring, PREMIUM, G-Barcode	484526	43865
K2E	K2EDTA	VACUETTE® TUBE 9 ml K2E K2EDTA 16x100 lavender cap-black ring, barcode label, non-ridged	485021	43865
K2E	K2EDTA	VACUETTE® TUBE 9 ml K2E K2EDTA 16x100 lavender cap-black ring, Bi-Barcode label, non-ridged	495045	43865
K2E	K2EDTA Sep	VACUETTE® TUBE 3.5 ml K2E K2EDTA Separator 13x75 lavender cap-yellow ring, non-ridged	454059	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 3.5 ml K2E K2EDTA Separator 13x75 lavender cap-yellow ring, PREMIUM	454235	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 3.5 ml K2E K2EDTA Separator 13x75 lavender cap-yellow ring, PREMIUM	454235P	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 3.5 ml K2E K2EDTA Separator 13x75 white cap-yellow ring, non-ridged	454509P	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 8 ml K2E K2EDTA Separator 16x100 white cap-yellow ring, non-ridged	455021	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 8 ml K2E K2EDTA Separator 16x100 white cap-yellow ring, non-ridged, USA	455021P	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 8 ml K2E K2EDTA Separator 16x100 lavender cap-yellow ring, non-ridged	455040	47587

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
K2E	K2EDTA	VACUETTE® TUBE 8 ml K2E K2EDTA Separator 16x100 lavender cap-yellow ring, non-ridged, USA	455040P	43865
K2E	K2EDTA Sep	VACUETTE® TUBE 8 ml K2E K2EDTA Separator 16x100 lavender cap-yellow ring, non-ridged, China	455040CN	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 lavender cap-yellow ring, PREMIUM	456011	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 lavender cap-yellow ring, PREMIUM, China	456011CN	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 4.5 ml K2E K2EDTA Separator 13x100 lavender cap-yellow ring, non-ridged, USA	456024P	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 white cap-yellow ring, non-ridged	456058	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 white cap-yellow ring, non-ridged	456058P	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 white cap-yellow ring, special label, non-ridged	456058CN	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 white cap-yellow ring, transparent label, non-ridged	456205	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 white cap-yellow ring, transparent label, non-ridged	456205P	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 white cap-yellow ring, PREMIUM	456211	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 white cap-yellow ring, PREMIUM	456211P	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 white cap-yellow ring, special label, non-ridged	456225	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 lavender cap-yellow ring, non-ridged	456282	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 lavender cap-yellow ring, non-ridged	456282P	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 lavender cap-yellow ring, special label, non-ridged, Japan	456282JP	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 pink cap-yellow ring, non-ridged	456350	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 pink cap-yellow ring, non-ridged	456350P	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 white cap-yellow ring, barcode label, PREMIUM	486506	47587
K2E	K2EDTA Sep	VACUETTE® TUBE 5 ml K2E K2EDTA Separator 13x100 white cap-yellow ring, barcode label, PREMIUM	486506P	47587
K3E	Crossmatch K3EDTA	VACUETTE® TUBE 6 ml K3E Crossmatch K3EDTA 13x100 pink cap-black ring, special label, non-ridged	456052	47588
K3E	Crossmatch K3EDTA	VACUETTE® TUBE 6 ml K3E Crossmatch K3EDTA 13x100 pink cap-black ring, special label, PREMIUM	456061	47588
K3E	Crossmatch K3EDTA	VACUETTE® TUBE 6 ml K3E Crossmatch K3EDTA 13x100 pink cap-black ring, special label, PREMIUM	456093	47588
K3E	Crossmatch K3EDTA	VACUETTE® TUBE 6 ml K3E Crossmatch K3EDTA 13x100 pink cap-black ring, special label, PREMIUM, China	456093CN	47588
K3E	Crossmatch K3EDTA	VACUETTE® TUBE 6 ml K3E Crossmatch K3EDTA 13x100 pink cap-black ring, special label, PREMIUM, Japan	456093JP	47588

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
K3E	Crossmatch K3EDTA	VACUETTE® TUBE 6 ml K3E Crossmatch K3EDTA 13x100 pink cap-black ring, special label, non-ridged	456242	47588
K3E	Crossmatch K3EDTA	VACUETTE® TUBE 6 ml K3E Crossmatch K3EDTA 13x100 pink cap-black ring, special label, non-ridged	456252	47588
K3E	Crossmatch K3EDTA	VACUETTE® TUBE 6 ml K3E Crossmatch K3EDTA 13x100 yellow cap-black ring, special label, PREMIUM	456255	47588
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 lavender cap-red ring, PREMIUM	454003	47588
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 pink cap-black ring, non-ridged	454012	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-red ring, non-ridged	454015	47588
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 lavender cap-black ring, non-ridged	454021	47588
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 lavender cap-black ring, non-ridged, High Altitude	454021H	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 red cap-black ring, non-ridged	454022	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 orange cap-black ring, non-ridged	454025	47588
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 lavender cap-black ring, special label, non-ridged	454026	47588
K3E	K3EDTA	VACUETTE® TUBE 1 ml K3E K3EDTA 13x75 lavender cap-white ring, PREMIUM	454034	47588
K3E	K3EDTA	VACUETTE® TUBE 2.5 ml K3E K3EDTA Liquid 13x75 lavender cap-black ring, PREMIUM	454035	47588
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 lavender cap-black ring, PREMIUM	454036	47588
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 lavender cap-black ring, PREMIUM, China	454036CN	47588
K3E	K3EDTA	VACUETTE® TUBE 4.5 ml K3E K3EDTA 13x75 lavender cap-black ring, PREMIUM	454039	47588
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 yellow cap-black ring, non-ridged	454040	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA Liquid 13x75 lavender cap-black ring, PREMIUM	454041	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA Liquid 13x75 lavender cap-red ring, non-ridged	454042	47588
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 lavender cap-black ring, PREMIUM, 100/Box	454048	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, PREMIUM	454086	47588
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, PREMIUM	454087	47588
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, PREMIUM, China	454087CN	47588
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 ping cap-black ring, pink label, PREMIUM	454099	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, non-ridged	454217	47588

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, non-ridged, High Altitude	454217H	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, non-ridged, Korea	454217KR	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 black cap-black ring, non-ridged	454220	47588
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, non-ridged	454222	47588
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, non-ridged, High Altitude	454222H	47588
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, non-ridged, China	454222CN	47588
K3E	K3EDTA	VACUETTE® TUBE 4.5 ml K3E K3EDTA 13x75 lavender cap-black ring, non-ridged	454223	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 royal blue cap-black ring, non-ridged	454227	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-white ring, PREMIUM	454232	47588
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 black cap-black ring, PREMIUM	454245	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, transparent label, non-ridged	454287	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, transparent label, non-ridged, Korea	454287KR	47588
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, transparent label, non-ridged	454317	47588
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, transparent special label, non-ridged	454317CN	47588
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 lavender cap-black ring, non-ridged 100/Box	454408	47588
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 lavender cap-black ring, non-ridged, transparent label	454425	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, special label, non-ridged	454452	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-red ring, transparent label, non-ridged	454515	47588
K3E	K3EDTA	VACUETTE® TUBE 9 ml K3E K3EDTA 16x100 black cap-black ring, black label, non-ridged	455030	47588
K3E	K3EDTA	VACUETTE® TUBE 9 ml K3E K3EDTA 16x100 lavender cap-black ring, non-ridged	455036	47588
K3E	K3EDTA	VACUETTE® TUBE 9 ml K3E K3EDTA 16x100 lavender cap-black ring, non-ridged, China	455036CN	47588
K3E	K3EDTA	VACUETTE® TUBE 9 ml K3E K3EDTA 16x100 pink cap-black ring, pink label, non-ridged	455038	47588
K3E	K3EDTA	VACUETTE® TUBE 9 ml K3E K3EDTA 16x100 gold cap-gold ring, non-ridged	455235	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 pink cap-black ring, non-ridged	456003	47588

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 pink cap-black ring, special label, non-ridged	456003CA	47588
K3E	K3EDTA	VACUETTE® TUBE 5 ml K3E K3EDTA 13x100 lavender cap-black ring, PREMIUM	456034	47588
K3E	K3EDTA	VACUETTE® TUBE 5 ml K3E K3EDTA 13x100 lavender cap-black ring, PREMIUM, China	456034CN	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 lavender cap-black ring, PREMIUM	456036	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA Liquid 13x100 lavender cap-black ring, PREMIUM	456037	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 lavender cap-black ring, non-ridged	456038	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 lavender cap-black ring, special label, non-ridged, SGH	456038SGH	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 yellow cap-black ring, PREMIUM	456039	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 royal blue cap-black ring, non-ridged	456042	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 lavender cap-black ring, transparent label, non-ridged	456053	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 pink cap-black ring, transparent label, non-ridged	456067	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 white cap-black ring, non-ridged	456074	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 pink cap-black ring, pink label, PREMIUM	456099	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 lavender cap-red ring, non-ridged	456209	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 red cap-black ring, non-ridged	456235	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 green cap-black ring, non-ridged	456236	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 yellow cap-black ring, non-ridged	456238	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 lavender cap-red ring, special label, non-ridged	456251	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 lavender cap-black ring, special label, PREMIUM	456256	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 pink cap-red ring, special label, non-ridged	456257	47588
K3E	K3EDTA	VACUETTE® TUBE 5 ml K3E K3EDTA Liquid 13x100 pink cap-black ring, special label, non-ridged	456258	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 pink cap-yellow ring, special label, non-ridged	456259	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 yellow cap-red ring, special label, non-ridged	456260	47588
K3E	K3EDTA	VACUETTE® TUBE 5 ml K3E K3EDTA Liquid 13x100 yellow cap-black ring, special label, non-ridged	456261	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 yellow cap-yellow ring, special label, non-ridged	456262	47588

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 red cap-red ring, special label, non-ridged	456263	47588
K3E	K3EDTA	VACUETTE® TUBE 5 ml K3E K3EDTA Liquid 13x100 red cap-black ring, special label, non-ridged	456265	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 red cap-yellow ring, special label, non-ridged	456266	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 lavender cap-black ring, special label, non-ridged	456276	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, barcode label, non-ridged	464217U	47588
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 lavender cap-black ring, transparent label, PREMIUM	474036	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, transparent label, PREMIUM	474086	47588
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, transparent label, PREMIUM	474087	47588
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, transparent label, PREMIUM, China	474087CN	47588
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, transparent label, PREMIUM, Japan	474087JP	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-white ring, transparent label, PREMIUM	474232	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 orange cap-black ring, barcode label, non-ridged	484025	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 red cap-black ring, barcode label, non-ridged	484026	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 pink cap-black ring, barcode label, non-ridged	484027	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, barcode label, non-ridged	484217	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-red ring, barcode label, non-ridged	484218	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 gold cap-black ring, barcode label, non-ridged	484404	47588
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, barcode label, non-ridged	484481	47588
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, barcode label, non-ridged	484500	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, barcode label, non-ridged	484501	47588
K3E	K3EDTA	VACUETTE® TUBE 1 ml K3E K3EDTA 13x75 lavender cap-white ring, barcode label, G-Code, PREMIUM	484509	47588
K3E	K3EDTA	VACUETTE® TUBE 4 ml K3E K3EDTA 13x75 lavender cap-black ring, barcode label, G-Code, PREMIUM	484510	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, barcode label, G-Code, PREMIUM	484511	47588

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
K3E	K3EDTA	VACUETTE® TUBE 2 ml K3E K3EDTA 13x75 lavender cap-white ring, barcode label, G-Code, PREMIUM	484512	47588
K3E	K3EDTA	VACUETTE® TUBE 9 ml K3E K3EDTA 16x100 lavender cap-black ring, barcode label, G-Code, non- ridged	485501	47588
K3E	K3EDTA	VACUETTE® TUBE 9 ml K3E K3EDTA 16x100 orange cap-black ring, barcode label, G-Code, non- ridged	485512	47588
K3E	K3EDTA	VACUETTE® TUBE 6 ml K3E K3EDTA 13x100 lavender cap-black ring, barcode label, G-Code, PREMIUM	486502	47588
K3E	K3EDTA	VACUETTE® TUBE 3 ml K3E K3EDTA 13x75 lavender cap-black ring, bi-barcode label, non-ridged	494217	47588
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 13x75 green cap-black ring, non-ridged	454029	47589
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 13x75 green cap-black ring, non-ridged, High Altitude	454029H	47589
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 13x75 green cap-black ring, non-ridged, China	454029CN	47589
LH	Lithium Heparin	VACUETTE® TUBE 4.5 ml LH Lithium Heparin 13x75 green cap-black ring, PREMIUM	454049	47589
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 13x75 green cap-black ring, transparent label, non-ridged	454056	47589
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 13x75 green cap-black ring, transparent special label, non- ridged, Japan	454056JP	47589
LH	Lithium Heparin	VACUETTE® TUBE 1 ml LH Lithium Heparin 13x75 green cap-white ring, PREMIUM	454081	47589
LH	Lithium Heparin	VACUETTE® TUBE 3 ml LH Lithium Heparin 13x75 green cap-black ring, PREMIUM	454082	47589
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 13x75 green cap-black ring, PREMIUM	454084	47589
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 13x75 green cap-black ring, special label, PREMIUM	454084CN	47589
LH	Lithium Heparin	VACUETTE® TUBE 2 ml LH Lithium Heparin 13x75 green cap-white ring, PREMIUM	454089	47589
LH	Lithium Heparin	VACUETTE® TUBE 2 ml LH Lithium Heparin 13x75 green cap-white ring, PREMIUM, China	454089CN	47589
LH	Lithium Heparin	VACUETTE® TUBE 3 ml LH Lithium Heparin 13x75 green cap-white ring, PREMIUM	454233	47589
LH	Lithium Heparin	VACUETTE® TUBE 2 ml LH Lithium Heparin 13x75 green cap-white ring, non-ridged	454237	47589
LH	Lithium Heparin	VACUETTE® TUBE 3 ml LH Lithium Heparin 13x75 green cap-black ring, non-ridged	454244	47589
LH	Lithium Heparin	VACUETTE® TUBE 3 ml LH Lithium Heparin 13x75 green cap-black ring, non-ridged, China	454244CN	47589
LH	Lithium Heparin	VACUETTE® TUBE 3 ml LH Lithium Heparin 13x75 green cap-white ring, non-ridged	454312	47589
LH	Lithium Heparin	VACUETTE® TUBE 3 ml LH Lithium Heparin 13x75 green cap-black ring, transparent label, non-ridged	454501	47589

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
LH	Lithium Heparin	VACUETTE® TUBE 9 ml LH Lithium Heparin 16x100 green cap-black ring, transparent label, non-ridged	455073	47589
LH	Lithium Heparin	VACUETTE® TUBE 9.5 ml LH Lithium Heparin 16x100 green cap-black ring, non-ridged	455080	47589
LH	Lithium Heparin	VACUETTE® TUBE 9 ml LH Lithium Heparin 16x100 green cap-black ring, non-ridged	455084	47589
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 16x100 green cap-black ring, non-ridged	455222	47589
LH	Lithium Heparin	VACUETTE® TUBE 6 ml LH Lithium Heparin 13x100 green cap-black ring, PREMIUM	456084	47589
LH	Lithium Heparin	VACUETTE® TUBE 6 ml LH Lithium Heparin 13x100 green cap-black ring, special label, PREMIUM, China	456084CN	47589
LH	Lithium Heparin	VACUETTE® TUBE 6 ml LH Lithium Heparin 13x100 green cap-black ring, non-ridged	456088	47589
LH	Lithium Heparin	VACUETTE® TUBE 6 ml LH Lithium Heparin 13x100 green cap-black ring, special label, non-ridged, China	456088CN	47589
LH	Lithium Heparin	VACUETTE® TUBE 6 ml LH Lithium Heparin 13x100 green cap-black ring, transparent label, non-ridged	456271	47589
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 13x75 green cap-black ring, transparent label, PREMIUM	474084	47589
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 13x75 green cap-black ring, transparent label, PREMIUM, China	474084CN	47589
LH	Lithium Heparin	VACUETTE® TUBE 6 ml LH Lithium Heparin 13x100 green cap-black ring, transparent label, PREMIUM	476084	47589
LH	Lithium Heparin	VACUETTE® TUBE 6 ml LH Lithium Heparin 13x100 green cap-black ring, transparent special label, PREMIUM, Japan	476084JP	47589
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 13x75 green cap-black ring, barcode label, non-ridged	484029	47589
LH	Lithium Heparin	VACUETTE® TUBE 4 ml LH Lithium Heparin 13x75 green cap-black ring, barcode label, G-Code, PREMIUM	484522	47589
LH	Lithium Heparin	VACUETTE® TUBE 9 ml LH Lithium Heparin 16x100 green cap-black ring, barcode label, G-Code, non-ridged	485508	47589
LH	Lithium Heparin	VACUETTE® TUBE 6 ml LH Lithium Heparin 13x100 green cap-black ring, barcode lable, G-Code, PREMIUM	486508	47589
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, non-ridged	454008	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, non-ridged	454008P	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, non-ridged, China	454008CN	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 2.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, PREMIUM	454046	46054

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
LH	Lithium Heparin Sep	VACUETTE® TUBE 2.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, PREMIUM	454046P	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-black ring, non-ridged	454057	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, PREMIUM	454083	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, non-ridged	454247	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, non-ridged	454247P	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, transparent label, non-ridged	454285	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3 ml LH Lithium Heparin Separator 13x75 gold cap-white ring, non-ridged	454446	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3 ml LH Lithium Heparin Separator 13x75 gold cap-white ring, non-ridged, USA	454446P	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, PREMIUM	454454	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, special label, non-ridged	454483	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 mint green cap-yellow ring, transparent label, non-ridged	454489	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3 ml LH Lithium Heparin Separator 13x75 mint green cap-yellow ring, non-ridged, special label	454490PCA	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, transparent label, non-ridged	454491	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, transparent label, non-ridged	454491P	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3 ml LH Lithium Heparin Separator 13x75 mint green cap-yellow ring, transparent label, non-ridged	454493	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, transparent label, non-ridged	454507	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 8 ml LH Lithium Heparin Separator 16x100 green cap-yellow ring, non-ridged	455083	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 8 ml LH Lithium Heparin Separator 16x100 green cap-yellow ring, non-ridged, USA	455083P	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 4 ml LH Lithium Heparin Separator 16x100 green cap-yellow ring, non-ridged	455229	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, PREMIUM	456083	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, PREMIUM	456083P	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, PREMIUM, China	456083CN	46054

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, non-ridged	456087	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, non-ridged	456087P	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, special label, non-ridged, USA	456087CA	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, non-ridged, Korea	456087KR	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, non-ridged, Canada	456087PCA	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, transparent label, non-ridged	456216	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, transparent special label, non-ridged, Japan	456216JP	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, transparent label, non-ridged, Korea	456216KR	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, non-ridged	456287	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, non-ridged, USA	456287P	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 mint green cap-yellow ring, special label, non-ridged	456289CA	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 mint green cap-yellow ring, special label, non-ridged	456289PCA	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-white ring, non-ridged	456303	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 mint green cap-yellow ring, non-ridged	456305	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, transparent label, PREMIUM	474080	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, transparent label, PREMIUM	474083	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, transparent label, PREMIUM	476083	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, transparent label, PREMIUM, Japan	476083JP	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-red ring, barcode label, non-ridged	484030	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, barcode label, non-ridged	484031	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 3.5 ml LH Lithium Heparin Separator 13x75 green cap-yellow ring, barcode label, G-Code, non-ridged	484523	46054

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
LH	Lithium Heparin Sep	VACUETTE® TUBE 8 ml LH Lithium Heparin Separator 16x100 green cap-yellow ring, barcode label, non-ridged	485500	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 5 ml LH Lithium Heparin Separator 13x100 green cap-yellow ring, barcode label, G-Code, non-ridged	486507	46054
LH	Lithium Heparin Sep	VACUETTE® TUBE 8 ml LH Lithium Heparin Separator 16x100 green cap-yellow ring, Bi-Barcode label, non-ridged	495083	46054
LH	Lithium Heparin/	VACUETTE® TUBE 4.5 ml Lithium Heparin / Iodoacetate 13x75 grey cap-black ring, PREMIUM	454080	45343
LH	Lithium Heparin/	VACUETTE® TUBE 4.5 ml Lithium Heparin / Iodoacetate 13x75 grey cap-black ring, PREMIUM	454248	45343
NH	Sodium Heparin	VACUETTE® TUBE 4 ml NH Sodium Heparin 13x75 green cap-green ring, non-ridged	454030	47592
NH	Sodium Heparin	VACUETTE® TUBE 4 ml NH Sodium Heparin 13x75 green cap-green ring, non-ridged, High Altitude	454030H	47592
NH	Sodium Heparin	VACUETTE® TUBE 4 ml NH Sodium Heparin 13x75 green cap-green ring, PREMIUM	454051	47592
NH	Sodium Heparin	VACUETTE® TUBE 4 ml NH Sodium Heparin 13x75 green cap-green ring, special label, PREMIUM	454051CN	47592
NH	Sodium Heparin	VACUETTE® TUBE 2 ml NH Sodium Heparin 13x75 green cap-green ring, non-ridged	454302	47592
NH	Sodium Heparin	VACUETTE® TUBE 4 ml NH Sodium Heparin 13x75 royal blue cap-black ring, transparent label, non-ridged	454502	47592
NH	Sodium Heparin	VACUETTE® TUBE 9 ml NH Sodium Heparin 16x100 green cap-green ring, non-ridged	455051	47592
NH	Sodium Heparin	VACUETTE® TUBE 9 ml NH Sodium Heparin 16x100 green cap-green ring, transparent label, non-ridged	455332	47592
NH	Sodium Heparin	VACUETTE® TUBE 6 ml NH Sodium Heparin 13x100 green cap-green ring, non-ridged	456028	47592
NH	Sodium Heparin	VACUETTE® TUBE 6 ml NH Sodium Heparin 13x100 green cap-green ring, special label, non-ridged, Japan	456028JP	47592
NH	Sodium Heparin	VACUETTE® TUBE 6 ml NH Sodium Heparin 13x100 green cap-green ring, PREMIUM	456051	47592
NH	Sodium Heparin	VACUETTE® TUBE 9 ml NH Sodium Heparin 16x100 green cap-black ring, barcode label, G-Code, non-ridged	485505	47592
NH	Trace Elements Sodium Heparin	VACUETTE® TUBE 6ml NH Trace Elements Sodium Heparin 13x100 royal blue cap-black ring, PREMIUM	456080	47592
NH	Trace Elements Sodium Heparin	VACUETTE® TUBE 6 ml NH Trace Elements Sodium Heparin 13x100 royalblue cap-black ring, non-ridged	456275	47592
Others	ACD-A	VACUETTE® TUBE 9 ml ACD-A 16x100 yellow cap-black ring, transparent label, non-ridged	455032	42583
Others	ACD-A	VACUETTE® TUBE 9 ml ACD-A 16x100 yellow cap-black ring, non-ridged	455055	42583
Others	ACD-A	VACUETTE® TUBE 6 ml ACD-A 13x100 yellow cap-black ring, PREMIUM	456055	42583
Others	ACD-B	VACUETTE® TUBE 9 ml ACD-B 16x100 yellow cap-black ring, non-ridged	455094	42370
Others	ACD-B	VACUETTE® TUBE 6 ml ACD-B 13x100 yellow cap-black ring, PREMIUM	456094	42370

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
Others	ACD-B	VACUETTE® TUBE 6 ml ACD-B 13x100 yellow cap-black ring, non-ridged	456204	42370
Others	CPDA	VACUETTE® TUBE 4 ml CPDA 13x75 yellow cap-black ring, non-ridged	454280	42914
Others	CPDA	VACUETTE® TUBE 9 ml CPDA 16x100 yellow cap-black ring, non-ridged	455056	42914
Others	CPDA	VACUETTE® TUBE 6 ml CPDA 13x100 yellow cap-black ring, PREMIUM	456056	42914
Others	CPDA	VACUETTE® TUBE 6 ml CPDA 13x100 yellow cap-black ring, non-ridged	456057	42914
Others	CTAD	VACUETTE® TUBE 3.5 ml CTAD 13x75 blue cap-yellow ring, sandwich tube, PREMIUM	454064	42569
Others	CTAD	VACUETTE® TUBE 3 ml CTAD 13x75 blue cap-yellow ring, sandwich tube, transparent label, non-ridged	454441	42569
Others	CTAD	VACUETTE® TUBE 3.5 ml CTAD 13x75 blue cap-yellow ring, sandwich tube, transparent label, non-ridged	454462	42569
Others	CTAD	VACUETTE® TUBE 3.5 ml CTAD 13x75 blue cap-yellow ring, sandwich tube, transparent label, PREMIUM	474064	42569
Others	CTAD	VACUETTE® TUBE 2 ml CTAD 13x75 blue cap-yellow ring, sandwich tube, transparent label, PREMIUM	474065	42569
Others	CTAD	VACUETTE® TUBE 3 ml CTAD 13x75 blue cap-yellow ring, sandwich tube, transparent label, PREMIUM	474304	42569
Others	EDTA Aprotinin	VACUETTE® TUBE 4 ml EDTA Aprotinin 13x75 pink cap-black ring, PREMIUM	454261	58135
Others	FC-Mix	VACUETTE® FC Mix TUBE 2 ml 13x75 pink cap-white ring, transparent label, non-ridged	454510	44208
Others	FC-Mix	VACUETTE® FC Mix TUBE 2 ml 13x75 grey cap-white ring, transparent label, non-ridged	454511	44208
Others	FC-Mix	VACUETTE® FC Mix TUBE 3 ml 13x75 pink cap-black ring, transparent label, non-ridged	454513	44208
Others	FC-Mix	VACUETTE® FC Mix TUBE 3 ml 13x75 grey cap-black ring, transparent label, non-ridged	454514	44208
Others	Homocysteine Detection Tube	VACUETTE® TUBE 2 ml Homocysteine Detection 13x75 white cap-red ring, sandwich tube, PREMIUM	454421	42585
Z	Discard Tube	VACUETTE® TUBE 5 ml Z Discard Tube 13x100 black cap-red ring, non-ridged	456202	47590
Z	Discard Tube	VACUETTE® TUBE 5 ml Z Discard Tube 13x100 white cap-black ring, non-ridged	456217	47590
Z	No Additive	VACUETTE® TUBE 4 ml Z No Additive 13x75 white cap-black ring, PREMIUM	454001	47590
Z	No Additive	VACUETTE® TUBE 4 ml Z No Additive 13x75 white cap-black ring, PREMIUM China	454001CN	47590
Z	No Additive	VACUETTE® TUBE 3.5 ml Z No Additive 13x75 white cap-black ring, non-ridged	454045	47590
Z	No Additive	VACUETTE® TUBE 4 ml Z Silicone coated 13x75 red cap-whiter ring, PREMIUM	454072	47590

VACUETTE® Tube Group	Additive	Product name - detailed product description	Item No.	GMDN
Z	No Additive	VACUETTE® TUBE 2 ml Z No Additive 13x75 white cap-white ring, PREMIUM	454088	47590
Z	No Additive	VACUETTE® TUBE 3 ml Z No Additive 13x75 red cap-white ring, transparent label, non-ridged	454224	47590
Z	No Additive	VACUETTE® TUBE 2 ml Z No Additive 13x75 blue cap-black ring, non-ridged	454240	47590
Z	No Additive	VACUETTE® TUBE 3 ml Z No Additive 13x75 white cap-black ring, non-ridged	454241	47590
Z	No Additive	VACUETTE® TUBE 3 ml Z No Additive 13x75 white cap-yellow ring, non-ridged, High Altitude	454241H	47590
Z	No Additive	VACUETTE® TUBE 3 ml Z No Additive 13x75 gold cap-gold ring, red label, non-ridged	454256	47590
Z	No Additive	VACUETTE® TUBE 2 ml Z No Additive 13x75 black cap-black ring, non-ridged	454318	47590
Z	No Additive	VACUETTE® TUBE 2 ml Z No Additive 13x75 black cap-yellow ring, non-ridged, High Altitude	454318H	47590
Z	No Additive	VACUETTE® TUBE 2 ml Z No Additive 13x75 lavender cap-black ring, non-ridged	454370	47590
Z	No Additive	VACUETTE® TUBE 2 ml Z No Additive 13x75 lavender cap-yellow ring, non-ridged, High Altitude	454370H	47590
Z	No Additive	VACUETTE® TUBE 4 ml Z No Additive 13x75 red cap-black ring, non-ridged	454418	47590
Z	No Additive	VACUETTE® TUBE 4 ml Z No Additive 13x75 white cap-black ring, non-ridged	454484	47590
Z	No Additive	VACUETTE® TUBE 9 ml Z No Additive 16x100 white cap-black ring, non-ridged	455001	47590
Z	No Additive	VACUETTE® TUBE 9 ml Z No Additive 16x100 white cap-black ring, non-ridged, Japan	455001JP	47590
Z	No Additive	VACUETTE® TUBE 9 ml Z No Additive 16x100 red cap-black ring, non-ridged	455331	47590
Z	No Additive	VACUETTE® TUBE 6 ml Z No Additive 13x100 white cap-black ring, PREMIUM	456001	47590
Z	No Additive	VACUETTE® TUBE 6 ml Z No Additive 13x100 white cap-black ring, non-ridged	456085	47590
Z	No Additive	VACUETTE® TUBE 6 ml Z No Additive 13x100 white cap-black ring, non-ridged	456240	47590
Z	No Additive	VACUETTE® TUBE 6.2 ml Z No Additive 13x100 white cap-black ring	456283	47590
Z	No Additive	VACUETTE® TUBE 5 ml Z No Additive 13x100 white cap-black ring	456285	47590
Z	No Additive	VACUETTE® TUBE 2 ml Z No Additive 13x75 white cap-white ring, Barcode label, PREMIUM	484088	47590
Z	No Additive	VACUETTE® TUBE 9 ml Z No Additive 16x100 white cap-black ring, barcode label, G-Code, non-ridged	485506	47590
Z	Trace Elements No Additive	VACUETTE® TUBE 6 ml Z Trace Elements No Additive 13x100 royal blue cap-black ring, non-ridged	456016	47590