

Standard Analysis Certificate

PRODUCT :DICHLOROMETANE STAB. ETHANOL RS SPECTROSOL
CODE :463025
METHOD :20066

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Colour (APHA)	-	<=10
Identification (I.R.)	-	Positive
Density at 20° C	-	1.322 ÷ 1.328
Refractive index at 20°C	-	1.4214 ÷ 1.4274
Boiling point	° C	39.6 ÷ 40.1
Water (K.F.)	ppm	<=100
Residue on evaporation	ppm	<=5
Acidity	meq/g	<=0.0005
Alcalinity	meq/g	<=0.0002
Assay (GLC)	%	>=99.9
Fluorescence	-	-
at 254 nm	ppb	<=2
at 365 nm	ppb	<=2
U.V. Transmittance	-	-
at 235 nm	%	>=35
at 240 nm	%	>=70
at 245 nm	%	>=90
at 250 nm	%	>=95
at 260 nm	%	>=99
Stabilizer (Ethanol)	%	0.1 ÷ 0.4

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

B. COULANGE (VDR)