

10.005
10.008**Shoecovers with 80 µm anti-slip strong sole****FEATURES & BENEFITS**

- Shoecovers with reinforced polyethylene sole (80 µm).
- Ideal for use in clean environments thanks to their low particle shedding – Tested using the Helmke drum method.
- Embossed soles provide good floor grip and extended durability against abrasion on uneven surfaces.
- They fit all types of footwear.

Suggested applications

Pharmaceutical industries / Research laboratories / Nuclear, electronics, and aerospace industries / Food industry

PRODUCT DETAILS

Product name :	Kolmi [®] Isoguard
Product type :	À usage unique, non stérile
Material :	Spunbond polypropylene Embossed polyethylene film
Spunbond grammage :	30 g/m ²
Sole thickness :	80 µm
Confection :	Gluing
Size :	165 x 400 mm (XL Size)
Unit weight:	12 g
Origin :	France

CERTIFICATIONS & NORMES

Manufacturing site : ISO 9001 and ISO 13485 certifications
Helmke drum tested according to IEST-RP-CC003.4

MICROBIOLOGICAL INFORMATION

Assessment of microbial contamination according to ISO 11737:2018+A1:2021 available upon request.
Additional microbiological tests (ASR, E. coli, staphylococci) available upon request.

**10.005**
WHITE
SOLE**10.008**
BLUE/GREY
SOLE

10.005
10.008

Shoecovers with 80 µm anti-slip strong sole

STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

LOGISTICS



**Outer case
with inside plastic bag**

Item	PCS qty per case	Size (mm)	Gross weight (kg)	QTY/pallet
10.005	200	300 x 400 x 400	3,4	30 (5 layers of 6 cases)
10.008	400	380 x 400 x 400	5,3	24 (4 layers of 6 cases)

EAN codes

Item	Outer case
10.005	3 700 381 500 157
10.008	3 700 381 500 188