

20.811

WHITE UNDERWEAR KIT 100% PP 1 TEE SHIRT, 1 BOXER SHORT, 1 PAIR OF SOCKS

FEATURES & BENEFITS

- Used in environments imposing airlocks and decontamination showers.
- Pleasant and light, polypropylene is very breathable.
- Tee-shirt with large form allowing the user to have a good comfort.
- One size fits all.
- Product destroyed after use to avoid any washing and contamination risk.

TECHNICAL DESCRIPTION

Product name: White Underwear Kit 100% PP
 Product type: Single use, non sterile
 Composition: PP
 Grammage: 30 g/m²
 Country of origin: China



Composition	Unit weight (g)	Length (cm)	Width (cm)
1 t-shirt	25 g	78	130
1 boxer short	12 g	115	41,5
1 pair of socks	13 g	57	31

INDUSTRIES

Industrial / Restaurant & catering / Community services / Hygiene

STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

LOGISTIC INFORMATION

OUTER CASE SPECIFICATIONS					SPECIFICATIONS CDT
Article	SIZE (mm)	WEIGHT (kg)	Qty / pallet	EAN CODE	CDT
20.811	450 x 280 x 450	4 – 5 kg	24 boxes		1 kit per polybag 100 polybags per box