

**40.301**  
**40.301R**  
**40.301V**

## CLASSIC HYGIENE PP GOWN 40g/m<sup>2</sup> COLLAR & SNAP FASTENER

### FEATURES & BENEFITS

- Single use hygiene gowns protect the working environment from particles coming from the body of the operator (dander, hair, sweating ...)
- They allow the operator to protect himself from the risks of soiling related to the activity he carries out.

### TECHNICAL DESCRIPTION

Product name: Classic Hygiene PP Gown  
Product type: Single use, non sterile  
Material: Polypropylene  
Confection: Sewed & elasticated (collar, wrists)  
Snap fastener  
Color: White / Green / Red  
Grammage: 40 g/m<sup>2</sup> (± 10%)  
Unit weight (Size M): 108 g  
Packaging: 50 individual polybags  
Origin: China



| Color                      | M         | L         | XL         | XXL        | XXXL       |
|----------------------------|-----------|-----------|------------|------------|------------|
| Blanc                      | 40.301/02 | 40.301/03 | 40.301/04  | 40.301/05  | 40.301/06  |
| Rouge                      | x         | x         | 40.301R/04 | 40.301R/05 | 40.301R/06 |
| Vert                       | x         | x         | 40.301V/04 | 40.301V/05 | x          |
| <b>Length<br/>(± 2 cm)</b> | 104 cm    | 104 cm    | 104 cm     | 104 cm     | 104 cm     |
| <b>Width<br/>(± 2 cm)</b>  | 130 cm    | 135 cm    | 140 cm     | 145 cm     | 150 cm     |

**40.301**  
**40.301R**  
**40.301V**

## CLASSIC HYGIENE PP GOWN 40g/m<sup>2</sup> COLLAR & SNAP FASTENER

### INDUSTRIES

Industrial / Restaurant & catering / Community services /  
Hygiene

### CERTIFICATIONS & STANDARDS

Manufacturing site : ISO 9001.

### STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

### LOGISTIC INFORMATIONS

| OUTERCASE SPECIFICATIONS |                 |             |               |                   | CDT SPECIFICATIONS            |                   |
|--------------------------|-----------------|-------------|---------------|-------------------|-------------------------------|-------------------|
| Article                  | SIZE (mm)       | WEIGHT (kg) | Qty / pallet  | EAN CODE          | CDT                           | EAN CODE          |
| 40.301/02                | 470 x 280 x 280 | 4,07        | 36 outercases | 3 700 381 505 510 | 1 x 50<br>individual polybags | 3 700 381 505 527 |
| 40.301/03                |                 | 4,5         |               | 3 700 381 500 799 |                               | 3 700 381 502 243 |
| 40.301/04                |                 | 4,9         |               | 3 700 381 500 805 |                               | 3 700 381 502 250 |
| 40.301/05                |                 | 5,3         |               | 3 700 381 500 812 |                               | 3 700 381 502 267 |
| 40.301/06                |                 | 5,8         |               | 3 700 381 506 999 |                               | 3 700 381 507 002 |
| 40.301R/04               |                 | 4,5         |               | 3 700 381 506 548 |                               | 3 700 381 506 562 |
| 40.301R/05               |                 | 4,4         |               | 3 700 381 506 555 |                               | 3 700 381 506 579 |
| 40.301R/06               |                 | 4,9         |               | 3 700 381 508 313 |                               | 3 700 381 508 320 |
| 40.301V/04               |                 | 4,4         |               | 3 700 381 505 558 |                               | 3 700 381 505 565 |
| 40.301V/05               |                 | 4,8         |               | 3 700 381 505 572 |                               | 3 700 381 505 589 |