

### Shoecovers with antislip sole

**10.007/XXL**

XXL size

Medical Device Class I

#### FEATURES & BENEFITS

- Hydrophobic spunbond material providing effective moisture protection.
- PE sole ensuring high mechanical resistance and strong ground adhesion.
- Anatomical shape designed to fit most types of footwear.
- Large size and extended height offering optimal compatibility with safety shoes.

#### Recommended uses

- In industrial environments on dry or wet floors
- In medical environments



#### PRODUCT DETAILS

Product Name:	Kolmi <sup>®</sup> Skidguard
Product type :	Single-use, non-sterile
Composition :	Spunbond polypropylene Embossed polyethylene film
Spunbond grammage :	35 g/m <sup>2</sup>
Sole thickness:	60 µm
Manufacturing :	Gluing
Size :	230 x 420 mm – XXL + very tall height
Unit weight :	14,4 g
Packaging :	Case of 300 pcs with protective bag
Origin :	France

Item	Color
10.007/XXL	<input type="checkbox"/> White

#### CERTIFICATIONS & STANDARDS

Meet requirements of Regulation (EU) 2017/745 on Medical Devices.  
Manufacturing site: ISO 9001 and ISO 13485 certified.

#### MICROBIOLOGICAL INFORMATION

Evaluation of initial microbial contamination according to standards ISO 11737:2018 + A1:2021 upon request.  
Additional microbiological tests: ASR, E. coli, staphylococci, available upon request.

**PRECAUTIONS**

The product must be used on intact skin only.

The reuse or extended use may lead to infection or cross-contamination.

Please refer to your national legislation on the subject of disposal of this product once it has been used.

**NORMAL STORAGE CONDITIONS**

Must not be exposed to moisture or sunlight, and should be stored in a cool or room temperature environment.

Product shelf life: 5 years.

**LOGISTICS**
**Outer case  
spec (with protective bag)**

Item	Size mm	Gross weight (kg)	QTY/ pallet
10.007/XXL	400 x 400 x 375	5	24 (4 layers of 6 cases)

 **EAN Code**

Item	Case
10.007/XXL	3 700 381 503 264