

**70.702**  
**70.703**

## PP ELASTICATED HOOD « TYPE HEAUME » WITH INTEGRATED MASK

### FEATURES & BENEFITS

- Covering ultra-effective protection of the product but also the user
- Integrated mask
- 1 item code = 2 Products
- Ability to color code tracking.
- Can also be used as an environmental control tool in HACCP and other industries.



### TECHNICAL DESCRIPTION

Product name : White PP elasticated hood « type heaume » with 2 pleats integrated PP mask

Product Type : Single use, non-sterile

Material : Polypropylene (PP)

Confection : Sewed and elasticated



Grammage : 25 g/m<sup>2</sup> (± 10%)

Dimensions : 18 x 37 cm (± 1 cm)

Unit weight : 8,2 g (± 10%)

Packaging : Polybags of 50 units

Origin : China

Item	Color
70.702	Blue 
70.703	White 

### INDUSTRIES

Industry / Catering / Food / Hygiene

### CERTIFICATIONS & STANDARDS

Manufacturing site : ISO 9001

### STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

### LOGISTIC INFORMATIONS

SPECIFICATIONS OUTERCASE					SPECIFICATIONS CDT	
Article	SIZE (mm)	WEIGHT (kg)	Qty / pallet	CODE EAN	CDT	CODE EAN
70.702	430 x 300 x 260	5,1	20 cases	3 700 381 503 240	10 ploybags of 50 units	3 700 381 503 257
70.703				3 700 381 504 803		3 700 381 504 810