

**20.100****20.200****20.300****PLASTIC COATED & ELASTICATED SLEEVE**  
**PE + FOOD CONTACT****FEATURES & BENEFITS**

- PE elasticated sleeves ensure forearm protection against various types of dust
- Exists in different colors : ability to use a color code for tracking
- Food contact



20.100

Blanc



20.200

Transparent



20.300

Bleu

**TECHNICAL DESCRIPTION**

Product name:	Plastic coated & elasticated sleeve
Product type :	Single use, non sterile
Material :	Polypropylene
Confection :	2 rubber elastics
Color :	White, transparent or blue
Grammage/Thickness :	20 µm
Dimensions :	40 x 20 cm
Weight :	2,8 g
Country of origin :	China

**INDUSTRIES**

Industrial / Restaurant &amp; catering / Community services / Hygiene

**STORAGE CONDITIONS**

Normal conditions of conservation &amp; storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

**CERTIFICATIONS & STANDARDS**

Manufacturing site : ISO 9001.

Meets the requirements of Reg.(EC)1935/2004 and Reg. (EU)10/2011. This product is suitable for food contact with all food types.

20.100

20.200

20.300

## PLASTIC COATED & ELASTICATED SLEEVE PE + FOOD CONTACT

### LOGISTIC INFORMATION

OUTERCASE SPECIFICATIONS					SPECIFICATIONS CDT	
Article	SIZE (mm)	WEIGHT (kg)	Qty / pallet	EAN CODE	CDT	CODE EAN
20.100	310 x 260 x 210	3,5	60 outercases (10 layers of 6 outercases)	3 700 381 500 256	10 polybags of 100 units	3 700 381 501 581
20.200				3 700 381 500 263		3 700 381 501 598
20.300				3 700 381 500 270		3 700 381 501 604