

70.011 70.012 PP BOUFFANT CAP - 12 g/m²
70.011S 70.013 SINGLE ELASTIC - Ø 53 CM

FEATURES & BENEFITS

- Lightweight for improved user comfort.
- Ability to use a color code for tracking.
- It adapts to different head sizes.
- Without latex.
- Can be used to control the environment for HACCP processes and various industries.



TECHNICAL DESCRIPTION

Product name : PP Bouffant Cap 12 g/m²
 Single elastic Ø 53 cm

Product type : Single use, non-sterile

Material : Polypropylene

Confection : Sewn and elasticated

Couleur : White / Blue / Green





Grammage : 12 g/m² (± 10%)

Dimensions : Diameter 53 cm (± 1 cm)

Unit weight : 2 g (± 10%)

Packaging : Bag of 100 units

Origin : China

| Item | Couleur | Conditionnement |
|---------|--|--------------------|
| 70.011 |  White | Box of 100 pcs |
| 70.011S |  White | Polybag of 100 pcs |
| 70.012 |  Blue | Box of 100 pcs |
| 70.013 |  Green | Box of 100 pcs |

INDUSTRIES

Industry / Catering / Food / Hygiene

CERTIFICATIONS & STANDARDS

Manufacturing site : ISO 9001

STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

LOGISTIC INFORMATIONS

| SPECIFICATIONS OUTERCASE | | | | | SPECIFICATIONS CDT | |
|--------------------------|-----------------|-------------|---------------|-------------------|--------------------------|-------------------|
| Article | SIZE (mm) | WEIGHT (kg) | Article | EAN CODE | CDT | EAN CODE |
| 70.011 | 460 x 250 x 250 | 2,2 | 54 outercases | 3 700 381 501 109 | 5 boxes of 100 units | 3 700 381 502 557 |
| 70.011S | 340 x 260 x 260 | 2,6 | | 3 700 381 500 881 | 10 polybags of 100 units | 3 700 381 503 097 |
| 70.012 | 460 x 250 x 250 | 2,2 | | 3 700 381 501 116 | 5 boxes of 100 units | 3 700 381 502 564 |
| 70.013 | | | | 3 700 381 501 123 | | 3 700 381 502 571 |