

10.006

10.005V

10.007

### Shoecovers with antislip sole

Medical Device Class I

#### CARACTÉRISTIQUES ET AVANTAGES

- The PE sole ensures good mechanical strength while providing excellent grip on the ground.
- Their anatomical shape allows them to adapt to all types of footwear, and to different body shapes when used without shoes, particularly in medical environments.






#### Recommended uses

- In industrial environments on dry or wet floors
- In medical environments

#### PRODUCT DETAILS

Product name:	Kolmi® Skidguard
Product type:	Single use, non sterile
Composition :	Spunbond polypropylene + polyethylene
Spunbond grammage :	30 g/m <sup>2</sup>
Sole thickness :	60 µm
Manufacturing :	Gluing
Size :	160 x 400 mm (Size XL)
Unit weight:	11 g
Packaging :	Case of 400 ou 200 pcs
Origin :	France

Item	Sole color
10.006	 Blue
10.005V	 Green
10.007	 White

#### CERTIFICATIONS & STANDARDS

Meet requirements of Regulation (EU) 2017/745 on Medical Devices.  
Manufacturing site: ISO 9001 and ISO 13485 certified.

#### MICROBIOLOGICAL INFORMATION

Evaluation of initial microbial contamination according to standards ISO 11737:2018 + A1:2021 upon request.  
Additional microbiological tests: ASR, E. coli, staphylococci, available upon request.

10.006  
10.005V  
10.007

## Shoecovers with antislip sole

Medical Device Class I

### NORMAL STORAGE CONDITIONS

Must not be exposed to moisture or sunlight, and should be stored in a cool or room temperature environment.

Product shelf life: 5 years.

### PRECAUTIONS

The device must only be used on intact skin.

Reuse or prolonged use may cause infection or cross-contamination.

After use, comply with the applicable national regulations for disposal of the device.

### CONDITIONS DE STOCKAGE

Normal conditions of conversation & storage: not to be exposed to moisture and sunlight, must be stored at room temperature.

Shelf life of the product: 5 years.

### LOGISTICS



Outer case  
spec.

Item	Packaging Type	Qty pcs	Size (mm)	Gross weight (kg)	QTY/pallet
10.006		400	380 x 400 x 400	4,8	24 (4 layers of 6 cases)
10.007	Case with plastic bag	400	380 x 400 x 400	4,8	24 (4 layers of 6 cases)
10.005V		200	300 x 400 x 400	3	30 (6 layers of 5 cases)

### EAN Code

Item	Case
10.006	3 700 381 500 164
10.007	3 700 381 500 171
10.005V	3 700 381 503 011