

# Technical Data Sheet

## PP CLIP CAP – 10g/m<sup>2</sup>

Medical Device Class I

### FEATURES & BENEFITS

- 100% Polypropylene spunbonded nonwoven
- Gluing and ultrasound assembly
- LATEX-FREE synthetic elastic under the hem for better comfort
- Soft and very resistant non-woven
- Very breathable for better comfort
- Extensible - suitable for all hair styles



### RECOMMENDED USE

- Can be used for nursing staff or for patients in the operating room and in the care departments
- Can also be used in a pharmaceutical, industrial, food industry environments

### PRODUCT DETAILS

Product name : Kolmi® – PP Clip Cap 10g/m<sup>2</sup>  
 Product type : Single use, non-sterile  
 Material : Polypropylene  
 Size : Length : 44,5 cm (± 0,5cm)  
 Width : 48 cm (± 1 cm)  
 Unit weight : 2,4g (±10%)  
 Manufacturing country : FRANCE 

| Item    |                                                                                       | Color |
|---------|---------------------------------------------------------------------------------------|-------|
| C12004  |  | White |
| C12103  |  | Blue  |
| C12513A |  | Blue  |

### MECHANICAL INFORMATION

Mesurements done by independent laboratories  
 (measures are available upon request)

| Test             | Standard   | Laboratory |                      | Average | Unit | Tolerance |
|------------------|------------|------------|----------------------|---------|------|-----------|
| Grammage         | ISO 9073-1 | Cetelor    |                      | 10      | g/m² | ±10%      |
| Tensile strength | ISO 9073-3 | Cetelor    | Tensile strength SM: | 21,3    | N    | ±10%      |
|                  |            |            | Tensile strength ST: | 8,3     | N    | ±10%      |
| Elongation       | ISO 9073-3 |            | Elongation SM:       | 36,8    | %    | ±10%      |
|                  |            |            | Elongation ST:       | 52,6    | %    | ±10%      |

# PP CLIP CAP - 10g/m<sup>2</sup>

## Medical Device Class I

### MICROBIOLOGICAL INFORMATION

Assessment of initial microbial contamination according to standards NF EN ISO 11737: 2018+A1:2021.  
Additional microbiological controls: ASR, E.coli, Staphylococci, available upon request.

### CERTIFICATIONS & STANDARDS

Meets the requirements of Regulation (EU) 2017/745 on medical devices.  
Manufacturing site : ISO 9001 and ISO 13485 certified.

### PRECAUTIONS

The product must be used on intact skin only.  
The reuse or extended use may lead to infection or cross-contamination.  
Please refer to your national legislation on the subject of disposal of this product once it has been used.

### STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be exposed to moisture and sunlight, must be stored in a cool environment or at room temperature.  
Shelf life of the product: 5 years.

### MARKING ON PACKAGING



\*Do not use if package is damaged

### LOGISTICS INFORMATION

| OUTERBOX SPECIFICATIONS |                 |             |                                        |                   | INNERCASE SPECIFICATIONS |                               |                   |
|-------------------------|-----------------|-------------|----------------------------------------|-------------------|--------------------------|-------------------------------|-------------------|
| Item                    | Size (mm)       | Weight (kg) | Qty/ pallet                            | Barcode           | Size (mm)                | Qty/carton                    | Barcode           |
| C12004                  | 280 x 260 x 190 | 0,5         | 96 cartons<br>(8 layers of 12 cartons) | 3 662 036 006 059 | /                        | Dispenser box<br>of 300 units | /                 |
| C12103                  | 280 x 260 x 190 | 0,5         | 96 cartons<br>(8 layers of 12 cartons) | 3 662 036 006 066 | /                        | Dispenser box<br>of 300 units | /                 |
| C12513A                 | 800 x 400 x 240 | 3,4         | 21 cartons<br>(7 layers of 3 cartons)  | 3 662 036 006 080 | /                        | 12 polybags of<br>150 unités  | 3 662 036 006 103 |