

# TECHNICAL DATA SHEET

## BLUE CPE SHOECOVERS

Item code : H10100



### FEATURES & BENEFITS

- Short wear overshoes
- Very economical
- Waterproof
- Widely used in mass catering and/or by visitors
- Suitable to HACCP procedures

### PRODUCT DETAILS

Product name : CPE shoe covers  
Product type : Single use, non sterile  
Material : CPE  
Confection : Welded elastic  
Color : Blue  
Dimensions : 150 x 360 mm  
Unit weight : 3,2 g  
Packaging : 20 polybags of 100 units  
Manufacturing country : China

### CERTIFICATIONS & STANDARDS

Manufactured under ISO 9001 certified quality system

### STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

### LOGISTICS INFORMATION

#### Outer case

Item	Size mm	Gross weight kg	QTY/ pallet
H10100	400 x 300 x 380	7	30

#### Polybag

QTY / case	Size mm
20 polybags of 100 units	220 x 190 x 60

### EAN CODES

Item	Outer case	Polybag
H10100	3 662 036 019 615	3 662 036 019 622