

# Technical Data Sheet

## PP CLIP CAP - 15g/m<sup>2</sup> Medical Device Class I

### FEATURES & BENEFITS

- 100% Polypropylene spunbonded nonwoven
- Gluing and ultrasound assembly
- LATEX-FREE synthetic elastic under the hem for better comfort
- Soft and very resistant non-woven
- Very breathable for better comfort
- Extensible - suitable for all hair styles



### RECOMMENDED USE

- Can be used for nursing staff or for patients in the operating room and in the care departments
- Can also be used in a pharmaceutical, industrial, food industry environments

### PRODUCT DETAILS

Product name : Kolmi® – PP Clip Cap 15g/m<sup>2</sup>  
 Product type : Single use, non-sterile  
 Material : Polypropylene  
 Size : Length : 44,5 cm (± 0,5cm)  
 Width : 48 cm (±1 cm)  
 Unit weight : 3,4g (±10%)  
 Manufacturing country : FRANCE

Item	Color
C12001	White
C12101	Blue
C12201	Green

### MECHANICAL INFORMATION

Mesurements done by independent laboratories  
(measures are available upon request)

Test	Standard	Laboratory		Average	Unit	Tolerance
Grammage	ISO 9073-1	Cetelor		15	g/m <sup>2</sup>	±10%
Tensile strength	ISO 9073-3	Cetelor	Tensile strength SM:	34,1	N	±10%
			Tensile strength ST:	15,6	N	±10%
Elongation	ISO 9073-3		Elongation SM:	44,5	%	±10%
			Elongation ST:	51,6	%	±10%

# PP CLIP CAP - 15g/m<sup>2</sup>

## Medical Device Class I

### MICROBIOLOGICAL INFORMATION

Assessment of initial microbial contamination according to standards NF EN ISO 11737: 2018+A1:2021.  
Additional microbiological controls: ASR, E.coli, Staphylococci, available upon request.

### CERTIFICATIONS & STANDARDS

Meets the requirements of Regulation (EU) 2017/745 on medical devices.  
Manufacturing site : ISO 9001 and ISO 13485 certified.

### PRECAUTIONS

The product must be used on intact skin only.  
The reuse or extended use may lead to infection or cross-contamination.  
Please refer to your national legislation on the subject of disposal of this product once it has been used.

### STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be exposed to moisture and sunlight, must be stored in a cool environment or at room temperature.  
Shelf life of the product: 5 years.

### MARKING ON PACKAGING



\*Do not use if package is damaged

### LOGISTICS INFORMATION

#### OUTERBOX SPECIFICATIONS

Item	Size (mm)	Weight (kg)	Qty / pallet	Barcode
C12001	295 x 260 x 190	0,9	108 cartons (9 couches de 12 cartons)	3 662 036 000 897
C12101				3 662 036 000 910
C12201				3 662 036 000 934