

70.061
PP BOUFFANT CAP - 12 g/m²
70.062
SINGLE ELASTIC - Ø 61 CM

FEATURES & BENEFITS

- Lightweight for improved user comfort.
- Ability to use a color code for tracking.
- It adapts to different head sizes.
- Without latex.
- Can be used to control the environment for HACCP processes and various industries.



TECHNICAL DESCRIPTION

Product name : PP Bouffant Cap 12 g/m²
 Single elastic Ø 61 cm

Product type : Single use, non-sterile

Material : Polypropylene

Confection : Sewn and elasticated

Coulour : White / Blue

Grammage : 12 g/m² (± 10%)

Dimensions : Diameter 61 cm (± 1 cm)

Unit weight : 2,8 g (± 10%)

Packaging : Bag of 100 units

Origin : China

Item	Colour
70.061	 White
70.062	 Blue

INDUSTRIES

Industry / Catering / Food / Hygiene

CERTIFICATIONS & STANDARDS

Manufacturing site : ISO 9001

STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

LOGISTIC INFORMATIONS

SPECIFICATIONS OUTER CASE

Article	SIZE (mm)	WEIGHT (kg)	Qty / pallet	CODE EAN
70.061	350 x 270 x 270	3,2	54 outercases	3 700 381 503 707
70.062				3 700 381 507 644

SPECIFICATIONS CDT

CDT	CODE EAN
10 polybags of 100 units	3 700 381 503 691
	3 700 381 507 651