

**70.801**  
**70.802**

**PP BEARD COVER WITH ELASTICS**



**FEATURES & BENEFITS**

- The beard cover covers the lower part of the face.
- It guarantees a protection adapted to the control of the environment against the falls of beard hair
- Can also be used as an environmental control tool in HACCP and other industries.



**TECHNICAL DESCRIPTION**

Product name : PP beard cover with elastics  
 Product Type : Single use, non-sterile  
 Material : Polypropylene (PP)  
 Confection : Sewed and elasticated  
 Grammage : 10 g/m<sup>2</sup> (± 10%)  
 Dimensions : Flat length : 45 cm (± 1 cm) / Width : 24 cm (± 1 cm)  
 Unit weight : 1,2 g (± 10%)  
 Packaging : Polybags of 100 units  
 Origin : China

Item	Color
70.801	Blanc 
70.802	Bleu 

**INDUSTRIES**

Industry / Catering / Food / Hygiene

**CERTIFICATIONS & STANDARDS**

Manufacturing site : ISO 9001

**STORAGE CONDITIONS**

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

**LOGISTIC INFORMATIONS**

SPECIFICATIONS OUTERCASE					SPECIFICATIONS CDT	
Article	SIZE (mm)	WEIGHT (kg)	Qty / pallet	CODE EAN	CDT	CODE EAN
70.801	420 x 200 x 260	1,5	60 cases	3 700 381 501 376	10 polybags of 100 units	3 700 381 502 823
70.802				3 700 381 504 827		