

A69955

WATERPROOF PP ISOLATION GOWN

With PE coating

FEATURES & BENEFITS

- Premium protection to meet all your infection control needs
- Comfort and resistant
- Large opening for the head
- Patented « Flexneck » technology with a flexible elastic design
- Long sleeves
- Thumbloop for extra protection

Recommended use:

- Disposable Impervious Isolation Gown designed to protect health professionals and patients from cross infection



PRODUCT DETAILS

Product name:	Impervious coated PP Isolation Gown
Product type:	Single use, non-sterile
Material:	Polypropylene + Polyethylene
Confection:	Ultrasonically sealed, without collar, without pocket, Flexneck neck closure, long sleeves + thumb loops, elasticated wrists, Back closure with laces
Colour:	Blue
Grammage:	40 g/m ² (± 10%)
Unit weight:	118 g (± 10%)
Packaging:	Polybag of 10 units
Origin:	China

INDUSTRIES

Health / Medical

Item	Size	Length (± 1 cm)	Underarm width (± 1 cm)	Gown width with sleeves (includes thumb loop) (± 1 cm)	Open gown width (± 1 cm)
A69955	XL	129 cm	79 cm	235 cm	158 cm

CERTIFICATION & STANDARDS

In accordance with ANSI/AAMI PB 70 level II barrier performance for surgical gown and tested for chemotherapy drug permeation as per ASTM D6978.

STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

LOGISTICAL INFORMATION

OUTERCASE SPECIFICATION					CDT SPECIFICATIONS	
Article	SIZE (mm)	WEIGHT (kg)	Qty / pallet	EAN CODE	CDT	CODE EAN
A69955	540 x 325 x 400	13	20 outercases (5 layers of 4 outercases)	206 868 640 518 97	10 polybags of 10 units	106 868 640 518 90