

TECHNICAL DATA SHEET

PP TROUSERS & TUNIC FOR HEALTHCARE WORKERS

J40480



FEATURES & BENEFITS

- Comfortable and light, the material is highly breathable
- Its ample shape allows the user to have a good ease of movement. Without collar.
- Trousers with elasticated waistband
- Individually wrapped for hygiene and controlled distribution.
- No pockets to prevent objects from falling out

PRODUCT DETAILS

Product name : Trousers & tunic for healthcare workers
Product type : Single use, non sterile
Composition : Polypropylene
Confection : Sewn, short sleeves, no collar, no pocket, elasticated trousers at the waist
Color : Dark blue
Grammage : 48 g/m² (+/-10%)
Packaging : 60 pcs individually wrapped
Manufacturing country : China



Item	Size	Tunic length (± 1 cm)	Tunic width (± 1 cm)	Tunic width with sleeves (± 1 cm)	Trousers length (± 1 cm)	Trousers width (± 1 cm)	Ankle width (± 1 cm)	Waist width (± 1 cm)	Leg length (± 1 cm)	Unit weight (± 10%)
J40480/S	S	74 cm	57,5 cm	81 cm	103 cm	52 cm	25 cm	45 cm	67 cm	93 g
J40480/M	M	75 cm	61 cm	83 cm	106 cm	56 cm	25 cm	49 cm	69 cm	101 g
J40480/L	L	76 cm	64,5 cm	85 cm	109,5 cm	60 cm	25,5 cm	53 cm	71,5 cm	107 g
J40480/XL	XL	77 cm	68 cm	88 cm	111 cm	64 cm	26 cm	57 cm	73 cm	111 g
J40480/XXL	XXL	78 cm	71,5 cm	90 cm	114 cm	68 cm	26,5 cm	60 cm	75 cm	119 g

CERTIFICATIONS & STANDARDS

Manufactured under ISO 9001 certified quality system

STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

INFORMATIONS LOGISTIQUES

OUTERCASE

POLYBAG

Item	Size mm	Gross weight kg	QTY/pallet	CODE EAN	QTY / case	CODE EAN
J40480/S		6,5		3 662 036 019 394		3 662 036 019 394
J40480/M	460 x 280 x 450	7,2	24 cartons (4 layers of 6 cartons)	3 662 036 019 417	60 individual polybags	3 662 036 019 417
J40480/L		7,4		3 662 036 019 431		3 662 036 019 431
J40480/XL		7,7		3 662 036 019 455		3 662 036 019 455
J40480/XXL		8,2		3 662 036 019 479		3 662 036 019 479

DOMINIQUE DUTSCHER SAS