

40.170/04 PP ELASTICATED TROUSER

FEATURES & BENEFITS

- Comfort and light, the PP is very breathable.
- Wide shape allows the user to have a good comfort.
- Easy to put.
- The resistance of trouser allows the user to work on their knees without it does not get damaged.
- Belt and ankles elasticated.
- Single size.
- Single packaging for hygiene and distribution control.

TECHNICAL DESCRIPTION

Product name: PP Elasticated Trouser
 Product type: Single use, non sterile
 Material: Polypropylene
 Confection: Sewn in double points
 Hankles & waist elasticated
 Color: White
 Thickness: 45 g/m2 (±5)
 Dimensions: 117 x 132 cm
 Weight: 70 g
 Packaging: Individual polybag
 Country of origin: China



INDUSTRIES

Industrial / Restaurant & catering / Community services / Hygiene

CERTIFICATIONS & STANDARDS

Manufacturing site : ISO 9001.

STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

LOGISTIC INFORMATION

OUTER CASE SPECIFICATIONS					SPECIFICATIONS CDT	
Article	SIZE (mm)	WEIGHT (kg)	Qty / pallet	EAN CODE	CDT	CODE EAN
40.170/04	470 x 280 x 280	4,5	36 (6 layers of 6 cases)	3 700 381 503 721	50 individual polybags	3 700 381 503 738